















































































## 18 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
<b>A-G</b>	<b>Standards</b>	in <b>500 µL</b> Dilution Buffer <b>VP</b>	-
<b>KS1</b>	<b>Control Serum 1</b>	in <b>250 µL</b> Dilution Buffer <b>VP</b>	<b>≥ 1:16</b> with <b>VP</b>
<b>KS2</b>	<b>Control Serum 2</b>	in <b>250 µL</b> Dilution Buffer <b>VP</b>	<b>≥ 1:16</b> with <b>VP</b>
<b>WP</b>	<b>Washing Buffer</b>	-	<b>1:20</b> with <b>Aqua dest</b>
<b>Sample dilution:</b> with Dilution Buffer <b>VP ≥ 1:16</b>			
Before assay procedure bring all reagents to room temperature <b>20-25°C</b> .			
<b>Assay Procedure in Double Determination:</b>			
Pipette	Reagents	Position	
50 µL	Antibody Conjugate <b>AK</b>	Pipette in <b>all</b> required number of wells	
50 µL	Standard <b>A (0 ng/mL)</b>	A1/A2	
50 µL	Standard <b>B (0.1 ng/mL)</b>	B1/B2	
50 µL	Standard <b>C (0.5 ng/mL)</b>	C1/C2	
50 µL	Standard <b>D (1 ng/mL)</b>	D1/D2	
50 µL	Standard <b>E (2 ng/mL)</b>	E1/E2	
50 µL	Standard <b>F (4 ng/mL)</b>	F1/F2	
50 µL	Standard <b>G (8 ng/mL)</b>	G1/G2	
50 µL	Control Serum <b>KS 1</b> (≥ 1:16 diluted)	H1/G2	
50 µL	Control Serum <b>KS 2</b> (≥ 1:16 diluted)	A3/A4	
50 µL	Sample (≥ 1:16 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
<b>Sample Incubation with Shaking: 1 h at 20°C - 25°C, 350 rpm</b>			
5 x 300 µL	Aspirate the contents of the wells and <b>wash</b> 5 x with 300 µL each Washing Buffer <b>WP/ well</b>	In each well	
100 µL	Enzyme Conjugate <b>EK</b>	In each well	
Cover the wells with the sealing tape.			
<b>Incubation with Shaking: 30 Minutes at 20°C - 25°C, 350 rpm</b>			
5 x 300 µL	Aspirate the contents of the wells and <b>wash</b> 5 x with 300 µL each Washing Buffer <b>WP/ well</b>	In each well	
100 µL	Substrate Solution <b>S</b>	In each well	
<b>Incubation: 15 Minutes in the Dark at 20°C - 25°C</b>			
100 µL	Stopping Solution <b>SL</b>	In each well	
Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with ≥ 590 nm as reference wavelength.			