hGH Sensitiv/-e ELISA

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| | |
|--|-----------|
| TABLE OF CONTENTS | 20 |
| ENGLISH Instructions for use | 21 |
| 1 INTENDED USE | 21 |
| 2 INTRODUCTION | 21 |
| 3 ASSAY PRINCIPLE | 22 |
| 4 WARNINGS AND PRECAUTIONS | 23 |
| 5 SAMPLES | 24 |
| 6 MATERIALS | 25 |
| 7 TECHNICAL NOTES | 26 |
| 8 ASSAY PROCEDURE | 27 |
| 9 QUALITY CONTROL | 28 |
| 10 EVALUATION OF RESULTS | 28 |
| 11 PERFORMANCE CHARACTERISTICS | 31 |
| 12 COMPARISON STUDIES | 34 |
| Instructions for use for scientific application | 35 |
| 13 SCIENTIFIC APPLICATION | 35 |
| 14 LITERATUR / REFERENCES | 36 |
| Internationale Assay Description | 40 |

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| hGH ELISA | 96 Determinations |
|-------------------------------|--|
| ○ | DE/CA40/00809/15-1 |
| Principle of the test | Sandwich ELISA |
| Duration (incubation period) | 3.25 h |
| Antibody Conjugate | ready for use |
| Enzyme Conjugate | ready for use |
| Buffer and Substrate | ready for use |
| Reference material | 2. Internationaler Standard WHO/ NIBSC 98/574, recombinant hGH |
| Standards | 5 single standards: 0.05 - 1 ng/mL, lyophilized, recombinant hGH |
| Assay Range | 0.0115 - 26 ng/mL |
| Control | 1 control serum, freeze-dried |
| Sample | human serum / plasma |
| Required sample volume | 10 µL |
| Sample dilution | 1:26 |
| Analytical sensitivity | 0.0115 ng/mL |
| Intra- / Inter-Assay Variance | < 10% |

1 INTENDED USE

The ELISA E022 is intended to be used for quantitative measurement of human growth hormone in human serum and plasma samples.

2 INTRODUCTION

The endocrine system of human Growth Hormone (hGH), also named Somatropin, is characterized by an extreme complexity. hGH is the product of the GH-1 gene located on chromosome 17 and expressed in pituitary cells. 80% of the hGH is a non-glycosylated 22 kDa protein consisting of 191 amino acids. About 20% is a variant form of 20 kDa resulting from alternative splicing. Additionally, several smaller variants can be found in circulation as well as translational modified proteins and different degrees of protein aggregation. Bioactivity of Growth Hormone is regulated by a specific binding protein (GHBP) formed by the extra cellular part of the cellular transmembran GH-receptor. These modifications allow a tight control of the half-life period hGH and of its bioactivity.

Not only synthesis and posttranslational modification but also secretion of hGH is tightly regulated. Spontaneous pulsatile secretion takes place with a single pulse every three hours and a maximal secretion during night's sleep. Several different stimuli as physiologic stress or hypoglycaemia result in additional hGH secretion, induced by the hypothalamic hormones Somatostatin and GH-Releasing Hormone (GHRH). Age, sexual steroids, nutritional status, illness and emotions influence the amount of secreted hGH. Because of the multitude of influencing factors the normal quantitative secretion is not known (see References 2-5).

Physiological functions are partially exerted by Insulin-like Growth Factors (IGFs). In children and adolescent the hGH system is the main regulator of growth. If the hGH system fails totally, human growth will end at 120 cm. Beside regulation of growth hGH exerts an anabolic effect on muscle and connective tissue as well as on bone and different other organs (heart, intestine). Further hGH was proved to have a lipolytic effect.

Growth Hormone pathology is characterized by extreme high or extreme low hGH secretion. During childhood it is the Growth Hormone deficiency congenital or acquired, which leads to microsomia. For diagnosis of Growth Hormone deficiency an hGH stimulation test has to be done or the spontaneous excretion must be investigated. The therapy consists of substitution of endogenous Growth Hormone by recombinant hGH resulting in normalization of growth.

In adulthood hGH deficiency is mostly caused by pituitary adenoma (and their surgical excision). hGH deficiency is associated with adipositas, muscle dystrophy, arteriosclerosis, osteoporosis, adynamia. Substitutional therapy is a well-known, approved and efficient therapy of severe Growth Hormone deficiency in adulthood. Therapeutic success is directly as well as indirectly proved by measurement of IGF in serum.

Excessive hGH secretion, mostly caused by pituitary adenoma, results in childhood in gigantism, in adulthood in acromegaly, leading to enlarged extremities, diabetes, heart insufficiency and tumor growth. Surgical excision of the adenoma is the therapy of choice. If tumor excision is not possible or incomplete, a medicinal therapy with somatostatin preparation will be conducted, resulting in inhibition of hGH production. Alternatively hGH analoga (e.g. Pegvisomat) are used to block the hGH receptor and thereby inhibit action of endogenous hGH. Measurement of human Growth Hormone (hGH, Somatropin) is done for diagnostic of Growth Hormone deficiency or Growth Hormone excess (arcomegaly). During medicinal and/or after surgical therapy of arcomegaly Growth Hormone (and IGF-I) measurement is used for therapy control.

3 ASSAY PRINCIPLE

The Mediagnost hGH SENSITIVE ELISA E022 is a so-called sandwich-assay. It utilizes a specific, high affinity polyclonal rabbit antiserum coated on the wells of a microtiter plate. The hGH in the samples binds quantitatively to the immobilized antiserum. In the following step, the biotinylated antibody in turn binds hGH. After washing, a streptavidin-peroxidase-enzyme conjugate will be added, which will bind highly specific to the biotin of the antibody and will catalyze the substrate to change the color quantitatively depending on the hGH level of the sample.

4 WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Control Serum KS**

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

Reagents A-E, AK, EK, VP, WP

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

| | |
|-----------|---|
| H317 | May cause an allergic skin reaction. |
| P280 | Wear protective gloves/ protective clothing/ eye protection/ face protection. |
| P272 | Contaminated work clothing should not be allowed out of the workplace. |
| P261 | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. |
| P333+P313 | If skin irritation or rash occurs: Get medical advice/ attention. |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water. |
| P501 | Dispose of contents/ container in accordance with local/ regional/ national/ international regulations. |

Substrate Solution (S)

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine (<0.05%)

| | |
|------------|---|
| H315 | Causes skin irritation. |
| H319 | Causes serious eye irritation. |
| H335 | May cause respiratory irritation. |
| P261 | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes. |
| P338 | Remove contact lenses, if present and easy to do. Continue rinsing. |

Stopping Solution (SL)

The Stopping solution contains 0.2 M acid sulphur acid (H₂SO₄)

| | |
|------------|---|
| H290 | May be corrosive to metals. |
| H314 | Causes severe skin burns and eye damage. |
| P280 | Wear protective gloves/ protective clothing/ eye protection/ face protection. |
| P301+P330+ | IF SWALLOWED: rinse mouth. |
| P331 | Do NOT induce vomiting. |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes. |
| P338 | Remove contact lenses, if present and easy to do. Continue rinsing. |
| P309+P310 | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |

4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

5.1 Sample type

Serum and Plasma

Serum and Heparin/EDTA Plasma yield comparable values. The hGH levels are reduced in citrate plasma samples, because of the relatively high amount of anticoagulant.

5.2 Specimen collection

Human GH is secreted pulsatile during the day/night. Therefore, in clinical application stimulation test are used to measure peak GH concentrations.

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

5.3 Required sample volume: 10 µL

5.4 Sample stability

In firmly closable sample vials

- Storage at 20-25°C: 3 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles max. 5

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized. 5 Freezing-Thawing showed no effect on samples.

5.5 Interference

Triglyceride and bilirubin in the sample do not interfere to a concentration of 100 mg/mL and 200 µg/mL. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

5.6 Sample dilution

- Dilution: **1:26** with Dilution Buffer **VP**
- Pipette **250 µL** Dilution Buffer **VP** in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series); add **10 µL sample** (dilution 1:26). After mixing use 2 x 100 µL of this dilution in the assay.
- Sample stability after dilution of the sample: maximum 1 hour at 20-25°C.

6 MATERIALS

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the standard curve.

| | | |
|---|--|---------------------|
| MTP | Microtiter plate , ready for use, coated with rabbit-anti-hGH-antibody. Wells are separately breakable. | (8x12) wells |
| A-E | Standards , lyophilized, (recombinant human hGH), concentrations are given on vial labels and on quality certificate in ng/mL. | 5 x 750 µL |
| KS | Control Serum , lyophilised, (human serum), concentration is given on quality certificate in ng/mL. | 1 x 500 µL |
| AK | Antibody Conjugate , ready for use, contains rabbit biotinylated anti-hGH antibody. | 1 x 12 mL |
| EK | Enzyme Conjugate , ready for use, contains HRP (Horseradish-Peroxidase)-labelled Streptavidin. | 1 x 12 mL |
| VP | Dilution Buffer , ready for use | 1 x 120 mL |
| WP | Washing Buffer , 20-fold concentrated solution | 1 x 50 mL |
| S | Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised H ₂ O ₂ Tetramethylbencidine. | 1 x 12 mL |
| SL | Stopping Solution , ready for use, 0.2 M sulphuric acid. | 1 x 12 mL |
| - | Sealing Tape , for covering the microtiter plate . | 2 x |
|  | Instructions for use | 1 x |
| -- | Quality Certificate | 1 x |

6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WP (A. dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** standards **A-E** and Control Serum **KS** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WP** is 4 weeks stable at 2-8°C

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Standards **A – E** and Control **KS** are reconstituted with the Dilution Buffer **VP**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Control **KS** with the Dilution Buffer **VP** in the same ratio (1:26) as the sample. The required volume of Washing Buffer **WP** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

Assay Procedure

When performing the assay, Blank, Standards **A-E**, Control **KS** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **AK** and the Enzyme Conjugate **EK** as well as the succeeding Substrate Solution **S** should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution **SL** should be added to the plate in the same order as Substrate Solution **S**. All determinations (Blank, Standards **A-E**, Control **KS** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate Solution **S**, stabilised H₂O₂-Tetramethylbenzidine, is photosensitive—store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WP** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8 ASSAY PROCEDURE

| Preparation of reagents | | Reconstitution: | Dilution |
|--|--|---|----------------------|
| A-E | Standards | in 750 µL Dilution Buffer VP | - |
| KS | Control Serum | in 500 µL Dilution Buffer VP | 1:26 with VP |
| WP | Washing Buffer | - | 1:20 with Aqua dest. |
| Sample dilution: with Dilution Buffer VP 1:26 | | | |
| Before assay procedure bring all reagents to room temperature 20-25°C. | | | |
| Assay Procedure in Double Determination: | | | |
| Pipette | Reagents | Position | |
| 100 µL | Dilution Buffer VP (Blank) | A1/A2 | |
| 100 µL | Standard A (0.05 ng/mL) | B1/B2 | |
| 100 µL | Standard B (0.15 ng/mL) | C1/C2 | |
| 100 µL | Standard C (0.30 ng/mL) | D1/D2 | |
| 100 µL | Standard D (0.6 ng/mL) | E1/E2 | |
| 100 µL | Standard E (1.0 ng/mL) | F1/F2 | |
| 100 µL | Control Serum KS (1:26 diluted) | G1/G2 | |
| 100 µL | Sample (1:26 diluted) | in the rest of the wells according the requirements | |
| Cover the wells with the sealing tape. | | | |
| Sample Incubation: 2 h at 20-25°C, 350 rpm | | | |
| 5 x 300 µL | Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well | In each well | |
| 100 µL | Antibody Conjugate AK | In each well | |
| Cover the wells with the sealing tape. | | | |
| Incubation: 30 Minutes at 20-25°C, 350 rpm | | | |
| 100 µL | Enzyme Conjugate EK, without washing the wells (!) – add to the previously pipetted Antibody Conjugate AK-solution <u>thereto</u> mix shortly through cautious tapping on the side of the MTP. Attention: high filled volume of the wells! | In each well | |
| Cover the wells with the sealing tape. | | | |
| Incubation: 30 Minutes at 20-25°C, without shaking | | | |
| 5 x 300 µL | Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well | In each well | |
| 100 µL | Substrate Solution S | In each well | |
| Incubation: 15 Minutes in the Dark at 20-25°C | | | |
| 100 µL | Stopping Solution SL | In each well | |
| | Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength. | | |

9 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All standards and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

9.1 Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of standard E should be above 1.00.

Samples, which yield higher absorbance values than Standard E, should be re-tested with a higher dilution.

10 EVALUATION OF RESULTS

10.1 Establishing of the standard curve

The 2nd International Standard for hGH, NIBSC Code 98/574 (6), was used as standard material. This was defined in an international study in the year 2001 with 3 International units per mg Protein (3 IU/mg). The exclusive application of this standard material is recommended in line with the current standardisation efforts for hGH Immunoassays. (7,8)

| Standard | A | B | C | D | E |
|-------------|------|------|------|------|------|
| ng/mL | 0.05 | 0.15 | 0.30 | 0.60 | 1.0 |
| pg/mL | 50 | 150 | 300 | 600 | 1000 |
| μ IU/mL | 0.15 | 0.45 | 0.9 | 1.8 | 3.0 |

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples and standards
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The hGH concentration in ng/mL (or pg/mL, or μ IU/mL, according the chosen unit for the standards) of the samples can be calculated by **multiplication** with the respective **dilution factor**.

10.2 Example of a typical standard curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

| | Blank | A | B | C | D | E |
|----------------------------|-------|--------|--------|--------|--------|-------|
| ng/mL | 0.0 | 0.05 | 0.15 | 0.3 | 0.6 | 1.0 |
| OD _(450-620 nm) | 0.08 | 0.2955 | 0.5885 | 1.1045 | 1.9935 | 2.761 |


Figure 1: Exemplary standard curve

The exemplary shown standard curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a standard curve for each test you conduct!

10.3 Exemplary calculation of GH concentrations

Sample dilution: 1:26

| | |
|------------------------------------|------|
| Measured extinction of your sample | 0.25 |
| Measured extinction of the blank | 0.08 |

Your measurement programm will calculate the hGH concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 3rd degree).

In this exemplary case the following equation is solved by the program to calculate the hGH concentration in the sample:

$$0.17 = -1.2868x^3 + 0.8619x^2 + 3.0523x + 0.0532$$

$$0.035 = x$$

If the dilution factor (**1:26**) is taken into account the hGH concentration of the undiluted sample is

$$0.035 \text{ ng/mL} \times 26 = 0.91 \text{ ng/mL}$$

10.4 Interpretation of results

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

Usually the cut-off value is determined as a maximal peak of growth hormone secretion in at least 2 independent stimulation assays (e.g. insulin or arginine stimulation). If the used test system is calibrated to WHO standard 98/574 a secretion peak of less than **8 ng/mL** indicates a possible growth hormone deficiency. But as growth hormone secretion is continuous between normal and pathological any cut-off is only a non-binding benchmark.

It is recommended to consider the international and national guidelines for diagnosis and treatment of growth hormone deficiency / acromegaly.

10.5 Limitation of procedure

The Mediagnost sensitive human Growth-Hormone ELISA, E022 is based on polyclonal rabbit antibodies. Generally, this technique is sensible to heterophilic antibodies in the sample. The influence of heterophilic antibodies is reduced by assay design, but cannot be excluded completely.

Interference of several physiological and pharmaceutical substances has been tested for the indicated concentrations. Higher concentrations or other substances may interfere with the measurement.

If a cut-off of 8 µg/L is used in the growth hormone stimulation test for diagnosis of growth hormone deficiency and the test system is calibrated against the WHO 98/574 standard, literature states a diagnostic sensitivity and specificity about 80%, if all clinical preconditions for testing are fulfilled (10).

10.6 Reference values

As growth hormone is secreted pulsatile mainly enduring the night sleep valid normal values can hardly be determined. Standard procedures are arginine or insulin stimulation tests, after injection of stimulating substance growth hormone concentration is measured over a period of time. We investigated hGH serum concentration of 104 healthy blood donors in the age of 18-69 years without any stimulation and undefined sampling time.

Table 1: Reference values hGH serum concentration of 104 healthy blood donors

| | female | male |
|--------------------------------------|--------|------|
| number | 54 | 50 |
| median [ng/mL] | 0.81 | 0.28 |
| minimal concentration [ng/mL] | 0.19 | 0.15 |
| maximal concentration [ng/mL] | 10.45 | 4.34 |

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests. Furthermore, we recommend that each laboratory determine its own range for the population tested.

11 PERFORMANCE CHARACTERISTICS

11.1 Sensitivity

Sensitivity was assessed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. The analytical sensitivity of the E022 is 0.0115 µg/L.

11.2 Specificity

Cross reactivity with recombinant human Prolactin has been tested and no significant signal was measured in an enriched serum sample containing 200 µg/L Prolactin.

Further, Pegvisomant (trade name Somavert), a growth hormone analogue and drug used in acromegaly therapy, has been tested for cross-reactivity in assay buffer in different concentrations. Here no significant influence of Pegvisomant was detected (Table 2).

Table 2: To determine the **specificity** Pegvisomant (trade name Somavert) was diluted in the Dilution Buffer (VP) at the indicated concentrations and used as a sample in E022. The cross-reactivity of Pegvisomant is presented.

| Pegvisomant [mg/L] | Concentration measured in E022 [mg/L] | % Cross-reactivity |
|--------------------|---------------------------------------|--------------------|
| 100 | 0.0114 | 0.0114 |
| 10 | 0.00845 | 0.0845 |
| 1 | 0.00436 | 0.436 |
| 0.1 | 0.00103 | 1.03 |
| 0.01 | 0.00017 | 1.7 |
| 0.001 | 0.00006 | 6 |
| 0.0001 | 0.00004 | 40 |

However, a study (9) showed that, after the addition of Pegvisomant (100 mg/L) to serum samples with enriched with hGH (2.6 and 10.1 µg/L), the measured concentrations of hGH were 154% and 108% of the expected values.

11.3 Precision

Intra-Assay-Variation

One sample has been measured 14 times in the same assay. The results are shown in Table 3. The measured coefficient of variation (CV) is 5.46%. Intra assay variance has also been evaluated externally (9), two serum samples with 0.45 and 5.94 µg/L hGH were measured 10 times within the same assay. The resulting coefficients of variation were 3.65% and 2.16%

Table 3: Intra-Assay Variation

| | Number of determinations | Mean value (µg/L) | Standard deviation (µg/L) | VC (%) |
|----------|--------------------------|-------------------|---------------------------|--------|
| Sample 1 | 16 | 2.41 | 0.19 | 7.99 |
| Sample 2 | 16 | 5.84 | 0.27 | 4.70 |
| Sample 3 | 16 | 14.98 | 0.55 | 3.70 |

Inter Assay Variance

Serum samples were measured in independent assays. On average the coefficient of variation was 4.34%. Results are shown in detail in table 4. Here also externally acquired data are available: The mean coefficient of variation for inter-assay variance at 2.39; 5.37 and 14.33 µg/L hGH was 5.98%; 3.93% and 3.12%, respectively (9).

Table 4: Inter-Assay Variation

| | Number of single determinations | Mean value (µg/L) | Standard deviation (µg/L) | VC (%) |
|----------|---------------------------------|-------------------|---------------------------|--------|
| Sample 1 | 14 | 5.37 | 0.21 | 3.93 |
| Sample 2 | 10 | 2.39 | 0.14 | 5.98 |
| Sample 3 | 11 | 14.33 | 0.45 | 3.12 |

11.4 Linearity

Linearity of the E022 was tested by dilution of 2 different serum samples. The samples were diluted in the range of 1:10 to 1:76800. Linearity of sample dilution was shown by linear regression in the dilution range of 1:10 - 1:9600 (Figure 2).


Figure 2: Linearity sample dilution. Two samples with different hGH concentrations were diluted a) 1:600 - 1:9600 b) 1:10 - 1:320. The hGH concentrations were recalculated.

11.5 Recovery and Accuracy

Recombinant human Growth Hormone (NIBSC 98/574) was added in different amounts to human serum. The hGH content of the so enriched samples was measured and recovery in comparison to enriched Dilution Buffer (VP) calculated. Results are shown in table 5.

Table 5: Recovery of recombinant human GH in Serum

| NIBSC Rec. hGH ng/mL | VP | | Serum 1 PAA539 | | Serum 2 PAA 574 | |
|----------------------------|-------|-------|-------------------|-------|--------------------|------|
| | ng/mL | % | ng/mL | % | ng/mL | % |
| 20 | 19.34 | 96.7 | 17.51 | 83.2 | 17.53 | 81.9 |
| 10 | 9.81 | 98.1 | 9.28 | 84.0 | 9.18 | 80.4 |
| 5 | 5.34 | 106.7 | 4.99 | 82.67 | 5.9 | 92.0 |
| 0 | 0 | - | 1.04 | - | 1.41 | - |

11.6 Interference

Interference of bilirubin and triglycerides has been tested in [9]. Here neither bilirubin (up to 200 mg/L) nor triglycerides (up to 100 g/L) showed a significant interference with hGH measurement (Table 6). The authors also tested the influence of growth hormone binding protein up to 10 µg/L on hGH measurement and haven't seen a significant effect (mean recovery 98%). Mediagnost data are shown in Table 7.

Table 6: Interference of GHBP on GH measurement

| | hGH [µg/L] | | |
|------------------|------------|----|----|
| | 2 | 8 | 20 |
| GHBP [µg/L] | 1 | 1 | 1 |
| hGH Recovery [%] | 95 | 95 | 99 |
| GHBP [µg/L] | 5 | 5 | 5 |
| hGH Recovery [%] | 87 | 95 | 97 |
| GHBP [µg/L] | 10 | 10 | 10 |
| hGH Recovery [%] | 87 | 93 | 95 |

Table 7: Interference of Bilirubin and Triglycerides on GH measurement

| Bilirubin [mg/L] | hGH Recovery [%] | Triglycerides [g/L] | hGH Recovery [%] |
|---------------------|---------------------|------------------------|---------------------|
| 25 | 111 | 12.5 | 89 |
| 50 | 116 | 25 | 109 |
| 100 | 112 | 50 | 85 |
| 200 | 108 | 100 | 110 |

12 COMPARISON STUDIES

Langkamp et al Growth Horm IGF Res 18 (2008) 526-532 compared an in-house assay used for hGH measurements for years with the Mediagnost E022 (9). Results show a very good comparability and thus cut-off values established by the in house-assay can also be used with the Mediagnost E022 (see Figure 3).


Figure 3: Assay Comparison

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Instructions for use for scientific application

13 SCIENTIFIC APPLICATION

In addition to serum and plasma samples hGH can be determined in other human body fluids and in cell culture supernatants of various human cell lines for research purposes. Additionally this hGH Sensitive ELISA E022 can be used for measurement of hGH in filter paper samples, which has been validated externally (9).

13.1 Samples suitable for scientific application

Serum, plasma, saliva, whole blood spots on filter paper, cell culture supernatant of different human cell lines and **urine** (e.g. hGH urine–excretion test)

The recommended dilution for serum and plasma samples in Dilution Buffer **VP: (1:26)**.

In the other samples, the hGH levels can vary considerable, the optimal dilution must be found out by the customer.

13.2 Species Cross-Reactivity

Several commercially available animal sera have been tested as samples in different dilutions (1:5 or 1:26) in this assay.

No signal was detected in serum of the following species: donkey, dog, goat, guinea pig, horse, rat, mouse, rabbit, sheep, cat, chicken and horse.

Whether this obvious non-reactivity is species specific should be assessed individually by each customer. We remind each customer that hGH secretion is pulsatile and thus commercially available animal serum samples may not be taken at the ideal daytime. An external laboratory was able to find good measurable signals in bovine serum.

14 LITERATUR / REFERENCES

- 1) Hauffa BP, Lehmann N, Bettendorf M, Mehls O, Dörr H-G, Partsch C-J, Schwarz HP, Stahnke N, Steinkamp H, Said E, Sander S, Ranke MB and participating Members of the German KIGS/IGLU Study Group (2004); Central reassessment of growth hormone concentrations measured at local treatment centers in children with impaired growth: consequences for patient management. *European Journal of Endocrinology* 150: 291 – 297
- 2) Baumann G (1991) Growth hormone heterogeneity: Genes, Isohormones, variants and binding protein. *Endocr Rev* 12:424-449
- 3) Clemmons DR, Chihara K, Freda PU, Ho KKY, Klibanski A, Melmed S, Shalet S, Strasburger CJ, Trainer PJ, Thorner MO (2003) Optimizing control of acromegaly: Integrating a growth hormone receptor antagonist into the treatment algorithm. *J Clin Endocrinol Metab* 88:4759-4767
- 4) Ranke, MB, Örskov H, Bristow AF, Seth J, Baumann (1999) Consensus on how to measure growth hormone in serum. *Horm res* 51 (Suppl.1):27-29
- 5) Ranke MB (2003) Diagnosis of growth hormone deficiency and growth hormone stimulation tests. In: *Diagnostics of Endocrine Function in Children and Adolescents* (3rd ed., Michael B Ranke, Herausg.) Basel, Karger pp 107-128
- 6) Adresse NIBSC: National Institute for Biological Standards and Controls, Blanche Lane, South Mimms, Potters Bar, Hertfordshire EN6 3QG, Great Britain
- 7) Strasburger CJ (2004) Taking one step at a time *Clin. Endocrinology* 60, 540
- 8) Wieringa GE, Barth JH, Trainer PJ (2004) Growth Hormone assay standardization: a biased view ? *Clin. Endocrinology* 60, 538 – 539
- 9) Langkamp M, Weber K, Ranke MB. Human growth hormone measurement by means of a sensitive ELISA of whole blood spots on filter paper. *Growth Horm IGF Res.* 2008 Dec;18(6):526-32. Epub 2008 Jun 24
- 10) S3-Leitlinie. Diagnostik des Wachstumshormonmangels im Kindes- und Jugendalter *Endokrinologie Informationen* 33 (2009) 3

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