

12 SUMMARY - MEDIAGNOST MOUSE-/RAT-GH ELISA

Preparation of reagents:		Reconstitution:	Dilution:
A-G	Standards	in 1 mL Dilution Buffer VP	-
KS1	Control Serum 1	in 150 µL Dilution Buffer VP	1:5 with Dilution Buffer VP
KS2	Control Serum 2	in 150 µL Dilution Buffer VP	1:5 with Dilution Buffer VP
WP	Washing Buffer	-	1:20 with Aqua dest.
Sample and Control Sera KS1 and KS2: dilute 1:5 with Dilution Buffer VP, mix immediately, incubate max. 60 min. Use 100 µL for each well in the assay.			
Before assay procedure bring all reagents to room temperature 20-25°C .			
Assay Procedure in Double Determination:			
Pipette	Reagents		Position
100 µL	Dilution Buffer VP (Blank)		A1/A2
100 µL	Standard A (0.15 ng/mL)		B1/B2
100 µL	Standard B (0.45 ng/mL)		C1/C2
100 µL	Standard C (0.90 ng/mL)		D1/D2
100 µL	Standard D (1.8 ng/mL)		E1/E2
100 µL	Standard E (3.6 ng/mL)		F1/F2
100 µL	Standard F (6.0 ng/mL)		G1/G2
100 µL	Standard G (9.0 ng/mL)		H1/H2
100 µL	Control Serum KS1 (1:5 diluted)		A3/A4
100 µL	Control Serum KS2 (1:5 diluted)		B3/A4
100 µL	Sample (1:5 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
Sample Incubation: 1 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Antibody Conjugate AK		In each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Enzyme Conjugate EK		In each well
Mit Klebefolie die Vertiefungen dicht abdecken.			
Incubation: 0.5 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WP/ well.		In each well
100 µL	Substrate Solution S		In each well
Substrat S Incubation: 0.5 h in the Dark at 20-25°C			
100 µL	Stopping Solution SL		In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			