

16 ASSAY PROCEDURE

Preparation of reagents		
A-E	Standards (each 750 µL)	Bring the ready for use reagents to room temperature: 20°-25°C
KS1	Control Serum 1 (750 µL)	
KS2	Control Serum 2 (750 µL)	
AK	Antibody Conjugate (12 mL)	
EK	Enzyme Conjugate (12 mL)	
WP	Washing Buffer	Dilute 1:20 with Aqua dest.
Dilution of samples is generally not necessary; just use 20 µL per single determination.		
Before assay procedure bring all reagents to room temperature 20-25°C .		
Assay Procedure in Double Determination:		
Pipette	Reagents	Position
100 µL	Dilution Buffer VP	Pipette in all required number of wells
20 µL	Dilution Buffer VP (Blank)	A1/A2
20 µL	Standard A (1 ng/mL)	B1/B2
20 µL	Standard B (5 ng/mL)	C1/C2
20 µL	Standard C (10 ng/mL)	D1/D2
20 µL	Standard D (15 ng/mL)	E1/E2
20 µL	Standard E (25 ng/mL)	F1/F2
20 µL	Control Serum KS1	G1/G2
20 µL	Control Serum KS2	H1/H2
20 µL	Sample	Pipette sample in the rest of the wells according the requirements (Pipetting Control = colour turns red!)
Cover the wells with the sealing tape.		
Sample Incubation: 1 h at 20-25°C, 350 rpm		
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WP/ well	In each well
100 µL	Antibody Conjugate AK	In each well
Cover the wells with the sealing tape.		
Incubation: 30 Minutes at 20-25°C, 350 rpm		
100 µL	Enzyme Conjugate EK, without washing the wells (!) – add to the previously pipetted Antibody Conjugate AK -solution thereto mix shortly through cautious tapping on the side of the MTP . Attention: high filled volume of the wells!	In each well
Cover the wells with the sealing tape.		
Incubation: 30 Minutes at 20-25°C, without shaking		
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WP/ well	In each well
100 µL	Substrate Solution S	In each well
Incubation: 15 Minutes in the Dark at 20-25°C		
100 µL	Stopping Solution SL	In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		