anti- Pseudomonas aeruginosa IgG ELISA

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
Enzyme Immunoassay for qualitative and quantitative Determination of
human antibodies against **Pseudomonas aeruginosa**
English

Europäische Union / European Union*:
für In-vitro-Diagnostik / for in-vitro diagnostics
Nur zum Gebrauch durch Fachpersonal! / Only for professional use!
Rest of the world and USA***: for research use only!

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according to DIN EN 980 and EDMA recommendations Standard News 6 2001

Expiry date/ Verfallsdatum/ Date de péremption/ Data di scadenza/ Fecha de caducidad/ Data de validade/ Uiterste gebruiksdatum/ Udløbsdato/ Bäst före-datum/ Termin ważności/ Lejárati idő/ Čas expirácie/ Doba expirace/ Срок на годност/ Aegumiskuupäev/ Ημερομηνία λήξης/ Data de expirare/ Rok uporabe/ Viimeinen käyttöpäivä

Consider instructions for use/ Bitte Gebrauchsanweisung beachten/ Consultez la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consulte las instrucciones de uso/ Respeitar as instruções de utilização./ A.u.b de gebruiksaanwijzing volgen/ Se brugsanvisningen/ Läs anvisningarna före användning/ Proszę przeczytać instrukcję obsługi/ Vegye figyelembe a használati utasításban foglaltakat/ Postupujte podľa pokynov na použitie/ Dodržujte návod k použití/ Моля, спазвайте инструкцията за употреба/ Palun järgige kasutusjuhendit./ Λάβετε υπόψη σας τις οδηγίες χρήσης/ Vă rugăm să respectați instrucțiunile de utilizare/ Upoštečajte navodila za uporabo! / Lue käyttöohje huolellisesti!

In vitro Diagnostic Medical Device (for in Vitro Diagnostic Use)/ in vitro Diagnostikum (zur In-vitro-Diagnostik)/ Dispositif médical de diagnostic in vitro (Pour usage diagnostique in vitro)/ Dispositivo médico per diagnostica in vitro (per uso diagnóstico in vitro)/ Dispositivo médico de diagnóstico in vitro (para uso diagnóstico in vitro)/ Dispositivo Médico para Diagnóstico in vitro (Para Utilização de Diagnóstico in vitro)/ Medisch hulpmiddel voor diagnostiek in vitro (voor diagnostisch gebruik in vitro)/ Medicinsk udstyr til in vitro-diagnostik (Udelukkende til in vitro diagnostisk anvendelse)/ Medicinteknisk produkt avsedd för in vitro-diagnostik (för in vitro-diagnostiskt bruk)/ Wyrób medyczny do diagnostyki in vitro/ In vitro orvosdiagnosztikai termék (in vitro diagnosztikai használatához)/ In vitro diagnostický zdravotnický materiál (určený na diagnostiku „in vitro“)/ In vitro diagnostický zdravotnický materiál (určeno pro diagnostiku „in vitro“)/ Медицинско устройство за ин-витро диагностика (за ин-витро диагностика)/ in vitro diagnostikaseade (in vitro diagnostika tegemiseks)/ In vitro διαγνωστικό (για διάγνωση in vitro)/ Dispozitiv de diagnosticare in vitro (pentru diagnosticarea in-vitro)/ In vitro diagnostika (o in vitro diagnostiki)/ in vitro-diagnostiikkakäyttö

Lot-Batch Number/ Charge-Chargennummer/ Lot-Code du lot/ Lotto-Numero di lotto/ Lote-Código de lote/ Lote-Código do Lote/Lot-Partijnummer/ Lot-Batchcode/ Partisatskod/ Numer serii/ Tétel-sarzs szám/ Číslo šarže/ Číslo šarže/ Партиден номер/Partii – partii number/ Παρτίδα-αριθμός παρτίδας/ Lot-numår lot/ Številka serije/ erä

Manufactured by/ Hergestellert von/ Fabriqué par/ Prodotto da/Fabricado por/ Fabricado por/ Vervaardigd door/Fabrikation af /Tillverkad av/ Wyprodukowane przez/ Gyártotta / Vyrobené/ Vyrobeno v/ Производител/ Tootja/ Κατασκευάζεται από/ Produs de/ Proizvajalec/ Valmistaja

Catalogue Number/ Bestellnummer/ Numéro de référence/Numero di riferimento/ Número de referencia/ Número de Referência/ Referentienummer/ Referencenummer /Beställningsnummer/ Numer katalogowy/ Rendelési szám/Katalógovú číslo/ Objednací číslo/Каталожен номер/Теллімісnumber/ Αρ. παραγγελίας/Număr comandă/ Številka naročila/ viite tai tilausnumero

Store at between/ Lagerung bei zwischen/ Conserver à entre/ Conservare a tra/ Conservar a temp. Entre/ Armazena entre/ Bewaar bij tussen/ Opbevares mellem/ Förvaras vid/ Przechowywać w/ Tárolási tartomány/ Skladujte v rozsahu / Skladujte v rozmezí/ Температурно ограничение/ Säilittä temperaturid/ Φύλαξη σε θερμοκρασία/ Depozitare între/ Skladiščenje med/ Säilytys x-y Celsiusasteen lämpötilassa

Contains sufficient for x tests/ Inhalt ausreichend für x Tests/ Contient suffisant pour x tests/ Contenido suficiente per x test/ Contenido suficiente para x pruebas/ Conteúdo suficiente para x testes/ Bevat voldoende voor x bepalingen/ Indeholder tilstrækkeligt til x prøver/ Innehållet räcker till x analyser/ Zawartość na x testów/ Tartalma x teszt elvégzésére elegendő/ Obsahuje materiál pre x testov / Obsah dostačuje pro x testů/ Съдържание достатъчно за x тестове/ Sisust jätkub x katse jaoks/ Το περιεχόμενο επαρκεί για x δοκιμές/ Conținut suficient pentru x teste/ Vsebina zadostuje za x preizkusov/ Sisältö riittää x testille

Keep away from sunlight/ Nicht dem Sonnenlicht aussetzen/ Conserver à l'abri de la lumière/ Conservare al riparo della luce solare/ No exponer a la luz solar/ Proteger da luz solar/ Niet aan zonlicht blootstellen/ Må ikke udsættes for sollys/ Utsätt inte för solljus/ Nie wystawiać na słońce/ Napfénytől távol tartandó/ Nevystavovat' slnečnému svetlu/ Nevystavovat' slnečnému svetlu/ Да се предпазва от слънчева светлина/ Kaitsta otsese päikesekiirguse eest/ Κρατήστε το μακριά από την ηλιακή ακτινοβολία/ Τηνετή departe de lumina soarelui/ Ne izpostavljajte sončni svetlobi/ suojaa auringonvalolta

Incubation time/ Inkubationszeit/ Temps d'incubation/ Tempo d'incubazione/ Tiempo de incubación/ Tempo de incubação/ incubationstijd/ Inkubationstid/ inkubationstid/ Czas inkubacji/ Inkubációs idő/ Inkubačná lehota/ Inkubační doba/ Инкубационен период/ Inkubatsiooniaeg/ Χρόνος επώασης/ Timp de incubare/ Inkubacijska doba/ inkubaatioaika

incubate at / Inkubation bei/ Incuber à/ Incubare a/incubar a/Incubar a/ incubatietemperatuur/ Inkubation ved/ inkubation vid/ Inkubacja przy/ Inkubáció hőmérséklete/Inkubácia pri/ Inkubace při/Инкубира се при/ Inkubatsioon temperatuuril/ Επώαση στους/ Incubare la/ Inkubacija pri/ inkubaatiolämpötila

Shaking/ Schütteln/ Mélanger/ Agitare/ Agitar/ Agitação/ Schudden/ Ryster/ Skaka/ Wstrząsanie/ Rázás/ Pretrepat/ Protřepat/ Разклащане/ Raputada/ Ανακινήστε/ Vibrare/ Stresite/ Sekoita

Mikrotiterplate/ Mikrotiterplatte/ plaque de microtitrage/ Piastra di microtitolazione/ Placa de microtitulación/ Placa de Microtitulação/ microtiterplaat/ Mikrotiterplade/ mikrotiterplatta/ microtiterplaat/ Płytko microtiter/ Mikrotiter lap/ Mikrotitračná podložka/ Mikrotitrační podložka/ Микротитърна плака/ Mikrotiterplaat/ Τρυβλίο μικροτιτλοδότησης/ Microplacă/ Mikrotitrská plošča/ Mikrotitrauslevy

Sample/ Probe /Echantillon/ campione/ Muestra/ Amostra/ monster/ Prøve/ prov/ Próbka/ Minta/ Vzorka/ Vzorek/ Проба/ Proov/ Δείγμα/ Probă/ Vzorec/ Näyte

Conjugate Concentrate / Konjugatkonzentrat/ conjugué conc / Coniugato concentrato / Conjugado concentrado / Konjugado Concentrado / conjugaat geconcentreerd / konjugat koncentrat/ konjugat koncentrat/ konjugat koncentrat / Koniugat koncentrátum / Antitest és enzim páros/ konjugát/ конюгат Концентрат / konjugaat konsentraat /

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Package Insert

English

TECHNICAL FEATURES	15
INTRODUCTION	15
TEST PRINCIPLE	15
INTENDED USE	16
PERFORMANCE CHARACTERISTICS AND VALIDATION	16
SPECIMEN COLLECTION, PREPARATION AND STORAGE	17
REAGENTS PROVIDED	17
MATERIALS REQUIRED BUT NOT PROVIDED	18
TECHNICAL NOTES	18
STORAGE CONDITIONS	18
WARNINGS AND PRECAUTIONS	19
ASSAY PROCEDURE	20
CALCULATION OF RESULTS	20
Qualitative Calculation:	20
Semi-quantitative calculation	21
LITREATURE / LITERATUR	24
KURZANLEITUNG – Mediagnost P. aeruginosa EIA E15	25
SUMMARY – MEDIAGNOST P. aeruginosa IgG EIA E15	27

* please ask for package inserts in your national language

** please ask for special package insert

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PACKAGE INSERT ENGLISH

TECHNICAL FEATURES

- ◆ Detects Antibodies against three Antigens Alkaline Protease, Exotoxin A and Elastase
- ◆ high sensitivity and specificity
- ◆ low intra- and interassay variability

INTRODUCTION

Pseudomonas aeruginosa, a Gram-negative bacterium ubiquitously distributed in the moist environment, causes about 10% of all nosocomial infections. This opportunistic pathogen leads to acute and chronic types of infection within the various organs of susceptible patient groups. Chronic pulmonary infection in patients with cystic fibrosis (CF) is very frequent, may start early in infancy and determines the life expectancy of these patients. *P. aeruginosa* infection provokes a rapid production of antibodies to a large number of *P. aeruginosa* antigens in CF patients. **Mediagnost's** sensitive antibody detection system discriminates between infected and uninfected patient groups, detects a *P.aeruginosa* infection very early at the onset, when microbiological information is not yet readily available as seen with very young CF infants and reveals the chronicity of the infection when high specific antibody titres are determined. Thus, the assay may allow the commencement of anti-*Pseudomonas* chemotherapy early after onset of the infection and may also serve as a tool to control the efficiency of antibiotic therapy. Due to the use of three *P.aeruginosa* antigens, which are highly immunogenic and present in different parts from nearly all *P.aeruginosa* strains, together with the high sensitivity of the assay, false-negatives as well as false-positive results are almost impossible. Most patients immune systems react with at least one of the three antigens. The present EIA is a further development of a radioimmunoassay used successfully for more than 10 years in the diagnosis of *P.aeruginosa* infection in CF and other patient groups.

The semi-quantitative detection of human sera using the **mediagnost anti-*Pseudomonas aeruginosa* IgG EIA , E15**, is subdivided into the following categories:

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Depending on the *Pseudomonas aeruginosa* species and the patient's immune reaction, antibodies can be detected against a single, two or even all three antigens simultaneously. A patient is regarded as sero-positive when the serum is positive for one or more of the antigens. In acute and chronic *Pseudomonas aeruginosa* infected patients suffering from immunosuppression, assessment of an infection via the detection of antibodies may be inadequate.

TEST PRINCIPLE

The **Mediagnost anti-*Pseudomonas aeruginosa* IgG EIA, E15**, is a sandwich enzyme immunoassay. Serum or plasma samples are diluted and added to the wells of a microtiter plate, which have been previously coated with the *Pseudomonas aeruginosa* antigens alkaline protease, elastase or exotoxin A. Specific antibodies in the sample bind to the antigens present during an incubation of 2h at 37°C. After washing, the conjugate (anti-human IgG peroxidase-labelled immunoglobulin) is added and incubated again (for 2 h at 37°C). After a final washing step,

substrate is added and further incubated for 30 min at room temperature. The reaction is terminated on addition of stop solution accompanied by a change from blue to yellow. The absorbance of the coloured reaction product is measured on a microtiter plate reader. The colour intensity of the reaction corresponds to the concentration of antibodies in the sample.

INTENDED USE

The Mediagnost anti-*Pseudomonas aeruginosa* IgG EIA, E15, is an enzyme immunoassay for qualitative and semi-quantitative detection of IgG-Antibodies in human serum against the extracellular proteins: **Alkaline Protease (AP)**, **Elastase (Ela)** and **Exotoxin A (Exo)** of *Pseudomonas aeruginosa*.

PERFORMANCE CHARACTERISTICS AND VALIDATION

Table 1 : Linearity

Dilution	Sample		
	AP	Ela	Exo
1:750	1445		3623
1:1000	1429		3846
1:1500	1371		4123
1:3000	1162	10893	4230
1:6000		11135	4375
1:12000		10905	4590

Table 2 : Inter-Assay-Variation

	Number	Mean (Titer)			Standard deviation			Coefficient of Variation (%)		
		AP	ELA	Exo	AP	ELA	Exo	AP	ELA	Exo
Sample 1	74	1041	958	806	66	60	76	6	6	9
Sample 2	73	1516	1662	1167	77	135	74	5	8	6
Sample 3	72	3031	4129	4566	147	306	303	5	7	7

Table 3: Intra-Assay-Variation

	Number	Mean (Titer)			Standard deviation			Coefficient of Variation (%)		
		AP	ELA	Exo	AP	ELA	Exo	AP	ELA	Exo
Sample 1	16	3858	5521	5673	173	135	127	4.47	2.45	2.25
Sample 2	16	1453	1530	1381	88	44	142	6.08	2.86	10.29

Clinical validation was done by the *Hygiene Institute of the University of Tübingen* Prof. Dr. Döring. In comparison with the microbiological detection of *Pseudomonas aeruginosa* in sputum and throat swab specimens the Mediagnost E15 Antibody Assay shows the same results in 104 samples (52 positive and 52 negative in microbiological culture). Recently, sensitivity and specificity were investigated again and published in *Thorax (Kappler, M et al. Thorax published online 31 Jan 2006; doi:10.1136/thx.2005.049536: Diagnostic and Prognostic Value of Serum Antibodies against Pseudomonas aeruginosa in Cystic Fibrosis)*. The results clearly demonstrate the high quality of the Mediagnost Assay.

SPECIMEN COLLECTION, PREPARATION AND STORAGE

Serum samples are suitable. A special external sample preparation prior to assay is not required. Slight hemolysis of the samples doesn't disturb the determination.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination, store the undiluted samples frozen -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote).

In most determinations (e.g. Serum or Plasma samples and no extreme values expected) the dilution of **1:1000 with Dilution Buffer VP is suitable**.

If the extinction of a sample exceeds the extinction of the positive control **PK** measurement should be repeated with the more diluted sample. For valid titer results samples should be measured in three dilutions (1:1000, 10000 and 100000).

Suggested Dilution Protocol:

Pipette 90 µl **Dilution Buffer VP** in PE-/PP-tubes (application of a multi-stepper is recommended in larger series), add **10 µl Serum- or Plasma** (dilution 1:10) and mix each tube.

In a second dilution step the pre-diluted samples should be diluted serial as follows:

10µl (1:10) Dilution in 990µl VP	1:1000
100µl (1:1000) Dilution in 900µl VP	1:10000
100µl (1:10000) Dilution in 900µl VP	1:100000

After mixing use **100 µl** of this solution within 1 hour **per determination** in the assay.

REAGENTS PROVIDED

1)	<input type="checkbox"/> MTP AP	Microtiter plate , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with alkaline Protease and labelled red , packed in a laminate bag.
2)	<input type="checkbox"/> MTP ELA	Microtiter plate , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with Elastase and labelled blue , packed in a laminate bag.
3)	<input type="checkbox"/> MTP EXO A	Microtiter plate , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with Exotoxin A and labelled green , packed in a laminate bag.
4)	<input type="checkbox"/> BUF VP	Dilution Buffer VP , 3 × 50 ml, ready-to-use, please use for the dilution of the conjugate concentrate KK as well as for the dilution of samples
5)	<input type="checkbox"/> Control POS PK1	Positive Control PK1 , 1.5 ml, ready-to-use, calibrated control serum for alkaline Protease AP red labelled shows a titer of 1:2500
6)	<input type="checkbox"/> Control POS PK2	Positive Control PK2 , 1.5 ml, ready-to-use, calibrated control serum for Elastase Ela blue labelled shows a titer of 1:2500
7)	<input type="checkbox"/> Control POS PK3	Positive Control PK3 , 1.5 ml, ready-to-use, calibrated control serum for Exotoxin A Exo green labelled shows a titer of 1:2500
8)	<input type="checkbox"/> Control NEG NK	Negative Control NK , 2 × 1.5 ml, ready-to-use, contains human serum, not reactive for Pseudomonas aeruginosa antigens
9)	<input type="checkbox"/> Control	Control Serum KS , 2 × 1.5 ml, ready-to-use, contains human serum and has to be determined as borderline in the assay
10)	<input type="checkbox"/> CONJ	Conjugate Concentrate KK , 2 × 250 µl, 100fold Concentrate, contains POD (horseradish-Peroxidase)-labelled anti-human IgG. After dilution the Cojugate is only limited stable, dilute only according to requirements.
11)	<input type="checkbox"/> WASHBUF 20x	Washing Buffer (WP) , 120 ml, 20 X concentrated solution. Washing Buffer (WP) has to be diluted 1:20 with distilled or demineralised water before use. Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
12)	<input type="checkbox"/> SUBST	Substrate (S) , 33 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H ₂ O ₂ Tetramethylbencidine.
13)	<input type="checkbox"/> H ₂ SO ₄	Stopping Solution (SL) , 33 ml, ready for use, 0.2 M sulphuric acid, Caution acid!
14)		Sealing tape for covering of the microtiter plate, 6 x, adhesive.

MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes and multichannel pipettes with disposable plastic tips
Distilled or deionized water for dilution of the Washing Buffer (WP)
Vortex-mixer
Microtiter plate washer (recommended)
Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm
Polyethylen PE/Polypropylen PP tubes for dilution of samples

TECHNICAL NOTES

In conducting the assay, follow strictly the test protocol. Room temperature incubation means: Incubation at 20 - 25°C.

Reagents with different lot numbers should not be mixed. The microtiter plate and all reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

The shelf life of the components after opening is not affected, if used appropriately. Store the unused stripes of the microtiter plate together with the desiccant airtight at 2-8°C.

Use the **Dilution Buffer VP** for the dilution of **Conjugate Concentrate KK**.

The 1:20 diluted Washing Buffer **WP** is only 4 weeks stable. Please dilute only according to daily requirements. This applies to the 1:100 diluted **Conjugate Concentrate KK** too.

Before use, all kit components should be brought to room temperature. **Precipitates, possible in buffers, should be dissolved before use through mixing and warming.**

Incubation at room temperature means: Incubation at 20-25°C

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 μ l at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtiter plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

The **Substrate Solution S**, stabilised H_2O_2 -Tetramethylbencidine, is photosensitive – store and incubate in the dark.

When performing the assay, Controls **KS**, **PK**, **NK** and the samples should be pipetted as fast as possible (e.g., 15 minutes).

To avoid distortions due to differences in incubation times the diluted Conjugate Concentrate **KK** as well as the succeeding

Substrate Solution S should be added to the plate in the same order and in the same time interval as the samples. **Stop**

Solution SL should be added to the plate in the same order as the Substrate Solution **S**

STORAGE CONDITIONS

The microtiter plate wells and all undiluted reagents are stable until the expiry date if stored in the dark at 2-8°C.

Store the unused microtiter wells together with the desiccant at 2° to 8°C.

The shelf life of the components after opening is not affected, if used appropriately.

The Substrate Solution (S), stabilised H_2O_2 -Tetramethylbencidine, is photosensitive – store and incubate in the dark.

WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use only. For professional use only.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to **room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature WILL affect the absorbance** readings of the assay. However, Values for the patient samples will not be affected.

Do not mix reagents of different lots. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of human and/or animal origin. Source human serum for the Control Serum provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

Stop solution contains 0.2 M Sulfuric Acid (H₂SO₄)

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.

Following components contain **0.01% 2-Methyl-4-isothiazolin-3-one Solution** as preservative: **KK, VP**

R34 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S36/37 Wear suitable protective clothing and gloves

S45 In case of accident or if you feel unwell seek medical advice

Following components contain **0.01% (w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-isothiazol-3-one** as preservative: **KK, VP, WP**

R36/38 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed

R36/37/38 Irritating to eyes, respiratory system and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves

General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

ASSAY PROCEDURE

NOTES: All determinations (Controls and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Positive Control, Negative Control, Control Serum and the samples should be pipette as fast as possible (e.g., <15 minutes). To avoid distortions due to differences in incubation times, the diluted **conjugate concentrate KK** as well as the following **Substrate Solution S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution SL** should be added to the plate in the same order as the Substrate Solution.

- 1) pipette in positions* A1/2 **100 µl Negative Control NK each**
- 2) pipette in positions* B1/2 **100 µl Positive Control PK each** (PK1,PK2 OR PK3, respectively)
- 3) pipette in positions*C1/2 **100 µl Control Serum KS each**
* of Microtiter plates PL1/PL2/P3
- 4) Pipette **100 µl each** of the **diluted sample** (generally 1:1000 diluted in Dilution Buffer **VP**) in the rest of the wells, according to requirements
- 5) Cover the wells with the sealing tape and incubate the plate for **2 hours at 37°C**
- 6) After incubation aspirate the contents of the wells and wash the wells 3 times with **300 µl Washing Buffer WP**. Empty the wells thoroughly.
- 7) Following the last washing step, pipette **100 µl** of the **1:100 diluted Conjugate Concentrate KK** in each well.
- 8) Cover the wells with the sealing tape and incubate the plates for **2 hours at 37°C**.
- 9) After incubation wash the wells 3 times with **Washing Buffer WP** as described in step 6)
- 10) Pipette **100 µl of the TMB-Substrate solution S** in each well.
- 11) Incubate the plates for **30 Minutes in the dark at room temperature**.
- 12) After incubation pipette **100 µl Stop Solution SL** in each well.

Measure the absorbance **within 30 minutes at 450 nm (Reference filter ≥590 nm)**

CALCULATION OF RESULTS

Calculations should be performed for each antigen (each plate) separately! Calculate the average of all multiple values.

For the evaluation of the assay it is preconditioned that the absorbance values of the Negative Controls NK should be below 0.25. The difference between the extinctions of Negative Controls NK and the respective Positive Controls PK must be at least 0.6.

Qualitative Calculation:

The negative control average is subtracted from the controls and samples to obtain absolute values (Blank).

Cut-off calculation:

The cut-off is 20% of the positive control average.

This value corresponds to a titre of 1:500 of a 1:1.000 diluted serum. 1:1.000 diluted sera with extinction values less than the cut off are regarded as negative.

Depending on the *Pseudomonas aeruginosa* species and the patient's immune reaction, antibodies can be detected against a single, two or even all three antigens simultaneously. A patient is regarded as sero-positive when the serum is positive for one or more of the antigens. In acute and chronic *Pseudomonas aeruginosa* infected patients suffering from immunosuppression, assessment of an infection via the detection of antibodies may be inadequate.

Negative Control NK		Extinction
1.Value		0.041
2.Value		0.056
Mean		0.049

Positive Control PK		Extinction
1.Value		1.120
2.Value		1.136
Mean		1.128

Calculation

PK – NK	: 1.128 – 0.049	= 1.079
Cut-off (20% of PK-NK)	: 0.2 × 1.079	= 0.216
Borderline (50% of PK-NK)	: 0.5 × 1.079	= 0.540

All samples with an extinction < 0.216 are determined as negative for anti-*P. aeruginosa* IgG. Samples with an extinction of > 0.216 and < 0.540 are judged as borderline and samples showing an extinction of > 0.540 are determined as positive for anti-*P. aeruginosa* IgG.

Semi-quantitative calculation

The measured extinction values of the individual serum dilutions are plotted onto a graph or calculated using an appropriate computer programme.

The extinction values of the negative (NK) and positive (PK) controls are plotted on the y-axis of a linear plot. A so-called titre factor is in turn plotted on the x-axis, for NK = 0 and PK = 2.5. For the quantification of the sera values, a straight line is drawn through the NK and PK values and extended to a value of 3.5.

The titre of the individual serum is determined by reading the **titre factor** of the measured extinction value through the NK-PK axis, which is multiplied by the serum dilution factor.

Titre factors lower than 0.25 and higher than 3.5 (x-axis) are not taken into consideration in these calculations, except that the titre factors from a 1:1000 dilution below 0,5 are regarded as negative and titre factors greater than 3,5 are always positive and should be diluted further and retested for quantitative determination.

Categories:


The semi-quantitative detection of human sera using the **mediagnost anti-*Pseudomonas aeruginosa* IgG EIA , E15**, is subdivided into the following categories:

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LITREATURE / LITERATUR

1. West S, Zeng L, Lee L, Kosorok M, Laxova A, Rock M, Splaingard M, Farrell : Respiratory infections with *Pseudomonas aeruginosa* in children with cystic fibrosis, early detection by serology and assesment of risk factors. JAMA, 2002, Vol 287, No.: 22: 2958-2967
2. Burns J, Gibson I, McNamara, Yim D, Emerson J, Rosenfeld M, Hiatt P, McCoy K, Castille R, Smith A, Ramsey B: Longitudinal Assesment of *Pseudomonas aeruginosa* in young children with cystic fibrosis. J Infect Dis 2001 ; 183 :444-52
3. Döring G, Obernesser H-J, Botzenhart K: Extracellular toxins of *Pseudomonas aeruginosa*. III Radioimmunoassay for detection of alkaline protease. Zentralbl Bakteriol Mikrob. Hyg (A) 1982; 252:239-47
4. Döring G., Obernesser H-J., Botzenhart K., Flehmig B., Høiby N., Hofmann A.: Proteases of *Pseudomonas aeruginosa* in cystic fibrosis. J Infect Dis 1983; 147:744-50
5. Döring G.; Høiby N: Longitudinal study of immune response to *Pseudomonas aeruginosa* antigens in cystic fibrosis. Infect Immun 1983; 42:197-201
6. Obernesser H-J, Döring G: Extracellular toxins of *Pseudomonas aeruginosa*. IV. Radioimmunoassay for detection of elastase. Zentralbl Bakteriol Mikrobiol Hyg (A) 1982; 252:248-56
7. Brett M M, Ghoneim A T M, Littlewood J M: Serum antibodies to *Pseudomonas aeruginosa* in cystic fibrosis. Archives of Diseases in Childhood, 1986, 61, 1114-20
8. Brett M M, Ghoneim A T M, Littlewood J M; Losowsky M S: Development of enzyme linked immunosorbent assay (ELISA) to detect antibodies to *Pseudomonas aeruginosa* cell surface antigens in sera of patients with cystic fibrosis. J Clin Pathol 1986; 39: 1124-29
9. Brett M M, Ghoneim A T M, Littlewood J M: Serum IgG antibodies in patient with cystic fibrosis with early *Pseudomonas aeruginosa* infection: Archives of Disease in Childhood, 1987, 62, 357-361
10. Brett M M, Ghoneim A T M, Littlewood J M: Prediction and Diagnosis of Early *Pseudomonas aeruginosa* Infection in Cystic Fibrosis: a Follow-up Study: J of Clin Microb, Aug 1988 p 1565-70
11. Brett M M, Simmonds E J, Ghoneim A T M, Littlewood J M: The value of serum IgG titres against *Pseudomonas aeruginosa* in the management of early peudomonal infection in cystic fibrosis. Archives of Diseases in Childhood 1992; 67: 1086-1088
12. Danielsen L, Westh H, Balselv E, Rosdahl VT, Döring G: *Pseudomonas aeruginosa* exotoxin A antibodies in rapidly deteriorating chronic leg ulcers. The Lancet Jan 1996 Vol. 347 No.8996 p 265
13. Dasgupta MK; Lam J, Döring G, Harley FL, Zuberbuhler P, Lam K, Reichert A, Costerton JW, Dossetor JB: Prognostic implications of circulating immune complex and *Pseudomonas aeruginosa*-specific antibodies in cystic fibrosis. J Clin Lab. Immunol. 1987, 23, 25-30
14. Döring G, Dalhoff A, Vogerl O, Brunner H, Dröge U, Botzenhard K. In vivo activity of Proteases of *Pseudomonas aeruginosa* in a rat model. The J of Infectious Diseases, Apri 1994, 4/Vol 149 p 532-537
15. Döring G, Buhl V, Høiby N, Schiøth P O, Botzenhart K: Detection of *Pseudomonas aeruginosa* in immune complexes isolated from sputum of cystic fibrosis patients. Acta path microb. immunol. scand Sect. C, 92: 307-312, 1984
16. Döring G, Goldstein W, Röll A, Schiøtz P O, Høiby N, Botzenhard K: Role of *Pseudomonas aeruginosa* Exoenzymes in Lung Infections of Patients with Cystic Fibrosis. Infection and Immunity, Sept 1985, p. 557-562

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SUMMARY – MEDIAGNOST P. aeruginosa IgG EIA E15

Dilution of Reagents and Samples		
Conjugate concentrate KK	in Dilution Buffer VP	1:100
Washing buffer WP	in Aqua dest. (e.g. add 100 ml WP in to a graduated flask and fill with A.dest to 2000 ml)	1:20

Dilute samples **1:1000** with dilution buffer (qualitative test). For quantitative antibody assays it is recommended to perform dilutions of **1:1.000**, **1:10.000** and **1:100.000**.

Proposal of Assay Procedure for double determinations

Pipette	Reagents	Plate 1 AP	Plate 2 Ela	Plate 3 Exo
2 x 100 µl	Negative Control NK	A1/A2	A1/A2	A1/A2
2 x 100 µl	Positive Control PK1	B1/B2		
2 x 100 µl	Positive Control PK2		B1/B2	
2 x 100 µl	Positive Control PK3			B1/B2
2 x 100 µl	Control Serum KS	C1/C2	C1/C2	C1/C2
2 x 100 µl	Sample dilution	Pipette in the rest of the wells according to requirements.		

Seal the wells with the sealing tape

Incubation: 2 h at 37°C

3 x 300 µl	Aspirate the contents of the wells and wash 3x with 300 µl each Wash Buffer WP/ well	in each well
100 µl	1:100 diluted Conjugate concentrate KK	in each well

Seal the wells with the sealing tape

Incubation: 2 h at 37°C

3 x 300 µl	Aspirate the contents of the wells and wash 3x with 300 µl each Wash Buffer WP/ well	in each well
100 µl	Substrate Solution S	in each well

Incubation: 30 min in the dark at RT

100 µl	Stopping Solution SL	in each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		

International Test description

CONJ	KK	1:100	DILU	BUF	VP	
WASHBUF	20x	WP				1:20 DILU A. dest.

SPE	1:1000	DILU	BUF	VP	↔	
°C	20-25 °C					

		MTP AP	MTP ELA	MTP Exo
100 µl	Control NEG NK	A1/2	A1/2	A1/2
100 µl	Control POS PK1	B1/2		
100 µl	Control POS PK2		B1/2	
100 µl	Control POS PK3			B1/2
100µl	CONTROL KS	C1/2	C1/2	C1/2
100 µl	SPE 1:1000 DILU BUF VP			
TAPE				

 2 h  37 °C

3x 300 µl	3x WASHBUF WP
100 µl	CONJ KK 1:100 DILU VP
TAPE	

 2 h  37 °C

3x 300 µl	3x WASHBUF WP
100 µl	SUBST TMB s

 0.5 h  20-25 °C 

100 µl	H₂SO₄ SL
MEASURE	