IGF-II - ELISA

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Human Insulin-like Growth Factor-II

(IGFBP-blockiert)

Deutsch

Enzyme Immunoassay for Quantitative Determination of

Human Insulin-like Growth Factor-II

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
Rest of the world and USA**: for research use only!

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according to DIN EN 980 and EDMA recommendations Standard News 6 2001

Expiry date/ Verfallsdatum/ Date de péremption/ Data di scadenza/ Fecha de caducidad/ Data de validade/ Uiterste gebruiksdatum/ Udløbsdato/ Bäst före-datum/ Termin ważności/ Lejárati idő/ Čas expirácie/ Doba expirace/ Срок на годност/ Aegumiskuurpäev/ Ημερομηνία λήξης/ Data de expirare/ Rok uporabe/ Viimeinen käyttöpäivä

Consider instructions for use/ Bitte Gebrauchsanweisung beachten/ Consultez la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consulte las instrucciones de uso/ Respeitar as instruções de utilização./ A.u.b de gebruiksaanwijzing volgen/ Se brugsanvisningen/ Läs anvisningarna före användning/ Proszę przeczytać instrukcję obsługi/ Vegye figyelembe a használati utasításban foglaltakat/ Postupujte podľa pokynov na použitie/ Dodržujte návod k použití/ Моля, спазвайте инструкцията за употреба/ Palun järgige kasutusjuhendit./ Λάβετε υπόψη σας τις οδηγίες χρήσης/ Vá rugám sã respectați instrucțiunile de utilizare/ Upoštevacite navodila za uporabo!/ Lue käyttöohje huolellisesti!

Lot-Batch Number/ Charge-Chargennummer/ Lot-Code du lot/ Lotto-Numero di lotto/ Lote-Código de lote/ Lote-Código do Lote/Lot-Partijnummer/ Lot-Batchkode/ Partisatskod/ Numer serii/ Tétel-sarzs szám/ Číslo šarže/ Číslo šarže/ Партиден номер/Partii – partii number/ Παρτίδα-αριθμός παρτίδας/ Lot-număr lot/ Številka serije/ erä

Manufactured by/ Hergestellt von/ Fabriqué par/ Prodotto da/Fabricado por/ Fabricado por/ Vervaardigd doo/Fabrikation af /Tillverkad av/ Wyprodukowane przez/ Gyártotta / Vyrobené/ Vyrobeno v/ Производител/ Tootja/ Κατασκευάζεται από/ Produs de/ Proizvajalec/ Valmistaja

Catalogue Number/ Bestellnummer/ Numéro de référence/Numero di riferimento/ Número de referencia/ Número de Referència/ Referentienummer/ Referencenummer /Beställningsnummer/ Numer katalogowy/ Rendelési szám/Katalógové číslo/ Objednací číslo/Каталожен номер/Tellimisnumber/ Αρ. παραγγελίας/Număr comandă/ Številka naročila/ viite tai tilausnumero

Store at between/ Lagerung bei zwischen/ Conserver à entre/ Conservare a tra/ Conservar a temp. Entre/ Armazena entre/ Bewaar bij tussen/ Opbevares mellem/ Förvaras vid/ Przechowywać w/ Tárolási tartomány/ Skladujte v rozsahu / Skladujte v rozmezi/ Температурно ограничение/ Säilätada temperatuuridel/ Φύλαξη σε θερμοκρασία/ Depozitare între/ Skladiščenje med/ Säilytys x-y Celsiusasteen lämpötilassa

Contains sufficient for x tests/ Inhalt ausreichend für x Tests/ Contient suffisant pour x tests/ Contenuto sufficiente per x test/ Contenido suficiente para x pruebas/ Conteúdo suficiente para x testes/ Bevat voldoende voor x bepalingen/ Indeholder tilstrækkeligt til x prøver/ Innehållet räcker till x analyser/ Zawartość na x testów/ Tartalma x teszt elvégzésére elegendő/ Obsahuje materiál pre x testov / Obsah dostahuje pro x testů/ Съдържание достатъчно за x тестове/ Sisust jätkub x katse jaoks/ Το περιεχόμενο επαρκεί για x δοκιμές/ Conținut suficient pentru x teste/ Vsebina zadostuje za x preizkusov/ Sisältö riittää x testille

Keep away from sunlight/ Nicht dem Sonnenlicht aussetzen/ Conserver à l'abri de la lumière/ Conservare al riparo della luce solare/ No exponer a la luz solar/ Proteger da luz solar/ Niet aan zonlicht blootstellen/ Må ikke udsættes for sollys/ Utsätt inte för solljus/ Nie wystawiać na słońce/ Napfénytől távol tartandó/ Nevystavovat slnečnému svetlu/ Nevystavovat slnečnému svetlu/ Да се предпазва от слънчева светлина/ Kaitsta otsese päikesekiirguse eest/ Κρατήστε το μακριά από την ηλιακή ακτινοβολία/ Țineți departe de lumina soarelui/ Ne izpostavljajte sončni svetlobi/ suojaa auringonvalolta

Incubation time/ Inkubationszeit/ Temps d'incubation/ Tempo d'incubazione/ Tiempo de incubación/ Tempo de incubação/ incubatietijd/Inkubationstid/ inkubationstid/ Czas inkubacji/ Inkubációs idő/ Inkubačná lehota/ Inkubační doba/ Инкубационен период/ Inkubatsiooniaeg/ Χρόνος επώασης/ Timp de incubare/ Inkubacijska doba/ inkubaatioaika

incubate at / Inkubation bei/ Incuber à/ Incubare a/incubar a/Incubar a/ incubatietemperatuur/ Inkubation ved/ inkubation vid/ Inkubacja przy/ Inkubáció hőmérséklete/Inkubácia pri/ Inkubace při/Инкубира се при/ Inkubatsioon temperatuuril/ Επώαση στους/ Incubare la/ Inkubacija pri/ inkubaatiolämpötila

Shaking/ Schütteln/ Mélanger/ Agitare/ Agitar/ Agitação/ Schudden/ Ryster/ Skaka/ Wstrząsanie/ Rázás/ Pretrepat/ Profřepat/ Разклащане/ Raputada/ Ανακινήστε/ Vibrare/ Stresite/ Sekoita

Miktrotiterplate/ Mikrotiterplatte/ plaque de microtitrage/ Piastra di microtitolazione/ Placa de microtitulación/ Placa de Microtitulação/ microtiterplaat/ Mikrotiterplatte/ mikrotiterplatta/ microtiterplaat/ Płytko microtiter/ Mikrotiter lap/ Mikrotitračná podložka/ Mikrotitrační podložka/ Микротитърна плака/ Mikrotiterplaat/ Τρυβλίο μικροπιλοδότησης/ Microplacă/ Mikrotitrska plošča/ Mikrotitruslevy

Reconstitute in/ Rekonstituieren in/ Reconstituer dans/ Ricostituire nel/ Reconstituier en/ Reconstituier em/ reconstituieren in/ Rekonstituier i/ rekonstituera/ Rekonstytuować w/ Helyreállítás/ Znovu připravit za/ Znovu pripraviti za/ Разтваряне в/ Moodustada uuesti / Ανασυστήστε σε/ Reconstituire în/ Predelava v/ rekonstituoi

Sample/ Probe /Echantillon/ campione/ Muestra/ Amostra/ monster/ Prøve/ prov/ Próbká/ Minta/ Vzorka/ Vzorek/ Проба/ Proov/ Δείγμα/ Probă/ Vzorec/ Näyte

Antibody and Enzyme Conjugate/ Antikörper und Enzym Konjugat/ anticorps conjugué et conjugué enzymatique/ Coniugato di anticorpo ed enzima/ Conjugado de anticuerpos y enzimas/ Conjugado Anticorpo-Enzima/ antilichaamen enzymconjugaat/ Antistoffer og enzym-konjugat/ antikropps- och enzymkonjugat (antikropp och enzym, konjugat)/ Koniugat antyciał i enzymów/ Antitest és enzim páros/ Protílátkový a enzymatický konjugát/ Protílátkový a enzymatický konjugát/ Антигяло и ензим конюгат/ Antikehad ja ensüümi konjugaat/ Σύμπλοκο αντισώματος-ενζύμου/ Compuși din anticorpi și enzime/ Antitelesa in konjugat encima/ Vasta-aine ja entsymi konjugaatti

Dilute in Buffer X/ Verdünnen in buffer X/ Fortyndes i buffer X/ spädi i buffert X/ Rozcieńczanie w buforze X/ Hígítás X pufferben/ Riedit' v pufrí X/ Ředit v pufru X/ Разреждане в буфер X/ Lahjendada puhvris X/ Αραιώστε σε ρυθμιστικό διάλυμα X/ Diluati în tamponul X/ Razredčiti v pufru X/ Iaimennetaan x puskuriin

Standard X /Standard X/ Etalon X/ Standard X/ Estándar X/ Standard X/ Standard X/ standaard X/ Standard X/ standard X/ Standard X/ Standard X/ Štandard X/ Standard X/ Стандарт X/ Standard X/ Πρότυπο X/ Standard X/ Standardni X/ Kalibraattori X

Standard X /Standard X/ Etalon X/ Standard X/ Estándar X/ Standard X/ Standard X/ standaard X/ Standard X/ standard X/ Standard X/ Standard X/ Štandard X/ Standard X/ Стандарт X/ Standard X/ Πρότυπο X/ Standard X/ Standardni X/ Kalibraattori X

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Package Insert

English

TECHNICAL FEATURES+APPLICATIONS	14
INTRODUCTION	14
INTENDED USE	14
PRINCIPLE	14
PERFORMANCE CHARACTERISTICS AND VALIDATION	15
SPECIMEN COLLECTION, PREPARATION, AND STORAGE	16
REAGENTS PROVIDED	17
MATERIALS REQUIRED BUT NOT PROVIDED	17
Technical Notes	18
WARNINGS AND PRECAUTIONS	20
ASSAY PROCEDURE	21
CALCULATION OF RESULTS	22
Establishing the Standard Curve	22
EXPECTED VALUES	23
LITREATURE / LITERATUR	24
KURZANLEITUNG – Mediagnost IGF-II-ELISA E30	25
SUMMARY – MEDIAGNOST IGF-II ELISA E30	26
<u>REF</u> E30 International Test Description	28

* please ask for package inserts in your national language

** please ask for special package insert

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(Blum W., Schweizer R.,: Insulin-like growth factors and their binding proteins; in Ranke MB (ed): Diagnostics of endocrine function in children and adolescents. Basel, Karger, 2003, pp 166-199 (22).

PACKAGE INSERT ENGLISH

TECHNICAL FEATURES+APPLICATIONS

- Measurement of IGF-II in **non-extracted** serum and plasma samples
- Calibrated against the **International Standard: WHO NIBSC 96/538**
- No interference by IGF-binding proteins through excess IGF-I
- Precise measurement of very low IGF-II levels: **high sensitivity** of **0.02 ng/ml**
- Inter- and Intra-Assay variance: maximal 7.2 and 6.6%

INTRODUCTION

The insulin-like growth factors (IGF)-I and -II play a pivotal role in the regulation of proliferation and differentiation of several tissue types (1-3). IGF-I also called Somatomedin C (4) has a molecular weight of 7.469 kDa (5). Its expression is mainly regulated by Growth Hormone and nutrition (6). But several hormones and peptide factors are known to influence IGF-II synthesis in different tissues. Bioavailability of the IGFs is regulated by specific binding proteins (IGFBP). Beside the high affinity Insulin-like Growth Factor Binding Proteins 1-6, IGFs are also bound by IGFBP-related Proteins (7, 8, 22). These binding proteins bind IGF-I and IGF-II with the same affinity or prefer IGF-II (9, 10). Direct measurement of IGFs in serum samples without pretreatment results in false values because of the extremely slow dissociation of the IGF/IGFBP complexes during the assay incubation only a part of the IGF-II in the specimen can bind to the antibodies and be detected.

Therefore, various techniques were applied to physically separate IGF-II from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2, 12, 13) These techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries.

This assay is easy, fast and results do not depend on the binding protein concentration of the sample. It is based on the high specificity of the employed antibodies for IGF-II. There is virtually no cross-reactivity with IGF-I. This allows the separation of IGF-II from the binding proteins by acidification and blocking of the free binding proteins with IGF-I. Thus, the endogenous IGF-II is free in solution.

INTENDED USE

Scientific investigations in the field of neonatal hypertrophy (IGF-II is a foetal growth factor) and malignancies (IGF-II is a monogenic growth factor). Age dependent reference values are shown in Table 4.

IGF-II seems to be of use in differential diagnostics of malignancies. Thus, it is possible to differentiate by IGF-II between adrenocortical carcinomas and adenomas (24). Further tumor staging and differentiation between hyperplasia and carcinoma can be improved by IGF-II measurements in prostate tumors (25). The IGF-System seem to be of relevance in neurodegeneration as well, e.g. Alzheimer's and Parkinson's diseases (26).

PRINCIPLE

The Mediagnost ELISA for IGF-II E30 is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. The first antibody, immobilized on the microtiter

plate, and the added second biotinylated antibody are binding the IGF-II in the sample. The Streptavidin-Peroxidase Enzyme Conjugate subsequently binds to the complex. In the closing substrate reaction the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-II-level of the samples.

IGF-II-IGFBP complex is dissociated by dilution in an acidic buffer. IGFBPs are blocked by IGF-I excess, thus allowing the measurement of free IGF-II. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the low cross-reactivity of the IGF-II antibody with IGF-I, excess IGF-I does not disturb the interaction of the first antibody with IGF-II

PERFORMANCE CHARACTERISTICS AND VALIDATION

The IGF-II ELISA E30 is calibrated against the **International Standard: WHO NIBSC 96/538**

The standards of the ELISA E30 are **human IGF-II** in concentrations **0.45; 1.5; 3; 5.63 and 9 ng/ml**, respectively the assay range covers –at recommended normal sample dilution- the **range to 2400 ng/ml**. By varying the sample dilution this can be adapted to the special individual requirements.

Sensitivity

The **analytical sensitivity** of the ELISA E30 yields **0.02 ng/ml** (2 SD of zero standard in 20fold determination).

The **Inter-** and **Intra-Assay variation** coefficients are less than 7.2% and 6.6 % respectively. Exemplary determinations are represented in the Table 1 and 2.

Table 1: Inter-Assay-Variation (n=12)

	Mean value (ng/ml)	Standard deviation (ng/ml)	Variation Coefficient (%)
Sample 1	381.53	27.54	7.22
Sample 2	817.81	57.70	7.06
Sample 3	639.41	45.47	7.11

Table 2: Intra-Assay-Variation (n=16)

	Mean value (ng/ml)	Standard deviation (ng/ml)	Variation Coefficient (%)
Sample 1	666.38	20.45	3.07
Sample 2	875.22	57.89	6.61

Because of the high specificity of the used antibodies the endogenous and structurally very similar IGF-I existing in variable concentration in the samples does not affect the correctness of the IGF-II determination up to 1000 ng/ml (s. Table 3).

Table 3: Dependence of the determined IGF-II concentration on the IGF-I content of the samples

Addition IGF-I (ng/ml)	Serum Sample (ng/ml)	NIBSC 96/538 in Samplebuffer (ng/ml)
0	603	3,49
100	630	3,39
250	605	3,44
500	614	3,59
1000	604	3,68

SPECIMEN COLLECTION, PREPARATION AND STORAGE

Serum samples as well as Heparin-, EDTA- and Citrat-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Furthermore suitable samples are: urine, saliva (low concentration, at least 1:10 dilution), cerebrospinal fluid (at least 1:10 dilution) and cell culture medium (inc. 5% FCS, at least 1:5 dilution). IGF-II concentration in other body fluids or cell culture supernatants can deviate strongly from the serum values.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although IGF-II levels were found to be unaffected by few cycles in our experiments.

The IGF-II ELISA E30 makes the correct determination possible of the samples starting from a dilution of 1:401 over a far concentration range.

The high sensitivity of the assays allows IGF-II determinations in small sample volumes, which is limited by pipetting accuracy rather than the amount of IGF-II. Generally a serum or plasma dilution of **1:401 is very well suited**. However, the respective suitable dilution should be examined first.

Suggestion for dilution protocol:

Please pipette 2000 μl **Sample Buffer PP** in PE/PP-Tubes (application of a multi-stepper is recommended in larger series); subsequently add 5 μl serum or plasma samples (dilution 1:401). After mixing **use 50 μl of this dilution per determination**.

Because the pipetting accuracy can rise by the use of 10 μl sample, a 2-step dilution is alternatively possible, for this place 1000 μl Sample Buffer PP in in PE/PP-Tubes, add 10 μl sample (samples are 1:101 diluted), mix from this solution 50 μl with 150 μl Sample Buffer (samples are 1:404 diluted). After mixing **use 50 μl of this dilution per determination**.

REAGENTS PROVIDED

1)	MTP	Microtiter plate , ready for use: Microtiter plate with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with anti-human IGF-II Antibody, packed in a laminate bag.
2)	CAL	Standards A-E , lyophilised, contain human IGF-II. Standard values are between 0.45 and 9 ng/ml (0,45; 1,5; 3; 5,63; 9 ng/ml) IGF-II and have to be reconstituted in 500 µl (each) in Sample Buffer PP . After using store the reconstituted standards in the original flasks as soon as possible at -20°C . When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. 50 µl per well are used in the assay
4)	BUF PP	Sample Buffer PP , 125 ml, ready for use, please use for reconstitution of Standards and Controls and for dilution of samples and Controls
5)	Control	Control Serum KS1, KS2, 250 µl , lyophilised, contain human serum and has to be reconstituted in 250 µl Sample Buffer PP . The reconstituted Control Sera must be stored in the original flask as soon as possible at -20°C after using. When using anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. The IGF-II target value concentrations and the respective ranges are given on the vial labels. The dilution of the Control Sera should be according to the dilution of the respective samples.
6)	Ab	Antibody Conjugate AK , 6 ml, ready for use, contains the biotinylated anti-IGF-II antibody. Use 50 µl for each well in the assay.
7)	CONJ	Enzyme Conjugate EK , 12 ml, ready for use, contains horseradisch-peroxidase conjugate to streptavidin, use 100 µl for each well in the assay.
8)	WASHBUF 20x	Washing Buffer (WP) , 50 ml, 20X concentrated solution. Washing Buffer (WP) has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
9)	SUBST	Substrate (S) , 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H_2O_2 Tetramethylbencidine.
10)	H_2SO_4	Stopping Solution (SL) , 12 ml, ready for use, 0.2 M sulphuric acid, Caution acid!
11)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes and multichannel pipettes with disposable plastic tips
 Distilled or deionized water for dilution of the Washing Buffer (WP)
 Vortex-mixer
 Microtiter plate shaker (350 rpm)
 Microtiter plate washer (recommended)
 Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm
 Polyethylen PE/Polypropylen PP tubes for dilution of samples

Technical Notes

The assay has to be conducted strictly according the test protocol herein.

Reagents with different lot numbers cannot be mixed. The microtiterplate and reagents are stable until the indicated expiry, if stored unopened and protected from sunlight at 2 – 8°C.

The shelf life of the components after opening is not affected, if used appropriately.

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming.

Incubation at room temperature means: 20-25°C

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must become adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µl at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtiter plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Standards and Controls

For the reconstitution of the lyophilised components (Standards A - E and Control Sera KS1 & KS2) the kit Sample Buffer PP has to be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The reconstituted standards and controls can be stored for 3 months at –20°C. Repeated freeze/thaw cycles have to be avoided. When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

Washing Buffer

The required volume of washing buffer is prepared by 1:20 dilution of the provided 20-fold concentrate with deionised water. The diluted Washing Buffer is stable for max. 4 weeks at 2-8°C.

Substrate Solution

The **Substrate Solution S**, stabilised H₂O₂-Tetramethylbencidine, is photosensitive – store and incubate in the dark.

Microtiter plate

Store the once unused microtiter strips and wells together with the desiccant in the tightly closed clip lock bag at 2-8°C use in the frame provided. The labelled expiry is not influenced in case of proper storage.

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WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use only. For professional use only.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood. The Mediagnost GmbH is not liable for any loss or harm caused by non-observance of the instructions, as far as no law withstands.

Temperature WILL affect the absorbance readings of the assay. However, values for the patient samples will not be affected.

Do not use expired reagents.

Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

Caution: This kit contains material of human and/or animal origin.

Human Serum

Contained in following components: **Control Serum KS1 and KS2.**

The sources of human sera were tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibody. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

Stop solution contains 0.2 M Sulfuric Acid (H₂SO₄)

R36/38 Irritating to eyes and skin
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1 After contact with skin, wash immediately with plenty of water
S36/37 Wear suitable protective clothing and gloves

2-Methyl-4-Isothiazolin-3-one

contained in following components: **AK, EK, VP, PP**

< 0.01% 2-Methyl-4-isothiazolin-3-one Solution

R34 Irritating to eyes and skin
R43 Sensibilisation through skin contact possible
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S36/37 Wear suitable protective clothing and gloves
S45 In case of accident or if you feel unwell seek medical advice

5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-Isothiazol-3-one

contained in following components: **AK, EK, VP, WP, PP**

< 0.01% (w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-Isothiazol-3-one Solution

R36/38 Irritating to eyes and skin
R43 Sensibilisation through skin contact possible
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1 After contact with skin, wash immediately with plenty of water

TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed
R36/37/38 Irritating to eyes, respiratory system and skin
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S28.1 After contact with skin, wash immediately with plenty of water
S36/37 Wear suitable protective clothing and gloves

General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

ASSAY PROCEDURE

NOTES: All determinations (Standards, Control Sera and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Standards, Control Sera and the samples should be pipetted as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times **Enzyme Conjugate EK**, the **Substrate Solution S** as well as the **Stop Solution SL** should be added to the plate in the same order and in the same time interval each, respectively.

- 1) Add **50 µl Antibody Conjugate AK** in **all** wells used
- 2) Pipette in positions A1/2 **50 µl Sample Buffer PP**
- 3) Pipette in positions B1/2 **50 µl of the Standard A (0.45 ng/ml)**,
pipette in positions C1/2 **50 µl of the Standard B (1.5 ng/ml)**,
pipette in positions D1/2 **50 µl of the Standard C (3 ng/ml)**,
pipette in positions E1/2 **50 µl of the Standard D (5.63 ng/ml)**,
pipette in positions F1/2 **50 µl of the Standard E (9 ng/ml)**.

To control the correct accomplishment of the assay **50 µl** of the 1:401 (or in respective dilution ratio of the samples) in Sample Buffer diluted **Control Sera KS1/KS2** can be pipetted in positions G1/2 and H1/2 .

Pipette **50 µl** each of the diluted samples (e.g. dilute 1:401 with Sample Buffer **PP**) In the rest of wells, according to your requirements.

- 4) Cover the wells with sealing tape and incubate the plate for **2 hours at room temperature** (shake at 350 rpm)
- 5) After incubation aspirate the contents of the wells and wash the wells **5 times 300 µl Washing Buffer WP** / well.
- 6) Following the last washing step pipette **100 µl** of the **Enzyme Conjugate EK** in each well.
- 7) Cover the wells with sealing tape and incubate the plate for **30 Minutes at room temperature** (shake at 350 rpm).
- 8) After incubation wash the wells **5 times** with Washing Buffer **WP** as described in step 5.
- 9) Pipette **100 µl** of the **Substrate Solution S** in each well.
- 10) Incubate the microtiter plate for **30 minutes in the dark at room temperature**.
- 11) Stop the reaction by adding **100 µl Stopping Solution SL** to all wells.
- 12) Measure the absorbance within **30 minutes at 450 nm (Reference filter ≥ 590 nm)**

CALCULATION OF RESULTS

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of standard E should be above 1.00.

Samples, which yield higher absorbance values than **Standard E**, are beyond the standard curve, for reliable determinations such samples should be retested at a higher dilution.

Establishing the Standard Curve

The standards provided contain the following concentration of recombinant hIGF-II:

Standard	A	B	C	D	E
ng/ml	0.45	1.5	3	5.63	9

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The **IGF-II concentration in ng/ml** of the samples can be calculated **by multiplication** with the **respective dilution factor**.

EXPECTED VALUES

Table 4: Serum levels of IGF-II in ng/ml in healthy persons at various ages*

Age group	Percentile		
	5th	50th	95th
Newborns	158	284	516
1-4 weeks	350	486	673
1-6 months	348	551	871
6-12 months	388	582	876
1-3 years	384	596	926
3-5 years	397	617	920
5-7 years	419	638	973
7-9 years	433	656	997
9-11 years	442	662	994
11-13 years	448	671	1006
13-15 years	455	679	1014
15-17 years	452	686	1042
20-30 years	436	679	1058
30-40 years	442	680	1049
40-50 years	407	650	1039
50-60 years	396	644	1049
60-70 years	373	611	1000

* Measurement was performed after acid-ethanol extraction, and values were corrected for recovery (correction factor 1.2). Blum W., Schweizer R.,: Insulin-like growth factors and their binding proteins; in Ranke MB (ed): Diagnostics of endocrine function in children and adolescents. Basel, Karger, 2003, pp 166-199 (22).

LITREATURE / LITERATUR

- 1) Baxter RC. 1986 The somatomedins:insulin-like growth factors. *Adv Clin Chem.*25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors.* New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology.* 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab.* 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
- 15) Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med.* 301:1138-1142
- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest.* 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence.* Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab.* 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh).*118:374-380.
- 22) Blum W., Schweizer R.: Insulin-like growth factors and their binding proteins; in Ranke MB (ed): *Diagnostics of endocrine function in children and adolescents.* Basel, Karger, 2003, pp 166-199.
- 24) Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGFBPs. *Growth Regulation* 4 (Suppl. 1):11-19
- 25) Trojan L, Bode C, Weiss C, Mayer D, Grobholz R, Alken P, Michel MS: IGF-II serum levels increase discrimination between benign prostatic hyperplasia and prostate cancer and improve the predictive value of PSA in clinical staging. *Eur Urol* (2006), 49(2):286-292
- 26) Soon PS, Gill AJ, Benn DE, Clarkson A, Robinson BG, McDonald KC, Sidhu SB: Microarray gene expression and immunohistochemistry analyses of adrenocortical tumors identify IGF2 and Ki-67 as useful in differentiating carcinomas from adenomas. *Endocr Relat Cancer* (2009), 16(2):573-583
- 27) Freude S et al: Neuronal IGF-1 resistance reduces Abeta accumulation and protects against premature death in a model of Alzheimer's disease. *FASEB J* (2009) Jun1 Epub ahead of print.

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SUMMARY – MEDIAGNOST IGF-II ELISA E30

Reconstitution / Dilution of Reagents		
Standards A-E	Reconstitution in Sample Buffer PP	500 µl
Control Serum KS1	Reconstitution in Sample Buffer PP	250 µl
Control Serum KS2	Reconstitution in Sample Buffer PP	250 µl
Wash Buffer WP	dilute in A. dest. (eg. total volume of 50 ml in a graduated flask and fill up to 1000 ml)	1:20
Sample + Control Sera KS1 and KS2: dilute 1:401 in Sample Buffer PP, mix immediately, incubate at least for 15 min, max. 2h. Use 50 µl for each well in the assay.		
Before conducting the assay equilibrate all reagents to room temperature.		

Assay Procedure for Double Determinations:

Pipette	Reagent	Position
50 µl	Antibody Conjugate AK	in all wells used
50 µl	Sample Buffer PP (blank)	A1 and A2
50 µl	Standard A (0.45 ng/ml)	B1 and B2
50 µl	Standard B (1.5 ng/ml)	C1 and C2
50 µl	Standard C (3 ng/ml)	D1 and D2
50 µl	Standard D (5.63 ng/ml)	E1 and E2
50 µl	Standard E (9 ng/ml)	F1 and F2
50 µl	Control Serum KS1	G1 and G2
50 µl	Control Serum KS2	H1 and H2
50 µl	Samples	following wells
Cover the wells with the sealing tape.		

Incubation: 2 h at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Enzyme Conjugate EK	each well

Incubation: 30 min at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Substrate S	each well

Incubation: 30 min in the dark RT

100 µl	Stop Solution SL	each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		

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CAL A-E	A -E	Rec in 500 µl BUF PP	
Control	KS1&KS2	Rec in 250 µl BUF PP	
WASHBUF 20x	WP		1:20 DILU A. dest.

SPE + Control 1:401 DILU BUF PP ↔ ⌚ 0.25 h
°C 20-25 °C

50 µl	Ab	A1 - End
50 µl	BUF PP	A1/2
50 µl	CAL A Std A (0.45 ng/ml)	B1/2
50 µl	CAL B Std B (1.5 ng/ml)	C1/2
50 µl	CAL C Std C (3 ng/ml)	D1/2
50 µl	CAL D Std D (5.63 ng/ml)	E1/2
50 µl	CAL E Std E (9 ng/ml)	F1/2
50 µl	CONTROL KS 1	G1/2
50 µl	CONTROL KS 2	H1/2
50 µl	SPE	
TAPE		

⌚ 2 h °C 20-25 ↔ 350 rpm

5x 300 µl	5x WASHBUF WP
100 µl	CONJ
TAPE	

⌚ 0.5 h °C 20-25 ↔ 350 rpm

5 x 300 µl	5x WASHBUF WP
100 µl	SUBST TMB S

⌚ 30 min °C 20-25

100 µl	H₂SO₄ SL
MEASURE	