HAVrealIDTECT

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Table of Contents

Symbols	2
1. Intended use.....	9
2. Principle.....	9
3. Specificity and Sensitivity	9
4. KIT- Contents	10
5. Storage.....	10
6. Required Materials Not Provided.....	10
7. Precautions	10
8. Principle of the real time PCR	10
9. Procedure of the one-step real time RT-PCR.....	11
10. Thermocycling protocol	12

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1. Intended use

The HAVrealDETECT Kit is suited for the detection of viral RNA from Hepatitis A Virus by means of real time Reverse Transcriptase Polymerase-Chain-Reaction (RT-PCR). The Kit is ready to use and contains a Mastermix and an Enzyme-Mix for the specific amplification of HAV-RNA extracted by an approved method from the sample material (serum or stool), for use with LightCycler® (Roche) as well as for Taqman® Platforms, (ABI or Biorad). The Mastermix contains a FAM labelled Taqman probe and an internal reference dye for fluorescent signal normalization which is detected in ROX channel. For quantification of the viral RNA the kit contains the quantification standards: QS1 to QS4 with 1.3×10^2 copies/reaction to 1.3×10^5 copies/reaction.

2. Principle

Hepatitis A is caused by the Hepatitis A Virus (HAV) and is responsible for about 50% of all viral caused inflammations of the liver in Germany. HAV is mostly transmitted faecal-orally, very rarely parenterally. The virus is excreted in stools and transmitted through contaminated food, especially raw seafood, ice cream, fruit, salad or polluted drinking water.

3. Specificity and Sensitivity

The primer/probe system was selected according to data bank analysis of the corresponding HAV-region of 100 published HAV-sequences. This system detects all human relevant HAV-genotypes.

Cross reactivity with organisms listed in the table below was NOT detected.

HIV Genotype 1a
HIV Genotype 1b
HCV Genotype 1a
HCV Genotype 1b
HEV
Coxsackie A7
Coxsackie A9
Echovirus 4
Echovirus 15
Echovirus 20
Norwalk Genotype II
Rotavirus

For determination of the sensitivity dilution series (0.3 IU/mL to 1000 IU/mL) of the NIBSC Standard (WHO International Standards for Hepatitis A Virus (HAV) RNA Nucleic Acid Amplification (NAT) assays, NIBSC code 00/560) were prepared. The analysed detection limit of the kit is ≥ 30 IU/mL.

4. KIT- Contents

	Product Code HAV-RD50
Content	Number of Reactions 50
2 x 425 µL	HAV-Master
1 x 50 µL	Enzyme-Mix
2 x 50 µL	QS1 1.3 x10 ² copies/reaction
2 x 50 µL	QS2 1.3 x10 ³ copies/reaction
2 x 50 µL	QS3 1.3 x10 ⁴ copies/reaction
2 x 50 µL	QS4 1.3 x10 ⁵ copies/reaction

5. Storage

The components of HAVrealDETECT have to be stored at –20°C and are stable to the expiry date given on the label. Please avoid repeated thawing and freezing (>2x). The Enzyme-Mix must be stored at –20°C.

6. Required Materials Not Provided

Powder-free gloves
RNA isolation kit
Adjustable pipettes
Sterile pipette tips with filter
Vortex-mixer
Real time PCR-instrument

7. Precautions

RNA is very sensitive to degradation through nucleases. Therefore use of sterile disposable materials for RNA handling is highly recommended. It is recommended to use an RNA isolation kit offered by various manufacturers for extraction of the nucleic acid from the sample material. Regarding the respective protocols a certain quantity of the sample material is purified and used in the reaction.

Sample preparation, set-up of RT/PCR and detection of the amplicons should be spatially separated.

All components should be used and stored regarding to the instructions.

8. Principle of the real time PCR

Specific regions of the organism to be determined are amplified by means of the RT-PCR. These genome regions are defined by choice of the primer/probe-systems. The detection by means of real time RT-PCR is carried out by fluorescence measurement, i.e. the more amplicon is generated the more fluorescence is released and can be measured. The intensity of the measured fluorescence is directly proportional to the quantity of the HAV-RNA to be determined in the sample material.

9. Procedure of the one-step real time RT-PCR

Number of Samples	1	1. Preparation of the Mastermix
HAV-Master	17 μ L	
Enzyme-Mix	1 μ L	
Total volume	18 μ L	

Preparation of the Mastermix

17 μ L HAV-Master per reaction and 1 μ L Enzyme-Mix per reaction are mixed in a sterile reaction tube (suitable for PCR).

Per reaction		2. Set-up of the RT-PCR reaction
Mastermix	18 μ L	
Sample	2 μ L	
Total volume	20 μ L	

Set-up of the RT-PCR


18 μ L Mastermix are added to 2 μ L sample in reaction tubes suitable for real time PCR (capillaries, plate, strips etc.). The reaction tubes are closed and transferred to the real time PCR cycle for the RT-PCR.

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Fig.1: Serial dilution of the quantification standard QS1-QS4 from 1.3×10^2 copies/reaction to 1.3×10^5 copies/reaction and the resulting standard curve.

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