Adiponectin ELISA

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ENGLISH **Instructions for use**

Adiponectin ELISA E09	96 Determinations
CE	DE/CA40/00809/18
Principle of the test	Enzyme-linked Immunoassay
Duration (incubation period)	1.75 h
Antibodies	monoclonal antibodies
Buffer	Ready for use and 20fold concentrate
Calibrator	5 single calibrators: 2 - 100µg/L, native human Adiponectin
Assay Range	0.27 – 31000 µg/L
Control	2 controls freeze-dried
Sample	human serum / plasma
Required sample volume	10 µL
Sample dilution	1:310
Analytical sensitivity	$\varnothing \leq 0.27 \mu\text{g/L}$
Intra- / Interassay Variance	$\varnothing < 10 \%$
Reference values	Langkamp M, Pridzun P, Böttner A, Kiess W, Thiery J, Kratzsch J, Reference Intervals for Adiponectin Levels in Human Serum. 2005 DGKL Annual Meeting Jena.

1 INTENDED USE

The ELISA E09 is intended to be used for quantitative measurement of human Adiponectin in human serum and plasma samples.

1 INTRODUCTION

Adiponectin is a 30kDa protein which percentage in serum proteins is 0.01%. It is mainly synthesized by adipocytes, but also muscle cells and hepatocytes have the ability to synthesize Adiponectin. It consists of a Collagen-like N-terminal and a globular C-terminal domain [1]. In vivo Adiponectin appears with different oligomers. Beside the trimer and dimer also high molecular multimers exist [1-3]. Two different receptors are known, both receptors are ubiquitarily expressed, though the distribution in the tissues varies. The Adiponectin Receptor 1 (AdipoR1) is especially in muscle- and AdipoR2 in liver tissue synthesized [4].

Several studies show, that adiponectin correlates negatively with BMI and thus it could have relevance for the energy metabolism for example through the regulation of fatty acid oxidation. Beside the correlation with BMI, Adiponectin level is associated with the Insulin-Resistance [5-7] and so also linked with Type II Diabetes. Adiponectin is associated also with glucose- und lipometabolism [8, 9]. The formerly proposed diagnostic value of the high molecular weight form of adiponectin was not verified using a commercially available test system for the determination of HMW adiponectin [10]. Blueher et al. evaluated the Mediagnost E09 regarding its diagnostic value in diagnosis of insulin resistance. Results of ROC analysis showed an area under curve of 0.92, which indicates a diagnostic value [10].

Furthermore adiponectin is involved in inflammatory processes [11-15] and therewith it is of importance for appearance of arteriosclerosis [4, 5, 16] and coronaritis [17, 18], thus the determination of Adiponectin level in plasma could serve to estimate the risk of coronary

Exemplary version, do not use to perform assays

disease [19, 20]. Beside this Adiponectin influences further physiological processes as for example the angiogenesis [21, 22].

2 ASSAY PRINCIPLE

The Mediagnost ELISA for Adiponectin E09 is a so-called Sandwich-Assay using two specific and highly affine antibodies. The Adiponectin in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-Adiponectin-Antibody binds in turn to the immobilised Adiponectin. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction, the turn of the colour will be catalysed quantitatively depending on the Adiponectin-level of the samples.

3 WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided. A Material Safety Data sheet is available on request.

Do not use obviously damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Controls CTR1 and CTR2, and Calibrators A-E.**

Source human serum for the controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore, all components and patient's specimens should be treated as potentially infectious.

Reagents A-E, DET, DIL, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution STP

The Stop solution contains 0.2 M acid sulphur acid (H₂SO₄)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

3.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

4 SAMPLES

4.1 Sample type

Serum and Plasma

Serum and heparin plasma levels are comparable. In EDTA- and Citrate Plasma-samples levels were found approx. 18% lower, because of the relatively high amount of anticoagulant.

4.2 Specimen collection

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. Hemolytic reactions have to be avoided.

4.3 Required sample volume: 10 µL

4.4 Sample stability

In firmly closable sample vials

- Storage at 20-25°C: 2 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles max. 3

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized, 3 freeze-thaw cycles showed no effect on the measured adiponectin concentration.

4.5 Interference

Hemoglobin, triglyceride and bilirubin in the sample do not interfere to a concentration of **5 mg/mL**, **100 mg/mL** and **100 µg/mL** respectively. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

4.6 Sample dilution

- Dilution: **1:310** with Dilution Buffer **DIL**
- Dilute for example **300 µL** Dilution Buffer **DIL** in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add **10 µL** Serum- or Plasma (dilution: 1:31). Add **900 µL** Dilution Buffer **DIL** in another PE-/PP-tube and **100 µL** of the thoroughly mixed first dilution. After mixing, use 2×100 µL from this **1:310** diluted sample in the assay.
- Alternatively, if you have the necessary technical equipment a one-step dilution of **1:301** is possible: Add **5 µL** to **1.5 mL** dilution buffer **DIL**.
- Sample stability after dilution of the sample: maximum 1 hour at 20-25°C.

5 MATERIALS

5.1 Materials provided

The reagents listed below are sufficient for 96 wells including the calibration curve.

MTP	Microtiter plate , ready for use, coated with mouse-anti-Adiponectin-antibody. Wells are separately breakable.	(8x12) wells
A-E	Calibrators , lyophilized, (native human Adiponectin), concentrations are given on vial labels and on quality certificate in ng/mL.	5 x 750 µL
CTR1	Control 1 , lyophilised, (human serum), concentration is given on quality certificate in ng/mL.	1 x 500 µL
CTR2	Control 2 , lyophilised, (human serum), concentration is given on quality certificate in ng/mL.	1 x 500 µL
DET	Antibody-HRP-Conjugate , ready for use, mouse-anti-hAdiponectin-antibody biotinylated + streptavidin horseradish peroxidase conjugate	1 x 12 mL
DIL	Dilution Buffer , ready for use	1 x 125 mL
WB	Washing Buffer , 20-fold concentrated solution	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine.	1 x 12 mL
STP	Stop Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x
	Instructions for use	1 x
--	Quality Certificate	1 x

5.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

6 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at –20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** calibrators **A-E** and Controls **CTR1** and **CTR2** must be stored at –20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Calibrators **A – E** and Controls **CTR1** and **CTR2** are reconstituted with the Dilution Buffer **DIL**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Controls **CTR1** and **CTR2** with the Dilution Buffer **DIL** in the same ratio (1:310) as the sample. The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

Assay Procedure

When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1**, **CTR2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **DET** as well as the succeeding Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**.

All determinations (Blank, Calibrators **A-E**, Controls **CTR1**, **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate**S**, stabilised H₂O₂-Tetramethylbenzidine, is photosensitive—store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

7 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 750 µL Dilution Buffer DIL	-
CTR1 and CTR2	Control 1 Control 2	in 500 µL Dilution Buffer DIL	1:310 with DIL
WB	Washing Buffer concentrate	-	1:20 with Aqua dest. → WB 1:20
Sample dilution: with Dilution Buffer DIL 1:310			
Before assay procedure bring all reagents to room temperature 20-25°C .			
Assay Procedure in Double Determination:			
Pipette	Reagents	Position	
100 µL	Dilution Buffer DIL (Blank)	A1/A2	
100 µL	Calibrator A (2 ng/mL)	B1/B2	
100 µL	Calibrator B (10 ng/mL)	C1/C2	
100 µL	Calibrator C (30 ng/mL)	D1/D2	
100 µL	Calibrator D (70 ng/mL)	E1/E2	
100 µL	Calibrator E (100 ng/mL)	F1/F2	
100 µL	Control CTR1 (1:310 diluted)	G1/G2	
100 µL	Control CTR2 (1:310 diluted)	H1/H2	
100 µL	Sample SPE (1:310 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
Sample Incubation: 1 h at 20-25°C, 350 rpm			
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WB 1:20/ well	In each well	
100 µL	Antibody-HRP-Conjugate DET	In each well	
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at 20-25°C, 350 rpm			
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WB 1:20/ well	In each well	
100 µL	Substrate Solution S	In each well	
Incubation: 15 Minutes in the Dark at 20-25°C			
100 µL	Stop Solution STP	In each well	
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			

8 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All calibrators and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

8.1 Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of calibrator E should be above 1.00.

Samples, which yield higher absorbance values than calibrator E, should be re-tested with a higher dilution.

9 EVALUATION OF RESULTS

9.1 Establishing of the calibration curve

Calibrator	A	B	C	D	E
ng/mL	2	10	30	70	100

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The Adiponectin concentration in ng/mL of the samples and controls **CTR1** and **CTR2** can be calculated by **multiplication** with the respective **dilution factor**.

9.2 Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

	Blank	A	B	C	D	E
ng/mL	0.0	2	10	30	70	100
OD(450-620 nm)	0.008	0.071	0.357	1.022	2.053	2.817


Figure 1 Exemplary calibration curve

The exemplary shown calibration curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

9.3 Exemplary calculation of Adiponectin concentrations

Sample dilution: 1:310

Measured extinction of your sample	0.408
Measured extinction of the blank	0.008

Your measurement programm will calculate the Adiponectin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2nd degree).

In this exemplary case the following equation is solved by the program to calculate the adiponectin concentration in the sample:

$$0.400 = -6 \times 10^{-5}x^2 + 0.0344x$$
$$11.89 = x$$

If the dilution factor (1:310) is taken into account the adiponectin concentration of the undiluted sample is

$$11.89 \times 310 = 3685.9 \text{ ng/mL} = 3.69 \text{ } \mu\text{g/mL}$$

The exemplary shown calibration curve in Fig.1 cannot be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

9.4 Interpretation of results

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

9.5 Limitation of procedure

The Mediagnost sensitive human Adiponectin ELISA, E09 is based on monoclonal antibodies. Generally, this technique could be sensible to heterophilic antibodies or rheumatic factors in the sample. Their influence is reduced by assay design, but cannot be excluded completely.

10 REFERENCE VALUES

The expected values for serum adiponectin, which were determined with the Mediagnost ELISA E09 in healthy donors and analysed by Prof. Dr. J. Kratzsch, Department of Laboratory Medicine, University Hospital Leipzig, are given below (Tab. 1).

These data show significant correlation between Adiponectin-Serum values and age as well as gender of the probands, in turn the correlation between the respective BMI seems to be less significant. In the samples of neonatal cord blood very high values were found. Several different statistical analyses were performed to adapt for certain individual demands. The best suited data can be chosen respectively for interpretation of the own measurements.

Table 1a Expected values for adults, gender specific mean as well as median, 5. and 95. percentiles are given.

Sex	Number [n]	Mean [µg/ml]	Median [µg/ml]	Standard Deviation	5 th Percentile [µg/ml]	95 th Percentile [µg/ml]
Female	101	10.2	9.1	4.6	4.0	19.4
Male	125	6.8	6.1	4.1	2.0	13.9
total	226	8.3	7.5	4.6	2.4	19.3

Table 1b Expected values for children, gender specific mean as well as median, 5. and 95. percentiles are given.

Sex	Number [n]	Mean [µg/ml]	Median [µg/ml]	Standard Deviation	5 th Percentile [µg/ml]	95 th Percentile [µg/ml]
Female	131	8.71	8.18	4.32	3.05	15.6
Male	134	8.97	8.12	5.13	3.36	18.6
total	265	8.84	8.18	4.74	3.33	16.5

Table 1c Expected values for Adiponectin, age specific mean as well as median, 5. and 95. percentiles are given.

Age group [a]	Number [n]	Mean [µg/ml]	Median [µg/ml]	5. Percentile [µg/ml]	95. Percentile [µg/ml]
< 7.99	46	12.82	11.7	2.33	26.5
8 – 9.99	40	8	8.09	3.96	14.9
10-11.99	55	8.02	7.14	3.36	13.8
12 – 13.99	26	8.21	7.54	4.5	13.2
14 – 15.99	59	8.12	8.09	3.67	13.7
16 – 19.99	41	7.97	7.79	2.74	13.3
alle	267	8.88	8.18	3.33	16.7
20 – 29.99	47	6.72	6.38	2.50	12.25
30 – 39.99	38	7.38	6.69	1.98	19.29
40 – 49.99	55	8.42	8.20	2.41	17.85
50 – 59.99	47	9.61	8.55	2.15	19.85
> 60	32	9.52	8.57	3.00	21.10
alle	226	8.33	7.5	2.41	19.29

Table 1d Expected values for Adiponectin, age as well as gender specific mean and median, BMI and 25. and 75. percentiles are given.

Female			Adiponectin (µg/ml):			
Age (Years):	n:	BMI: AV ± SD	AV ± SD:	Median:	Percentile: 25.- 75.	Min. – Max.:
Newborn Cord blood	19	--	29.80 ± 12.49	26.1	19.5 -35.2	16.9 - 61.4
< 3.99	9	15.73 ± 0.79	14.43 ± 7.76	11.2	8.2 - 21.8	2.3 - 26.7
4.0 - 7.99	11	16.01 ± 1.94	8.46 ± 4.73	9.3	2.9 - 12.1	1.4 - 15.6
8.0 - 9.99	22	17.58 ± 3.84	7.92 ± 3.00	8.2	5.2 - 10.0	3.6 - 15.1
10.0 - 11.99	33	17.83 ± 1.86	7.66 ± 4.59	6.6	5.0 - 8.8	3.1 - 20.9
12.0 - 13.99	11	19.85 ± 2.31	8.22 ± 5.64	7.5	6.5 - 9.2	4.9 - 13.2
14.0 - 15.99	27	19.91 ± 1.72	8.83 ± 9.25	8.9	5.2 - 11.8	2.6 - 17.7
16.0 – 19.99	18	21.64 ± 2.64	9.00 ± 3.22	8.7	6.9 - 11.2	2.7 - 14.0
20.0 - 29.99	24	23.12 ± 5.01	7.39 ± 3.35	7.3	5.7 - 9.0	3.4 - 17.8
30.0 - 39.99	17	23.20 ± 2.86	9.19 ± 3.89	8.6	7.2 - 10.4	3.6 - 19.3
40.0 - 49.99	26	24.50 ± 4.11	9.93 ± 3.59	9.5	7.5 - 11.6	4.4 - 19.6
50.0 - 59.99	21	24.61 ± 3.31	11.5 ± 5.49	10.0	8.0 - 15.9	2.0 - 23.1
>60.0	8	24.63 ± 1.89	15.6 ± 4.64	15.3	11.4 - 18.2	11.2 - 24.1

Male			Adiponectin (µg/ml):			
Age (Years):	n:	BMI: AV ± SD	AV ± SD:	Median:	Percentile: 25.- 75.	Min. – Max.:
Newborn Cord blood	10	--	27.80 ± 7.68	26.7	22.2 -31.0	15.6 - 40.6
< 3.99	14	16.17 ± 1.81	16.57 ± 6.55	14.3	11.6 - 21.2	5.8 - 40.3
4.0 - 7.99	12	15.69 ± 1.05	11.24 ± 5.43	9.7	8.9 - 15.9	3.5 - 18.6
8.0 - 9.99	18	16.45 ± 1.76	8.11 ± 2.93	7.6	6.2 - 9.1	5.00 - 15.4
10.0 - 11.99	21	18.34 ± 2.18	8.43 ± 3.91	7.8	5.2 - 10.9	3.4 - 20.2
12.0 - 13.99	14	18.61 ± 2.11	7.59 ± 2.86	7.1	6.0 - 10.3	2.4 - 12.2
14.0 - 15.99	32	19.86 ± 2.00	7.53 ± 2.52	7.4	5.1 - 9.3	3.8 - 15.4
16.0 -19.99	23	22.03 ± 2.42	7.16 ± 3.53	6.9	4.2 - 9.6	2.0 - 13.9
20.0 - 29.99	23	23.43 ± 2.48	5.44 ± 2.29	5.8	4.0 - 6.9	1.3 - 10.3
30.0 - 39.99	21	23.33 ± 2.72	5.92 ± 4.60	4.4	2.7 - 6.7	1.9 - 20.6
40.0 - 49.99	22	23.79 ± 2.41	6.13 ± 2.92	5.5	3.8 - 8.3	2.1 - 11.6
50.0 - 59.99	23	26.68 ± 2.77	7.45 ± 4.50	6.7	5.0 - 8.8	1.4 - 19.6
>60.0	24	25.72 ± 2.12	7.48 ± 3.92	7.6	4.6 - 9.2	3.0 - 21.1

n= Number of Proband AV=Average Value BMI=Body Mass Index (kg/m²) SD=Standard Deviation

11 PERFORMANCE CHARACTERISTICS

11.1 Sensitivity

The analytical sensitivity of the ELISA E09 was measured by the variability of the signal of the blank. Based on the twofold standard deviation of the blank the mean analytical sensitivity is < 0.27 ng/mL (Range 0.094 to 0.59 ng/mL).

11.2 Specificity

Adiponectin exists in different oligomeric forms: the high, medium and low molecular weight form. Different numbers of the adiponectin monomer aggregate specifically to form a complex. In Figure 2a the five different forms of human adiponectin are shown schematically. In parallel the results of a size-exclusion chromatography of human serum measured with the Mediagnost E09 Adiponectin ELISA are shown.


Figure 2 Adiponectin Structure (A) proteinaceous structure of human adiponectin including posttranslational modifications and multimeric forms (taken from Nakano et al¹). (B) Results of a serum separation by size-exclusion chromatography. The sample was fractionated and the adiponectin content of each fraction was measured by Mediagnost E09 Adiponectin, for comparison the corresponding seize in kDa is indicated by black arrows.

The results shown in Figure 2b clearly demonstrate that the Mediagnost E09 Adiponectin ELISA detects all forms of Adiponectin present in human serum: the trimer at 65 kDa, the hexamer at 150 kDa and the high molecular weight forms of >280 kDa. The Mediagnost E09 Adiponectin ELISA therefore measures total Adiponectin.

11.3 Precision

Intra-Assay Variance & Accuracy

Intra assay variance and accuracy is exemplarily shown with two samples. The adiponectin concentration of these samples was repeatedly measured in one assay.

Table 2 Intra-Assay Variation. Recombinant adiponectin was diluted in dilution buffer and the adiponectin concentration of the dilution was measured repeatedly within one assay.

	Determinations [n]	Mean value [µg/L]	Standard deviation [µg/L]	VC [%]	Target Value [µg/L]
Sample 1	8	7.108	0.22	3.14	6
Sample 2	8	107.96	3.97	3.67	100

In both samples the variability is less than 5% and the deviation from the target value is <20%.

¹ Nakano Y, Tahima S, Yoshimi A, Akiyama H, Tsushima M, Tanioka T, Negoro T, Tomita M, Tobe T: A novel enzyme-linked immunosorbent assay specific for high-molecular-weight adiponectin. J Lipid Res 2006, 47(7):1572-1582
IFU E09 IVD D/E

Inter-Assay Variance

Serum samples were repeatedly measured in independent assays of different lots. On average the coefficient of variation was 7.5% (SD 1.6). The results of 5 samples are shown in table 3.

Table 3 Inter-Assay Variation

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
Mean [µg/mL]	4.72	5.25	8.36	5.45	22.29
CV [%]	8.16	8.14	6.93	8.05	7.30
n	99	68	62	174	62

11.4 Linearity

Linearity of sample dilution was tested by serial dilution (1:100 – 1:4000) of human serum samples and recalculation of the adiponectin content in comparison to the mean adiponectin concentration of all dilutions (Table 4), no diluted sample showed a deviation of >30%.

Table 4 Linearity. Human serum sample were diluted in DIL and adiponectin content was recalculated. Measurements results are shown in [mg/L]. No deviation of the mean >30% was detected.

µg/mL	Mean	1:100	1:200	1:400	1:800	1:1000	1:2000	1:4000
Sample 1	5.76	6.53	6.331	5.764	5.49	6.067	6.114	4.056
Sample 2	11.53	10.93	12.107	11.395	11.454	11.567	12.884	10.362
Sample 3	12.07	13.57	12.853	12.03	11.974	11.338	11.548	11.169
Sample 4	4.89	4.659	4.886	4.384	4.425	5.851	5.13	n/a

Additionally, dilutions of 1:1500 to 1:96000 were evaluated with two samples. In Figure 3 the results are shown and demonstrate that in the tested samples no effect of dilution could be detected on measured adiponectin concentrations. The deviation of the target concentration of each dilution was less than 30% in all dilutions.


Figure 3 Linearity. Adiponectin concentration was measured in two human samples diluted 1:1500 to 1:96000.

11.5 Recovery

Trueness and traceability of the Mediagnost Adiponectin ELISA E09 was evaluated by recovery of recombinant adiponectin in human serum. The recovery of recombinant Adiponectin yielded in a serum matrix on average 110%.

Table 5 Recovery of recombinant human Adiponectin in Serum. Recombinant Adiponectin was added in different amounts to human serum. The Adiponectin content of the so enriched samples was measured and recovery in comparison to enriched buffer calculated.

R&D recombinant Adiponectin lot 1022911	DIL	Serum	Recovery
ng/mL	ng/mL	ng/mL	%
0	0	0.00778	---
75	87.92	95.71	109
37.5	51.84	59.98	116
18.75	26.55	26.7	101
9.375	13.35	15.15	113

11.6 Interference

Interference of bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing adiponectin. For comparison the same amount of buffer without any substance was also added to the serum. Table 6 demonstrates that neither bilirubin nor triglycerides or haemoglobin exert any influence on the measurement of adiponectin in human serum.

Table 6 Interference. Serum samples were enriched with different amount of triglycerides, bilirubin or haemoglobin. The relative amount of adiponectin measured in comparison with native serum is shown here [%].

Triglyceride 100 mg/mL	Bilirubin 100 µg/mL	Hemoglobin 1 mg/mL	Hemoglobin 5 mg/mL
94	96	90	109
90	93	97	--
95	94	93	--

12 COMPARISON STUDIES

Mediagnost Adiponectin, E09 was compared with two different, commercially available test systems. In linear regression analysis both tests showed a good coefficient of determination ($R^2=0.896$ and $R^2=0.97$). Thus, the comparability of the results between the Mediagnost E09 and other test systems is possible. In dependence on the respective system the absolute deviation of measurement results is different, but because of the excellent correlation the results are comparable after applying a factor. The results of both studies are shown in figure 4.


Figure 4 Comparison of Mediagnost E09 with commercially available test systems (left) radioimmunoassay (n=14) and (right) ELISA (n=84)

Instructions for use for scientific application

13 SCIENTIFIC APPLICATION

In addition to serum and plasma samples Adiponectin can be determined in other human body fluids and in cell culture supernatants of human cell lines for research purposes.

13.1 Samples suitable for scientific application

Serum, plasma, saliva, urine, breast milk and cell culture supernatant of human cell lines

The recommended dilution for serum and plasma samples in Dilution Buffer **DIL (1:310)**.

In the other samples, the Adiponectin levels can vary considerably, the optimal dilution must be found out by the customer.

Table 7 Results of sample matrix tests . Adiponectin was added to the respectively diluted samples. Enriched samples were measured without further dilution. Shown is the relative recovery of added Adiponectin of the concentration measured in respectively enriched assay buffer.

Matrix Dilution	Urine	Saliva	Breast Milk	Cell Culture Medium 10% FCS	Cell Culture Medium
1:2	80	87	89	83	95
1:5	95	80	92	92	97
1:10	92	87	-	101	85
1:20	94	99	-	83	91

13.2 Species Cross-Reactivity

Serum of the cited species was diluted and used as sample in this assay system. No signal was detected in serum of the following species:

Horse, Cow, Chicken, Rabbit, Dog, Guinea pig, Sheep, Mouse, Goat, Donkey, Rat, Cat

Whether this obvious non-reactivity is species specific should be assessed individually by each customer.

14 LITERATUR / REFERENCES

Literature

1. Nakano, Y., et al., Isolation and characterization of GBP28, a novel gelatin-binding protein purified from human plasma. *J Biochem (Tokyo)*, 1996. 120(4): p. 803-12.
2. Pajvani, U.B., et al., Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity. *J Biol Chem*, 2003. 278(11): p. 9073-85.
3. Tsao, T.S., et al., Role of disulfide bonds in Acrp30/adiponectin structure and signaling specificity. Different oligomers activate different signal transduction pathways. *J Biol Chem*, 2003. 278(50): p. 50810-7.
4. Shimada, K., T. Miyazaki, and H. Daida, Adiponectin and atherosclerotic disease. *Clin Chim Acta*, 2004. 344(1-2): p. 1-12.
5. Higashiura, K., et al., Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations. *Clin Endocrinol (Oxf)*, 2004. 61(6): p. 753-9.
6. Spranger, J., et al., Adiponectin is independently associated with insulin sensitivity in women with polycystic ovary syndrome. *Clin Endocrinol (Oxf)*, 2004. 61(6): p. 738-46.
7. Zoico, E., et al., Adipocytokines, fat distribution, and insulin resistance in elderly men and women. *J Gerontol A Biol Sci Med Sci*, 2004. 59(9): p. M935-9.
8. Ye, J.M., et al., PPARalpha /gamma ragaglitazar eliminates fatty liver and enhances insulin action in fat-fed rats in the absence of hepatomegaly. *Am J Physiol Endocrinol Metab*, 2003. 284(3): p. E531-40.
9. Yamauchi, T., et al., Adiponectin stimulates glucose utilization and fatty-acid oxidation by activating AMP-activated protein kinase. *Nat Med*, 2002. 8(11): p. 1288-95.
10. Blüher (Blueher), M., et al., Total and high-molecular weight adiponectin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays. *Diabetes Care*, 2007. 30(2): p. 280-5.
11. Delaigle, A.M., et al., Induction of adiponectin in skeletal muscle by inflammatory cytokines: in vivo and in vitro studies. *Endocrinology*, 2004. 145(12): p. 5589-97.
12. Winzer, C., et al., Plasma adiponectin, insulin sensitivity, and subclinical inflammation in women with prior gestational diabetes mellitus. *Diabetes Care*, 2004. 27(7): p. 1721-7.
13. Xydakis, A.M., et al., Adiponectin, inflammation, and the expression of the metabolic syndrome in obese individuals: the impact of rapid weight loss through caloric restriction. *J Clin Endocrinol Metab*, 2004. 89(6): p. 2697-703.
14. Motoshima, H., et al., Adiponectin suppresses proliferation and superoxide generation and enhances eNOS activity in endothelial cells treated with oxidized LDL. *Biochem Biophys Res Commun*, 2004. 315(2): p. 264-71.
15. Wolf, A.M., et al., Adiponectin induces the anti-inflammatory cytokines IL-10 and IL-1RA in human leukocytes. *Biochem Biophys Res Commun*, 2004. 323(2): p. 630-5.
16. Okamoto, Y., et al., Adiponectin reduces atherosclerosis in apolipoprotein E-deficient mice. *Circulation*, 2002. 106(22): p. 2767-70.
17. Schlegel, A., Adiponectin and risk of coronary heart disease. *Jama*, 2004. 292(1): p. 40; author reply 40.
18. Choi, K.M., et al., Inflammation, Insulin Resistance, and Glucose Intolerance in Acute Myocardial Infarction Patients without a Previous Diagnosis of Diabetes Mellitus. *J Clin Endocrinol Metab*, 2004.
19. Nakamura, Y., et al., Implications of plasma concentrations of adiponectin in patients with coronary artery disease. *Heart*, 2004. 90(5): p. 528-33.
20. Pischon, T., et al., Plasma adiponectin levels and risk of myocardial infarction in men. *Jama*, 2004. 291(14): p. 1730-7.
21. Shibata, R., et al., Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling. *J Biol Chem*, 2004. 279(27): p. 28670-4.
22. Fernandez-Real, J.M., et al., Adiponectin is associated with vascular function independent of insulin sensitivity. *Diabetes Care*, 2004. 27(3): p. 739-45.

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