# Anti-HAV ELISA

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## Table of Contents

1.	INTENDED USE	20
2.	INTRODUCTION	20
3.	ASSAY PRINCIPLE	21
4.	WARNINGS AND PRECAUTIONS	22
4.1.	General first aid procedures:	22
5.	SAMPLES	23
5.1.	Sample type	23
5.2.	Specimen collections	23
5.3.	Required sample volume:	23
5.4.	Sample stability	23
5.5.	Interference	23
5.6.	Sample dilution	23
6.	MATERIALS	24
6.1.	Materials provided	24
6.2.	Materials required, but not provided	24
7.	TECHNICAL NOTES	25
8.	ASSAY PROCEDURE	26
9.	QUALITY CONTROL	27
9.1.	Quality criteria	27
10.	EVALUATION OF RESULTS	27
10.1.	Qualitative test analysis	27
10.2.	Semi-Quantitative test analysis	28
10.2.1.	Antibody titre calculation	28
10.2.2.	Analysis with Mediagnost standards	28
10.2.3.	Analysis with reference preparation	30
11.	LIMITATION OF PROCEDURE	31
12.	PERFORMANCE CHARACTERISTICS	31
12.1.	Assay Calibration	31
12.2.	Sensitivity	31
12.3.	Specificity	32
12.4.	Precision Data	32
12.5.	Linearity	33
13.	LITERATUR / LITERATURE	34
14.	International Test description for QUALITATIVE Assay	36

## ENGLISH Instructions for use

Anti-HAV ELISA E10	96 Determinations
Regulatory Status	For Research Use Only. Not for diagnostic purposes
Principle of the test	Sandwich ELISA
Duration (incubation period)	3.5 h
Conjugate Concentrate	100-fold concentrated
Buffer and Substrate	Ready for Use
Reference material	WHO NIBSC 97/646 anti-hepatitis A immunoglobulin, WHO-Standard
Standard	3 Single Standards: ready for use, 50 mIU/mL, 30 mIU/mL, 10 mIU/mL
Assay Range	0 – 50 IU/L
Control	Ready for use
Sample	humanes Serum / Plasma
Required sample volume	25 µL
Sample dilution	1:10
Analytical sensitivity	1.10 mIU/mL
Inter-Assay Variance	< 30 %

### 1. INTENDED USE

Mediagnost anti-HAV ELISA, E10 is an enzyme immunoassay for the detection of antibodies against Hepatitis A virus in Human Serum **for research use**.

### 2. INTRODUCTION

After infection with the Hepatitis A-Virus, neutralising antibodies appear at the same time of Anti-HAV of IgG-Class formation (1). The Anti-HAV titers 3 to 6 months after naturally occurring infections are very high, within the range of 100 to more than 300 International Units per ml (IU/ml). Even after more than 10 years the titers usually remain at more than 1 to 10 IU/ml (2,3).

A value of 10 - 30 (mIU/ml) milli International Units can be considered the minimal protective level (2, 4). Since the incidence of HAV infection in children has diminished in Northern Europe in recent years, children and juveniles are predominantly Anti-HAV negative. In an adult population could be observed, that the percentage of anti-HAV positive individuals raises with increasing age: only 3.9 % of the young people (18-24 years) show antibody against HAV, where as in older population (≥50 years) HAV-Antibody was found in 40.3 % of the population (12).

In countries with less favourable sanitary conditions around the Mediterranean area, Africa or Asia, the incidence is very high. (5). Higher risk is also given for clinical staff, especially in the paediatric (6).

Since a vaccine against Hepatitis A virus infections is available, vaccinations are recommended for people travelling to countries where a high risk of HAV infections exists and for health care employees.

### 3. ASSAY PRINCIPLE

Mediagnost anti-HAV ELISA, E10 is a pseudo-competitive enzyme immunoassay. Serum or plasma samples are added to the wells of a microtiter plate, which have been previously coated with inactivated HAV antigen, and incubated for 2 hours at 37 °C. Anti-HAV antibodies bind to the antigen. The conjugate (peroxidase labeled anti-HAV) is added and incubated again for 1 h at 37 °C. Free binding sites of the antigen are bound with conjugate. Excess conjugate is washed of the plate and the substrate is added and incubated for 30 min at room temperature. The bound conjugate changes the colour of the substrate to blue. The reaction is terminated by adding the stopping solution. The colour turns yellow. The absorbance of the coloured reaction product is measured on a microtiter plate reader. The extinction is reciprocal to the anti-HAV titer.

For **semi-quantitative determination** use the included serum standards.

The **preparation of titration curve** e.g. for calibration of sera by means of standard reagents is also possible.

## 4. WARNINGS AND PRECAUTIONS

### For In Vitro Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro use and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

### Human Serum

Following components contain human material: **PK, NK, KK, STD 1-3**

The test plate MTP is coated with inactivated Antigen, which was negative in the rest infectivity test.

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

### Reagents KK, VP, WP

Contain as preservative 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

### Substrate Solution (S)

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

### Stopping Solution (SL)

The Stopping solution contains 0.2 M acid sulphur acid (H<sub>2</sub>SO<sub>4</sub>)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

#### 4.1. General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

## 5. SAMPLES

### 5.1. Sample type

Serum and Plasma

Serum and Heparin / EDTA-Plasma yield comparable values.

### 5.2. Specimen collections

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

### 5.3. Required sample volume:

25 µL

### 5.4. Sample stability

In firmly closable sample vials

- Storage at 4°C: max. 3 days
- Freeze-thaw cycles: max. 3

Freezing and thawing of samples should be minimized. 3 Freezing-thawing showed no effect on samples.

### 5.5. Interference

Anticoagulants like EDTA and heparin in the usual concentrations do not influence the test. Bilirubin concentrations do not interfere up to 200 µg/mL and triglyceride concentrations up to 5 mg/mL. Strongly hemolyzed samples (hemoglobin > 10 mg/mL) can result in false negative or false low values. Strongly hemolyzed samples should not be used.

### 5.6. Sample dilution

- Qualitative Test: Dilution: **1:10** with Dilution Buffer **VP**
- Quantitative Test: Dilution **1:10** with Dilution Buffer **VP**, > **1:10** with **VPN**

Use a 10 % solution of **anti-HAV antibody negative serum** in **Dilution Buffer VP (VPN)** for all sample dilutions > **1:10**

Example: Dilution 1:10:

**25 µL Sample** is added to **225 µL Dilutionbuffer VP**.

Example Dilution 1:50:

**1 mL anti-HAV Antibody negative Serum** is mixed with **9 mL Dilutionbuffer VP (VPN)**.

**10 µL Sample** is added to **490 µL** to this Dilution Buffer including negative Serum (**VPN**).

## 6. MATERIALS

### 6.1. Materials provided

<b>MTP</b>	<b>Microtiter plate</b> , ready for use, coated with inactivated HAV-antigen. Wells are separately breakable.	<b>(8x12) wells</b>
<b>KK</b>	<b>Conjugate Concentrate KK, 100fold concentrated</b> (peroxidase labeled Anti-HAV Ig)	<b>1 x 100 µL</b>
<b>PK</b>	<b>Positive Control, PK</b> , anti-HAV positive control serum >500 mIU/mL, ready for use.	<b>1 x 1 mL</b>
<b>NK</b>	<b>Negative Control, NK</b> anti-HAV negative control serum, ready for use	<b>1 x 1 mL</b>
<b>VP</b>	<b>Dilution Buffer</b> , ready for use, <b>Please shake before use!</b>	<b>1x 120 mL</b>
<b>WP</b>	<b>Washing Buffer</b> , 20-fold concentrated solution	<b>1 x 50 mL</b>
<b>S</b>	<b>Substrate</b> , ready for use, horseradish-peroxidase-(HRP) substrate	<b>1 x 12 mL</b>
<b>SL</b>	<b>Stopping Solution</b> , ready for use, 0.2 M sulphuric acid.	<b>1 x 12 mL</b>
-	<b>Sealing Tape</b> , for covering the <b>microtiter plate</b> .	<b>2 x</b>
	<b>Instructions for use</b>	<b>1 x</b>
--	<b>Quality Control Certificate</b>	<b>1 x</b>

#### Additional to Semi-Quantitative Test:

<b>STD1</b>	<b>Serum standard 1</b> (anti-HAV titer 50 mIU/mL) (STD1) ready for use.	<b>1 x 1 mL</b>
<b>STD2</b>	<b>Serum standard 2</b> (anti-HAV titer 30 mIU/mL) (STD 2) ready for use	<b>1 x 1 mL</b>
<b>STD3</b>	<b>Serum standard 3</b> (anti-HAV titer 10 mIU/mL) (STD 3) ready for use	<b>1 x 1 mL</b>

### 6.2. Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WP (A.dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- **Anti-HAV antibody negative serum** (for dilutions > 1:10 to the semi-quantitative measurement)
- Incubator (suitable for incubations at **37°C**)
- Vortex-mixer
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥590 nm

## 7. TECHNICAL NOTES

### Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. Avoid repeated thawing and freezing.

### Storage Life

The shelf life of the components after initial opening is warranted for 4 weeks, store the unused strips and microtiter wells airtight together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The 1:20 diluted Washing Buffer WP is 4 weeks stable at 2-8°C.

### Preparation of reagents

Bring all reagents to room temperature (20 – 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Dilution

Use the Dilution Buffer **VP** for the dilution of the Conjugate Concentrate **KK**. Prepare only the amount required. The diluted conjugate can be stored for at least one week at 4°C.

The required volume of Washing Buffer **WP** is prepared by 1:20 dilution of the provided 20fold concentrate with **Aqua dest.** Please dilute only according to daily requirements.

### Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate Solution **S**, is photosensitive—store and incubation in the dark.

### Assay Procedure

When performing the assay reagents and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, diluted Conjugate Concentrate **KK**, as well as the Substrate Solution **S** should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution **SL** should be added to the plate in the same order as Substrate Solution **S**. All determinations should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

### Washing

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WP** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter plate washer**, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

## 8. ASSAY PROCEDURE

Reagent Preparation			
Before use bring all reagents to room temperature: <b>20-25 °C</b>			
<b>KK</b>	Conjugate Concentrat	in Dilution Buffer <b>VP</b> (e.g. 10 µL KK + 990 µL <b>VP</b> )	<b>1:100</b>
<b>WP</b>	Washing Buffer	in Aqua dest. (e.g 50 mL WP + 950 mL <b>A.dest</b> )	<b>1:20</b>
-	Samples	In Dilution Buffer <b>VP</b> (e.g. 25 µL Sample + 225 µL <b>VP</b> )	<b>1:10</b>
<b>VPN</b>	<b>VP +10% negative Serum*</b>	In Dilution Buffer <b>VP</b> (z.B. 1 mL negative Serum + 9 mL <b>VP</b> )	-
-	Samples	in Dilution Buffer <b>VP</b> + Serum <b>VPN</b> (e.g. 10 µL Sample + 490 µL <b>VPN</b> )	<b>&gt;1:10</b>
Assay Procedure in Doble Determination			
Pipette	Reagents		Position
100 µL	-		A1/A2
100 µL	Positive Control		B1/B2
100 µL	Negative Control = Standard 4 0 mIU/mL		C1/C2
100 µL	Standard 1 50 mIU/mL (only in semi-quantiative measurement)		D1/D2
100 µL	Standard 2 30 mIU/mL (only in semi-quantiative measurement)		E1/E2
100 µL	Standard 3 10 mIU/mL (only in semi-quantiative measurement)		F1/F2
100 µL	Samples e.g. 1:10 diluted		according to demand
Cover the wells with the sealing tape.			
Incubation: 2 h at 37°C			
50 µL	1:100 diluted Conjugate Concentrate		from B1 on in each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash <b>3x</b> with <b>300 µL</b> Washing Buffer <b>WP</b> / well		in each well
100 µL	Substrate Solution <b>S</b>		in each well
Substrate S Incubation: 30 minutes in the dark at 20-25°C			
100 µL	Stopping Solution <b>SL</b>		in each well
Measure the absorbance within 30 min at <b>450 nm</b> (reference filter ≥ 590 nm).			

\*not included in the test

Instead of the standards contained in the test, an internal standard can be used. This may be the dilution of a reference preparation, or prepared by titration of serum against a reference preparation.

**Reference preparations** are available from the Paul-Ehrlich-Institute, (Postfach 17 40, 63207 Langen, Tel. 06103/770, Telefax: 06103/77123, anti-HAV Reference sample, "PEI-Standard") or from the Central laboratory of the Netherlands Red Cross (Plesmanlaan 125, 1006 CX Amsterdam, Telefax 31 20 5123474, WHO International Reference Preparations, 97/646 anti-hepatitis A immunoglobulin, WHO-Standard).

## 9. QUALITY CONTROL

GLP requires that controls be run with each calibration curve. A statistically significant number of controls should be frequently assayed to establish mean values and acceptable ranges to assure proper performance.

### 9.1. Quality criteria

In order to evaluate the results it should be ensured that the absorbances of the blank (**A1/A2**) **0.1** do not exceed the Positive Control **PK 0.15** OD<sub>450</sub> units. The signal of the Negative Control **NK** must be greater than **1.0** OD<sub>450</sub> units. The difference between the absorbance of the Negative Control **NK** and Positive Control **PK** must be at least **0.4** OD<sub>450</sub> units.

The absorbances of the samples should be within the standard curve (**NK - STD3**) in the semi-quantitative method. If the absorbances are outside this range, for reliable determinations, they should be determined again at higher dilutions in a second test.

## 10. EVALUATION OF RESULTS

First, the absorbance determined with the chosen reference wavelength (> 590 nm) is subtracted from each absorbance value at 450 nm, regardless of the chosen evaluation method. Following this, the absorbance of the blanks (A1/A2) is also subtracted from all values. The measured values thus obtained are the basis for all further analyzes.

### 10.1. Qualitative test analysis

For the qualitative determination of anti-HAV antibodies in human serum samples a method was created in the context of product development, which enables the distinctions of positive and negative probands with a sensitivity of 98.26 % (n = 801).

For this purpose a cut-off value is calculated from the signal of the positive control and negative control according to the following formula.

$$cut - off = \frac{\text{Extinction Negative Control} + \text{Extinction Positive Control}}{2}$$

Samples with mean absorbance values **higher than the cut-off value** are considered to be **negative**. Samples with **absorbance values less than** the cut-off are considered to be **positive**. Samples  $\pm 10\%$  around the cut-off should be determined again.

An exemplary test result is shown in Table 1. The absorbance difference between negative and positive control is greater than 0.4 (1.390 to 0.024). The calculation of the cut-off value gives the following result.

$$cut - off = \frac{NK + PK}{2} = \frac{1.390 + 0.024}{2} = 0.707$$

The cut off value is therefore 0.707. Thus, all samples whose signal is higher than this are classified as negative (containing no anti-HAV antibody) and all samples in which the signal is lower than 0.707 are assessed as positive (containing anti-HAV antibody).

The absorbance of serum sample 1 is greater than the cut-off value, it is thus negative. The absorbance of serum sample 2 is less than the cut-off value, it is positive. In serum sample 3, the absorbance value is in the range of 0.778 to 0.636 (= 0.707 + 10 %), the determination must be repeated in the above example.

**Table 1:** Examples of absorbances of anti-HAV determinations using the qualitative test methods are shown as the difference between the absorbances at 450 nm signal and the signal at 620 nm as a reference.

Sample:	Extinction	Mean
Negative Control	1.372	1.390
Negative Control	1.408	
Positive Control	0.024	1.387
Serum Sample 1	1.461	
Serum Sample 1	1.312	0.023
Serum Sample 2	0.025	
Serum Sample 2	0.021	0.739
Serum Sample 3	0.735	
Serum Sample 3	0.743	

## 10.2. Semi-Quantitative test analysis

### 10.2.1. Antibodytitre calculation

In addition to the calculation of the antibody content in mIU/ml by using a defined standard, a semi-quantitative analysis can be carried out by the determination of the antibody titer. For this the unknown sample is used in various dilutions in the Mediagnost E10 ELISA. Depending on the expected level of antibodies, the sample is diluted more or less. For these dilutions the Dilution Buffer including 10 % negative serum (VPN) is used. To determine the titer of the (cut-off) value the following formula is used

$$cut - off = \frac{\text{Extinction Negative Control} + \text{Extinction Positive Control}}{2}$$

The serum dilution with an absorbance directly below the cut-off value represents the antibody titer.


### 10.2.2. Analysis with Mediagnost standards

A more accurate determination of serum titer can be carried out in semi- quantitative test method. **For** this purpose, the test kit includes 3 standards, which contain a defined amount of anti-HAV antibodies

These are ready to use and can be tested directly in the assay.

These samples are used to estimate the anti-HAV titer based on a calibration curve. The absorbance values of each serum standard (STD1 - STD3) and the negative control (NK) as standard with the concentration 0 mIU/mL are given on the y-axis against the concentration (mIU/mL) of the anti-HAV antibody on the x-axis and the regression line is placed through the points.

For computational evaluation various regression methods are useful. The method with the best curve fit (at least R2 > 0.9) should be selected. Ideally, the calculation of the antibody content of the sample is performed by means of an evaluation program. An exemplary standard curve determined by 4-parameter logistics is shown in Figure 1.

**Figure 1:** Exemplary standard curve with negative control as standard with the concentration of 0 mIU/mL.

In the table 2, exemplary results of anti-HAV Mediagnost antibody ELISA are shown.

**Table 2:** An example of an extinction with the semi-quantitative method when using the Mediagnost standards including negative control as standard with 0 mIU/mL. The sample dilution >1:10 was carried out in dilution buffer with 10% negative serum (VPN). The calculation of the antibody content of the sample was performed using an evaluation program.

		Blank corrected OD <sub>450-620</sub>	Anti-HAV Titer in the diluted Probe (mIU/mL)	Anti-HAV Titer in the undiluted sample (mIU/mL)
Positive Control		0.01	-	-
Negative Control		1.631	0	-
Standard 4				
Standard 1		0.083	50	-
Standard 2		0.190	30	-
Standard 3		0.735	15	-
Sample	Dilution			
Serum 1	1:10	0.089	48.14	481
Serum 2	1:10	0.294	22.18	222
Serum 3	1:10	1.264	4.17	42
Serum 4	1:50	0.598	12.36	618
Serum 5	1:50	0.121	39.97	1998

### 10.2.3. Analysis with reference preparation

In addition it is also possible to prepare standards in the desired concentrations from reference material to quantify the anti-HAV antibody.

The reference material should be diluted with the Dilution Buffer with 10% anti-HAV antibody negative serum (VPN). An exemplary result is shown in Table 3.

The cut off value for the evaluation of the results shown here is 0.724 OD450 units  $((0.02 + 1.427)/2)$ . The dilution of the sample, which absorbance is directly below the cut off value, is chosen for the calculation of the antibody content.

In serum 1 the absorbance of the 1:400 diluted sample is immediately below the cut-off. The anti-HAV titre in the diluted sample is 8.8 mIU/mL. The anti-HAV titer of the serum is therefore  $8.8 \times 400 = 3\,520$  mIU/mL.

In serum 2, the absorbance of 1:1600 diluted sample is immediately below the cut-off. The anti-HAV titre in the diluted sample is 8.5 mIU/mL. The anti-HAV titer in the undiluted serum is therefore  $8.5 \times 1600 = 13\,600$  mIU/mL.

**Table 3:** Exemplary extinction measurement by the semi-quantitative method using the reference preparation (NIBSC 97/646): The reference preparation and the unknown samples were diluted in Dilution Buffer with 10 % negative serum (VPN) to the concentration as given in the table and used in the assay. The calculation of the antibody content of the sample was performed using an evaluation program.

	Extinction-Meanvalue	Anti-HAV Titer in the diluted samples (mIU/mL)
Positive Control:	0.020	-
Negative Control	1.427	-
Standard 1*	0.068	50
Standard 2*	0.157	30
Standard 3*	0.334	15
Standard 4*	0.542	10
Standard 5*	0.838	5
Standard 6*	0.953	3
Standard 7*	1.245	1
<b>Serum 1:</b>		
1:100 diluted	0.110	
1:200 diluted	0.262	
1:400 diluted	0.571	$8.8 \times 400 = 3520$
1:800 diluted	0.761	
1:1600 diluted	1.048	
<b>Serum 2:</b>		
1:100 diluted	0.013	
1:200 diluted	0.036	
1:400 diluted	0.088	
1:800 diluted	0.228	
1:1600 diluted	0.532	$8.5 \times 1600 = 13600$

\*Not included in the test kit

## 11. LIMITATION OF PROCEDURE

The influence of the heterophilic antibodies, rheuma factors and anti-species antibodies is reduced by the assay design, but can not be completely excluded. The interference of various physiological and pharmaceutical substances was tested for the concentrations indicated. Higher concentrations or other substances may interfere with the measurement.

## 12. PERFORMANCE CHARACTERISTICS

### 12.1. Assay Calibration

The assay is calibrated to the International Standard WHO 97/646 material. The traceability of the measured results is shown in Figure 2. Exemplary results for the International Standard NIBSC 97/646 are shown for 30 and 10 mIU/mL from 1997 to 2014.


Figure 2: Traceability of Mediagnost anti HAV Antibody ELISA, E10

### 12.2. Sensitivity

The sensitivity of the test system was evaluated by the International standard NIBSC 97/646. The standard material was diluted and the dilution was applied as sample in the Mediagnost anti-HAV antibody ELISA, E10. The measured signals ( $OD_{450}$ ) are shown in table 4.

Table 4: Analytical sensitivity. Shown are the measured signals ( $OD_{450}$ ).

NIBSC [mIU/mL]	0	2	5	10	20
Replicates-	1.577	1.327	1.136	0.925	0.564
	1.495	1.337	1.159	0.83	0.475
	1.587	1.297	1.108	0.898	0.557
	1.329	1.373	1.142	0,79	0,51
	1.562	1.071	1.111	0,819	0,505
	1.278	1.219	1.083	0,854	0,398
		1.225	1.093	0,812	0,464
		1.331	1.199	0,865	0,45
Mean	1.47	1.27	1.13	0,85	0,49
SD	0.13	0.10	0.04	0,05	0,06
CV [%]	9	8	3	5	11

The recalculated analytical sensitivity resulting from the 2fold standard deviation of the negative control (n=6) is 1.10 mIU/mL. In an additional experiment the theoretically calculated analytical sensitivity was proven by 1 mIU/mL diluted NIBSC 97/646. Here the negative control showed a signal of 2.5 (SD 0.18) and the diluted international standard had a signal of 2.12 (SD 0.2).

Furthermore, through the measurement of 801 human serum samples (of which 287 were positive and 514 negative for anti-HAV antibodies) can be shown that the Mediagnost E10 ELISA distinguishes very well between positive and negative samples, the sensitivity was 98.26 %.

### 12.3. Specificity

Assessment of the specificity of the E10 ELISA was done with 25 samples negative for anti-HAV-IgG antibody positive for IgG Antibodies against different viruses (Table 5).

**Table 5: Specificity of the Mediagnost anti-HAV antibody test**

Sample No	Tested Positive for	Antibody	E10 result
1	anti-HBc, anti-HBs	IgG	negative
2	anti-HBs, anti-HCV	IgG	negative
3	anti-HCV	IgG	negative
4	anti-HEV	IgG	negative
5	anti-HCV; anti-HEV	IgG	negative
6	anti-HCV	IgG	negative
7	anti-HCV	IgG	negative
8	anti-HBe; anti-HBc	IgG; IgG & IgM	negative
9	anti-HBs	IgG	negative
10	anti-HBs	IgG	negative
11	anti-HBs	IgG	negative
12	anti-HBs, anti-HBc	IgG; IgG & IgM	negative
13	anti-HCV	IgG	negative
14	anti-HEV	IgG	negative
15	anti-HBs	IgG	negative
16	anti-HBs	IgG	negative
17	anti-HBs	IgG	negative
18	anti-HBs	IgG	negative
19	anti-HBe; anti-HBc	IgG; IgG	negative
20	HBsAG; anti-HBs	IgG; IgG & IgM	negative
21	anti-HEV	IgG	negative
22	anti-HBs	IgG	negative
23	anti-HBs; anti-HBe; anti-HEV	IgG; IgG; IgG	negative
24	HBsAG, anti-HBe; anti-HBc;	IgG; IgG; IgG & IgM	negative
25	anti-HBc, anti-HBs	IgG; IgG	negative

### 12.4. Precision Data

Serum samples were measured in at least 13 independent assays and the variation of the measurement result was calculated (Table 6).

**Table 6:** Inter-Assay Variation. SD= Standard deviation, CV% = Coefficient of Variation in %, Number = Number of the independent determinations.

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6	Sample 7	Sample 8	Sample 9	Sample 10
Mean [mIU/mL]	19303	3077	297	378	46	1210	45	27	9	5
SD	2557	378	36	54	13	156	5	3	2	1
CV%	13	12	12	14	29	13	11	13	17	19
Number	20	21	16	19	19	18	20	21	20	13

## 12.5. Linearity

The linearity of the assay was exemplary by the dilution of two human sera with high anti-HAV antibody titer. The dilution of the samples was carried out until 1:10 in Dilution Buffer (VP) and for all dilutions > 1:10 in Dilution Buffer + 10 % anti -HAV antibody negative serum (VPN). In Figure 3 the dilutions of 1: 2.5 to 1:80 for serum 1 and 1:5 to 1:320 for serum 2 are shown. By means of linear regression can be shown, that the linearity of the test system is very high.


Figure 3 Linearity of the sample dilution.

### 13. LITERATUR / LITERATURE

- 1) Flehmig, B.: Hepatitis A. Baillière`s Clinical Gastroenterology, Vol.4 No. 3, p. 707, 1990.
- 2) Ambrosch, F. et al.: Comparison of HAV-Antibodies Induced by Vaccination, Passive Immunization, and Natural Infection. *Viral Hepatitis and Liver Disease*. S. 98 (1990). (Editors: Hollinger, Lemon, Margolis, erschienen bei Williams and Wilkins, Baltimore; ISBN 0-683-04120-7).
- 3) Delem, A. et al.: Characterization of the immune response of volunteers vaccinated with a killed vaccine against hepatitis A. *Vaccine*, Vol. 11, Issue 4, S. 479 (1993).
- 4) Heinricy, U. et al.: Schedule-dependent Immune Response to Hepatitis A Vaccination. *Viral Hepatitis and Liver Disease*. S. 108 (1990). (Editors: Hollinger, Lemon, Margolis, erschienen bei Williams and Wilkins, Baltimore; ISBN 0-683-04120-7).
- 5) Jilg, W.: Hepatitis-A-Impfung - Anmerkungen zur Diagnostik. *Klin. Lab.* 38, S.697 (1992).
- 6) Van Damme, P.; Cramm, M.: Hepatitis A vaccination for health care workers. *British Medical Journal*, Vol. 306, S1615, 1993).
- 7) Turner, P. C. et al: Screening Before Hepatitis-A Vaccination. *Lancet* 340/8828, S.1160 (1992).
- 8) Germanaud, J.; Causse, X.: Stratègie de la vaccination contre l`hepatite A. *La Presse Médicale*, 22, n° 21, S 1104, (1993).
- 9) Brede, H. D.: Hepatitis-A-Gefahr für Hong-Kong Touristen, *BIOforum* 6, S.222 (1993).
- 10) Euteneuer, B.: Hepatitis A und Immunprophylaxe. *LABOR AKTUELL* 5.93, Herausg. Bioscientia Ingelheim.
- 11) Karpinski, K. F. et al.: Statistical considerations in the quantitation of serum immunoglobulin levels using the enzyme-linked immunosorbent assay. *Journal of Immunological Methods*, 103, S.189 (1987).
- 12) Bölke, E.; Flehmig B.: New epidemiological patterns of Hepatitis A and B in Germany. *Zbl.Hyg.* 196, 511-514 (1995)

#### 14. International Test description for QUALITATIVE Assay

<b>KK</b>	<b>CONJ</b>	1:100 <b>DILU</b> <b>BUF</b> VP
<b>WP</b>	<b>WASHBUF</b> <b>20x</b>	1:20 <b>DILU</b> A. dest.
<b>SPE</b>		1:10 <b>DILU</b> <b>BUF</b> VP
<b>°C</b> 20-25 °C		

-	-	A1/A2
100 µL	<b>Control+</b> PK	B1/B2
100 µL	<b>Control-</b> NK	C1/C2
100 µL	<b>SPE</b>	
<b>TAPE</b>		
 2 h <b>°C</b> 37°C		
50 µL	<b>CONJ</b> KK	B1/B2 → <b>End</b>
<b>TAPE</b>		
 1 h <b>°C</b> 37°C		
3x 300 µL	3x <b>WASHBUF</b> WP	
100 µL	<b>SUBST</b> <b>TMB</b> S	
 0.5 h <b>°C</b> 20-25°C 		
100 µL	<b>H<sub>2</sub>SO<sub>4</sub></b> SL	
<b>MEASURE</b>		