# Anti- Pseudomonas aeruginosa IgG ELISA

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Exemplary version, do not use to perform assays

Daines et al.	J Cystic Fibrosis	2014	AP/ ELA/ Exo A Mediagnost	68	100	66	56
---------------	-------------------	------	---------------------------	----	-----	----	----

## ENGLISH INSTRUCTIONS FOR USE

### TABLE OF CONTENTS

Instructions for use.....	18
1 INTENDED USE.....	18
2 INTRODUCTION.....	18
3 PRINCIPLE.....	19
4 WARNINGS AND PRECAUTIONS.....	20
5 SAMPLES.....	21
6 MATERIALS.....	22
7 TECHNICAL NOTES.....	23
8 ASSAY PROCEDURE.....	24
9 QUALITY CONTROL.....	24
10 EVALUATION OF RESULTS.....	25
11 LIMITATION OF PROCEDURE.....	27
12 PERFORMANCE CHARACTERISTICS.....	27
13 LITERATUR / LITERATURE.....	30
14 INTERNATIONALE ASSAY DESCRIPTION.....	31

## Instructions for use

<b>Anti-Pseudomonas aeruginosa IgG ELISA E15</b>	<b>3x96 Determinations</b>	
CE	DE/CA40/00809/7	
Principle of the test	Enzyme-linked Immunoassay	
Duration (incubation period)	4.5 h	
Antigen	Alkaline Protease, Elastase, Exotoxin A	
Buffer	Ready-for-use/ 20x concentrate	
Positive and Negative Control	Human Serum	
Sample	Human Serum/ EDTA-Plasma	
Sample Volume	10 µL	
Sample Dilution	1:1000	
Intra- / Inter-Assay Variance AP/ELA/Exo A > 1250 [1/Titre]	ø < 20%	
Result evaluation	<b>Cut-Off Values</b>	<b>Test Result</b>
	< 1:500	Negative
	500 - 1250	Borderline
	> 1250	Positive
Diagnostic sensitivity [%] <sup>1</sup>	ø 84 (Range 66 – 96)	
Diagnostic specificity [%] <sup>1</sup>	ø 76 (Range 48 – 96)	
<sup>1</sup> Mean value from	Kappler et al., Thorax 2006 Tramper-Stranders et al., Thorax 2006 Ratjen et al. Pediatric Pulmonol 2007 Dogru et al. Turkish Journal of Pediatrics 2013 Anstead et al J Cystic Fibrosis 2012 Daines et al J Cystic Fibrosis 2014	

### 1 INTENDED USE

The Mediagnost Anti-Pseudomonas aeruginosa IgG ELISA, E15, is an enzyme immunoassay for qualitative and quantitative detection of IgG-Antibodies in human serum against the extracellular proteins: Alkaline Protease (AP), Elastase (ELA) and Exotoxin A (Exo A) of Pseudomonas aeruginosa.

### 2 INTRODUCTION

*Pseudomonas aeruginosa*, a Gram-negative bacterium ubiquitously distributed in the moist environment, causes about 10% of all nosocomial infections. This opportunistic pathogen leads to acute and chronic types of infection within various organs of susceptible patient groups. Chronic pulmonary infection in patients with cystic fibrosis (CF) is very frequent, may start early in infancy and determines the life expectancy of these patients. *P. aeruginosa* infection provokes a production of antibodies to a large number of *P. aeruginosa* antigens in CF patients. **Mediagnost's** sensitive antibody detection system is an aid to discriminate between infected and uninfected patient groups.

Thus, the test is a supplement of microbiological diagnostics, which can be influenced by sampling technique (e.g. throat swab).

Depending on the *Pseudomonas aeruginosa* species and the patient's immune reaction, antibodies can be detected against a single, two or even all three antigens simultaneously. A patient is regarded as sero-positive when the serum is positive for one or more of the antigens.

## 1 PRINCIPLE

The **Mediagnost Anti-*Pseudomonas aeruginosa* IgG ELISA, E15**, is a sandwich enzyme immunoassay. Serum or plasma samples are diluted and added to the wells of a microtiter plate, which have been previously coated with the *Pseudomonas aeruginosa* antigens alkaline protease, elastase or exotoxin A. Specific antibodies in the sample bind to the respective antigens during an incubation of 2 h at 37°C. After washing, the conjugate (anti-human IgG peroxidase-labelled immunoglobulin) is added and incubated again 2 h at 37°C. After a final washing step, substrate is added and incubated for 30 min at room temperature. The color forming reaction is terminated on addition of stop solution accompanied by a color change from blue to yellow. The absorbance of the coloured reaction product is measured on a microtiter plate reader. The colour intensity of the reaction corresponds to the concentration of antibodies in the sample.

## 2 WARNINGS AND PRECAUTIONS

### For In Vitro Diagnostic Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro use and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use expired or obviously damaged or microbial contaminated or spilled material.

**Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.**

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations. Material Safety Data Sheet is available on request.

### Human Serum

Following components contain human serum: **PC1, PC2, PC3, NC, CTR**

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

### Reagents DET, DIL, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (< 0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

### Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (< 0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

### Stop Solution STP

The Stop solution contains 0.2 M sulphur acid (H<sub>2</sub>SO<sub>4</sub>)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

### 2.1 General first aid procedures:

**Skin contact:** Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

**Eye contact:** In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

**Ingestion:** After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

### **3 SAMPLES**

#### **3.1 Sample type**

Serum and EDTA-Plasma samples

#### **3.2 Specimen collection**

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

#### **3.3 Required sample volume: 10 µL**

#### **3.4 Sample stability**

In firmly closed sample vials:

- Storage at 20-25°C: at least 2 days
- Storage at -20°C: at least 2 years
- Freeze-thaw cycles max. 3

After the separation of coagulated and / or corpuscular components, it is recommended to store the samples as soon as possible at -20°C and to avoid freeze-thaw cycles.

#### **3.5 Sample dilution**

Samples should be diluted in PE/PP-tubes. For larger sample numbers usage of a multi-stepper is recommended. If the signal of the sample exceeds the signal of the positive control (PC), the measurement must be repeated for an exact titer result with a higher diluted sample e.g. in three dilutions (1:1000; 1:10000 and 1:100000).

Recommended dilution: **1:1000** with Dilution Buffer **DIL**.

#### **Exemplary Dilution Protocol**

##### **1:10 Pre-Dilution**

**90 µL** Dilution Buffer **DIL** + **10 µL** Serum- oder EDTA-Plasma

The pre-diluted sample should then be diluted as follows:

##### **1:1000 Dilution**

Add **10 µL** of the 1:10 pre-diluted sample to **990 µL** Dilution Buffer **DIL**

##### **1:10000 Dilution**

Add **100 µL** of the 1:1000 dilution to **900 µL** Dilution Buffer **DIL**

##### **1:100000 Dilution**

Add **100 µL** of the 1:10000 dilution to **900 µL** Dilution Buffer **DIL**

After mixing use **100 µL** each of the needed sample dilution per well within 1 hour (Pipetting control = orange color of the solution in the microtiter plate well).

## 4 MATERIALS

### 4.1 Materials provided

The reagents listed below are sufficient for 3x96 tests. The microtiter plates are divided up in 12 strips with separately breakable wells.

AP	<b>Microtiter plate 1</b> , ready for use: coated with alkaline Protease, labelled red.	<b>(8x12) Wells</b>
ELA	<b>Microtiter plate 2</b> , ready for use: coated with Elastase, labelled blue.	<b>(8x12) Wells</b>
Exo A	<b>Microtiter plate 3</b> , ready for use: coated with Exotoxin A, labelled green.	<b>(8x12) Wells</b>
PC1	<b>Positive Control 1</b> , ready for use, calibrated control serum for alkaline Protease, AP, <b>red labelled</b> , shows a titer of 1:2500.	<b>1 x 1.5 mL</b>
PC2	<b>Positive Control 2</b> , ready for use, calibrated control serum for Elastase, ELA, <b>blue labelled</b> , shows a titer of 1:2500.	<b>1 x 1.5 mL</b>
PC3	<b>Positive Control 3</b> , ready for use, calibrated control serum for Exotoxin A, Exo A, <b>green labelled</b> , shows a titer of 1:2500.	<b>1 x 1.5 mL</b>
NC	<b>Negative Control</b> , ready for use, contains human serum, <b>not reactive</b> for Pseudomonas aeruginosa antigens.	<b>2 x 1.5 mL</b>
CTR	<b>Control</b> , ready for use, contains human serum and has to be determined as <b>borderline</b> for Pseudomonas aeruginosa antigens in the assay.	<b>2 x 1.5 mL</b>
DET	<b>Antibody-HRP-Conjugate, 100-fold concentrated</b> , contains horseradish peroxidase (HRP)-labelled anti-human IgG.	<b>2 x 250 µl</b>
DIL	<b>Dilution Buffer</b> , read for use	<b>3 x 50 mL</b>
WB	<b>Washing Buffer</b> , 20-fold <b>concentrated</b> solution.	<b>1 x 120 mL</b>
S	<b>Substrate</b> for horseradish peroxidase (HRP), ready for use, stabilised Tetramethylbencidine (TMB).	<b>1 x 33 mL</b>
STP	<b>Stop Solution</b> , ready for use, 0.2 M sulphuric acid.	<b>1 x 33 mL</b>
-	<b>Sealing Tape</b> , for covering the <b>microtiter plate</b> .	<b>6 x</b>
	<b>Instructions for use</b>	<b>1 x</b>
-	<b>Quality Certificate</b>	<b>1 x</b>

### 4.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, **2400 mL**.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm
- Incubator **37°C**

## 5 TECHNICAL NOTES

### Storage Conditions

Store the kit at **2-8°C** until its expiry date.

### Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided.

### Preparation of reagents

Bring all reagents to room temperature (**20°C - 25°C**) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Dilution

Use the **Dilution Buffer DIL** for the **1:100** dilution of **Antibody-HRP-Conjugate DET**. Please dilute only according to daily requirements.

The required volume of Washing Buffer WB is prepared by 1:20 dilution of the provided 20-fold concentrate with Aqua dest. The **1:20** diluted Washing Buffer **WB** is stable **4 weeks at 2-8°C**. Please dilute only according to daily requirements.

### Assay Procedure

When performing the assay, Controls **CTR**, **PC**, **NC** and the samples should be pipetted as fast as possible (e.g., 15 minutes). To avoid distortions due to differences in incubation times the diluted Antibody-HRP-Conjugate **DET** and the **Substrate S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution STP** should be added to the plate in the same order as the Substrate **S**.

### Incubation

**Sample Incubation: 2 h at 37°C, Conjugate Incubation 2 h at 37°C, Substrate Incubation: 30 min at room temperature 20°C - 25°C.** The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive: storage and incubation in the dark.

**Incubation at room temperature means: Incubation at 20°C - 25°C.**

### Washing

Proper washing is of importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

When using an automatic microtiter plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

## 6 ASSAY PROCEDURE

Preparation of Reagents		Dilution		
DET	Antibody-HRP-Conjugate Concentrate (100x)	1:100 in Dilution Buffer <b>DIL</b> → <b>DET 1:100</b>		
WB	Washing Buffer Concentrate (20x)	in Aqua dest. (z.B. 100 mL <b>WB</b> mit <b>A.dest</b> auf 2000 mL auffüllen) → <b>WB 1:20</b>		
Dilute samples <b>1:1000</b> with Dilution Buffer <b>DIL</b> (qualitative test). For quantitative antibody assays it is recommended to perform serial dilutions of <b>1:1000</b> , <b>1:10000</b> and <b>1:100000</b> with Dilution Buffer <b>DIL</b> .				
Before the assay procedure bring all reagents to room temperature (20°C -25°C).				
<b>Assay procedure in double determination</b>				
Pipette	Reagents	Plate 1 AP	Plate 2 ELA	Plate 3 Exo A
100 µL	Negative Control <b>NC</b>	A1/A2	A1/A2	A1/A2
100 µL	Positive Control <b>PC1</b>	B1/B2	--	--
100 µL	Positive Control <b>PC2</b>	--	B1/B2	--
100 µL	Positive Control <b>PC3</b>	--	--	B1/B2
100 µL	Control <b>CTR</b>	C1/C2	C1/C2	C1/C2
100 µL	Sample Dilution <b>SPE</b>	Pipette in the rest of the wells according to requirements.		
Seal the wells of the 3 microtiter plates with the sealing tape				
Incubation: <b>2 h at 37°C</b>				
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL each Washing Buffer <b>WB 1:20/ well</b>			in each well
100 µL	Diluted Antibody-HRP-Conjugate <b>DET 1:100</b>			in each well
Seal the wells of the 3 microtiter plates with the sealing tape				
Incubation: <b>2 h at 37°C</b>				
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL each Washing Buffer <b>WB 1:20/ well</b>			in each well
100 µL	Substrate <b>S</b>			in each well
Incubation: <b>30 min in the dark at RT (20°C -25°C).</b>				
100 µL	Stop Solution <b>STP</b>			in each well
Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with <b>≥ 590 nm</b> as reference wavelength.				

## 7 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The Kit-control must be found within the acceptable range, which has been stated on the QC Certificate. The test results are only valid, if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws.

## 7.1 Quality criteria

For the evaluation of the assay it is preconditioned that the absorbance values of the Negative Controls **NC** should be **below 0.25**. The difference between the extinctions of Negative Controls **NC** and the respective Positive Control **PC** must be **at least 0.6**.

## 8 EVALUATION OF RESULTS

### 8.1 Calculation of results

Calculations should be performed for each antigen (each plate) separately! Calculate the average of all multiple values.

### 8.2 Qualitative calculation:

The negative control (Blank) average is subtracted from the controls and samples to obtain absolute values.

The **cut-off value** is 20% of the absorbance of the Positive Control PC.

1: 1000 diluted samples with extinction values less than the cut off value are classified as **negative**. The antibody titre is below 1:500.

1: 1000 diluted samples, having an absorbance in the range of  $\geq 20\%$  but  $< 50\%$  of positive control, cannot be clearly classified and must therefore be assessed as **borderline**.

1: 1000 diluted samples, which extinctions are  $\geq 50\%$  of PC extinction, are to be classified as **positive**.

### Exemplary qualitative calculation:

Negative Control NC	Extinction
1.Value	0.041
2.Value	0.056
<b>Mean</b>	<b>0.049</b>

Positive Control PC	Extinction
1.Value	1.120
2.Value	1.136
<b>Mean</b>	<b>1.128</b>

### Calculation

$$\text{PC} - \text{NC} \quad : 1.128 - 0.049 \quad = 1.079$$

$$\text{Cut-off (20\% of PC-NC)} \quad : 0.2 \times 1.079 \quad = 0.216$$

$$\text{Borderline (50\% of PC-NC)} \quad : 0.5 \times 1.079 \quad = 0.540$$

All samples with an extinction  $< 0.216$  are determined as negative for Anti-P. aeruginosa IgG. Samples with an extinction of  $> 0.216$  and  $< 0.540$  are judged as borderline and samples showing an extinction of  $> 0.540$  are determined as positive for the content of Anti-P. aeruginosa IgG.

### 8.3 Quantitative calculation

The evaluation is carried out graphically or via an appropriate evaluation program.

The extinction values of the negative NC and positive PC controls are plotted on the y-axis in a double-linear coordinate system against a titer factor. NC = 0 and PC = 2.5 on the x-axis. For the quantification of the sera values, a straight line is drawn through the NC and PC values and extended to a titer factor of 3.5.


The titre of the individual serum is determined by reading the titre factor of the measured extinction value through the NC-PC axis, which is multiplied by the serum dilution factor.

**Titre factors lower than 0.25 and higher than 3.5 (x-axis) are not taken into consideration in these calculations.**

### Exemplary quantitative calculation:

Negative Control NC	Extinction
1.Value	0.041
2.Value	0.056
<b>Mean</b>	<b>0.049</b>

Positive Control PC	Extinction
1.Value	1.220
2.Value	1.176
<b>Mean</b>	<b>1.198</b>

**Figure 1:** Exemplary graphic for the quantitative evaluation for serum B.

Sample	Dilution of Sample	Extinction (average)	Titre factor (see graphic)	Titre	Interpretation
Serum B	1:1000	2.82	> 3.5	n.a.	positive
	1:10000	1.45	3.05	1:(10000 x 3.05) => 1:30 500	positive
	1:100000	0.210	0.3	1:(100000 x 0.3) => 1:30000	positive

If an electronic analysis program is used, NC and PC have to be set as standards with the titer values of 0 and 2500. The program calculates the titer of the unknown sample. The dilution factor of the sample has to be taken into account.

#### 8.4 Interpretation of results

This ELISA Kit is used as a complementary test to a classical microbiological detection of infection with *Pseudomonas aeruginosa*. If a sample is tested positive for antibodies against one of the three *Pseudomonas aeruginosa* exotoxins, this is a sign of corresponding infection. The evaluation of the measured antibody titres has to be done according to the scheme below. In the literature there are several studies describing the diagnostic performance of the assay (see table 5). These studies often do not differentiate between borderline and positive and apply only one cut-of value (1:500) to differentiate between negative and positive.

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For evaluation of the test results, the diagnostic sensitivity and specificity of the test system must be taken into account (see Section 12.7).

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other

diagnostic investigations. It is recommended to consider the international and national guidelines for diagnosis and treatment.

## 9 LIMITATION OF PROCEDURE

Basically, the result of immunological test systems can be affected by various sample components such as medications or lipids. Their influence is reduced by the assay design, but cannot be excluded completely. The evaluation of the results must be made in the context of diagnostic sensitivity and specificity of the assay (table 5). In acute and chronic *P. aeruginosa* infected patients suffering from immunosuppression, assessment of an infection via the detection of antibodies may be inadequate.

## 10 PERFORMANCE CHARACTERISTICS

### 10.1 Analytical Sensitivity

The analytical sensitivity as a measure for the minimal amount of specific antibody detectable by this test system was determined by the signal variability of the negative control. Therefore, the recalculated antibody titre of the negative control was evaluated and the 3-fold standard deviation was used to determine the limit of detection. The results are shown in table 1. The limit of detection of antibody titers for AP, ELA, Exo A in 1:1000 diluted samples is < 1:300.

**Table 1 Analytical sensitivity / Limit of Detection** based on 3-fold standard deviation of the blank. Shown is the recalculated antibody titre of three independent experiments.

Dilution	AP	ELA	Exo A
-	0.105	0.268	0.181
1:100	10	27	18
1:1000	105	268	181

### 10.2 Analytical Specificity

The analytical specificity was assessed by the evaluation of the microbiological status of serologically and microbiologically *Pseudomonas*-negative cystic fibrosis patients. Most frequently infections with *Staphylococcus aureus*, *Candida albicans* and *Hemophilus influenza* were detected. These infections do not result in antibodies cross reacting with *Pseudomonas aeruginosa* antigens.

### 10.3 Precision

#### Intra-Assay Variance

Serum samples were diluted 1:1000 and measured 16-fold within one assay. The results are shown in table 2. Mean variance was 4.95, 4.49 and 8.41% for alkaline Protease, Elastase and Exotoxin A, respectively.

**Table 2 Intra-Assay Variance**

	Sample 1			Sample 2			Sample 3		
	AP	ELA	Exo A	AP	ELA	Exo A	AP	ELA	Exo A
Mean Titre	806	529	715	3858	5521	5673	1453	1530	1381
SD	35	43	91	173	135	127	88	44	142
VC [%]	4.3	8.14	12.71	4.47	2.45	2.25	6.08	2.86	10.29
Number[n]	16	16	16	16	14	16	16	16	16

#### Inter-Assay Variance


Samples were portioned and stored at -20°C for up to 12 years and measured irregularly during this time. Table 3 shows the summarized results for 3 different samples. The variability measured as coefficient of variation does not exceed 20% for any antigen or sample.

**Table 3 Inter-Assay Variance**

	Sample 1			Sample 2			Sample 3		
	AP	ELA	Exo A	AP	ELA	Exo A	AP	ELA	Exo A
<b>Mean Titre</b>	<b>1041</b>	<b>958</b>	<b>804</b>	<b>1743</b>	<b>1519</b>	<b>1239</b>	<b>2971</b>	<b>4071</b>	<b>4343</b>
<b>SD</b>	67	60	70	228	231	176	164	309	392
<b>CV [%]</b>	6.39	6.31	8.69	13.09	15.22	14.22	5.51	7.59	9.03
<b>Number [n]</b>	73	73	72	295	295	293	112	112	112

### 10.4 Linearity

Linearity of sample dilution was evaluated by dilution dependent signal decrease in three serum samples with high, medium and low antibody titres. In figure 2 the measured signals are shown in dependence of the sample dilution. Linear regression analysis revealed coefficients of correlation of >0.9 for all tested samples. The recalculated antibody titres show a good linearity, too.

**Figure 2 Linearity of sample dilution.** Samples with high, medium and low antibody titre were diluted (1:500 – 1:24000) depending on the expected signal intensity and measured results were analysed by linear regression analysis.

### 10.5 Interference

Interference was not investigated, no information regarding the influence of triglycerides, hemoglobin or bilirubin is available. Therefore lipaemic, hemolytic or icteric samples should not be used for this assay.

### 10.6 Traceability / Assay Calibration

Neither an international standard preparation, purified human Anti-Pseudomonas aeruginosa antibodies nor an established reference reference method for measurement of these antibodies are available. Thus, the assay calibration is based on manufacturer positive controls with arbitrary antibody titres of 2500 (Mediagnost [1/Titre]). Selected human serum samples with defined antibody content are used for assay production and quality control. This panel allows the traceability of measurement results.

### 10.7 Diagnostic Quality

For evaluation of diagnostic sensitivity (true positive) and specificity (true negative) a comparison of E15 titre results in serum samples and the microbiological detection of Pseudomonas aeruginosa growth in throat swab samples by classical culture techniques has been conducted.

For the evaluation of the test results regarding Pseudomonas aeruginosa infection only a distinction between positive and negative was made.

## Positive

At least one microbiological detection of *Pseudomonas aeruginosa* within 12 months before the serum sample was collected.

## Negative

No microbiological detection of *Pseudomonas aeruginosa* within 12 months before the serum sample was collected

**Table 4 Diagnostic Specificity and Sensitivity** Here the absolute numbers of positive and negative results are shown in relation to the microbiological gold standard for different cut-off values (AP: alkaline Protease, ELA: Elastase; Exo A: Exotoxin A; Total: all three antigens). Based on the absolute numbers sensitivity (% true positive) and specificity (% true negative) are calculated).

Cut- off	1:100				1:500				1:1250			
	AP	ELA	Exo A	Total	AP	ELA	Exo A	Total	AP	ELA	Exo A	Total
True Positive <b>TP</b>	39	40	43	45	27	31	27	40	19	23	19	32
False Positive <b>FP</b>	10	8	11	12	2	2	3	5	1	0	0	1
True Negative <b>TN</b>	16	18	15	14	24	24	23	21	25	26	26	25
False Negative <b>FN</b>	12	11	8	6	24	20	24	11	32	28	32	19
<b>Sensitivity [%]</b> TP/(TP+FN)	76	78	84	88	53	61	53	78	37	45	37	63
<b>Specificity [%]</b> TN/(TN+FP)	62	69	58	54	92	92	88	81	96	100	100	96

A population of 77 cystic fibrosis patients was investigated, 51 patients were microbiologically positive for *Pseudomonas aeruginosa* before the serum sample was collected and 26 negative. These results conform to the data published by Kappler et al in 2006. Here the authors used a cut-off Titre of 500 and calculated an overall sensitivity of 86% and a specificity of 96%. In detail it was shown that sensitivity was 64; 53 and 73% for AP, ELA and Exo A, respectively. Specificity for AP, ELA and Exo A was calculated as 99; 100 and 96%. In comparison to our results Kappler et al. detected higher sensitivity and specificity values which probably are the consequence of an increased sample number (183 vs 77).

The Mediagnost test system has been evaluated by different independent research groups regarding its diagnostic sensitivity and specificity. In table 5 the results are summarized. Mean diagnostic specificity is 76% (48 - 96) and mean diagnostic sensitivity is 84% (66–96).

**Table 5 Summary of literature results** describing diagnostic sensitivity and specificity in dependence of the applied cut-off titre value.

Autor	Journal	Jahr	Mediagnost Assay	Patient [n]	Cut-off Titre	Diagnostic Sensitivity [%]	Diagnostic Specificity [%]
Kappler et al.	Thorax	2006	AP/ ELA/ Exo A	183	500	86	96
Tramper-Stranders et al.	Thorax	2006	AP/ ELA/ Exo A	220	1250	66	96
				220	500	79	89
					AP > 0 ELA > 35 Exo A > 0	96	79
Ratjen et al.	Pediatric Pulmonol	2007	AP/ ELA/ Exo A	375	AP > 285 ELA > 300 Exo A > 1000	92.7	93.3
Dogru et al.	Turkish Journal of Pediatrics	2013	AP/ ELA/ Exo A	90	500	87.5	70
Anstead et al.	J Cystic Fibrosis	2013	AP/ ELA/ Exo A	304	100	78	48
Daines et al.	J Cystic Fibrosis	2014	AP/ ELA/ Exo A	68	100	66	56

## 11 LITERATUR / LITERATURE

### Literature

1. West S, Zeng L, Lee L, Kosorok M, Laxova A, Rock M, Splaingard M, Farrell : Respiratory infections with *Pseudomonas aeruginosa* in children with cystic fibrosis, early detection by serology and assessment of risk factors. JAMA, 2002, Vol 287, No.: 22: 2958-2967
2. Burns J, Gibson I, McNamara, Yim D, Emerson J, Rosenfeld M, Hiatt P, McCoy K, Castille R, Smith A, Ramsey B: Longitudinal Assesment of *Pseudomonas aeruginosa* in young children with cystic fibrosis. J Infect Dis 2001 ; 183 : 444-52
3. Döring G, Obernesser H-J, Botzenhart K: Extracellular toxins of *Pseudomonas aeruginosa*. III Radioimmunoassay for detection of alkaline protease. Zentralbl Bakteriol Mikrob. Hyg (A) 1982; 252:239-47
4. Döring G., Obernesser H-J: Botzenhart K., Flehmig B: Høiby N., Hofmann A.: Proteases of *Pseudomonas aeruginosa* in cystic fibrosis. J Infect Dis 1983; 147:744-50
5. Döring G.; Høiby N: Longitudinal study of immune response to *Pseudomonas aeruginosa* antigens in cystic fibrosis. Infect Immun 1983; 42:197-201
6. Obernesser H-J, Döring G: Extracellular toxins of *Pseudomonas aeruginosa*. IV. Radioimmunoassay for detection of elastase. Zentralbl Bakteriol Mikrobiol Hyg (A) 1982; 252:248-56
7. Brett M M, Ghoneim A T M, Littlewood J M: Serum antibodies to *Pseudomonas aeruginosa* in cystic fibrosis. Archives of Diseases in Childhood, 1986, 61, 1114-20
8. Brett M M, Ghoneim A T M, Littlewood J M; Losowsky M S: Development of enzyme linked immunosorbent assay (ELISA) to detect antibodies to *Pseudomonas aeruginosa* cell surface antigens in sera of patients with cystic fibrosis. J Clin Pathol 1986; 39: 1124-29
9. Brett M M, Ghoneim A T M, Littlewood J M: Serum IgG antibodies in patient with cystic fibrosis with early *Pseudomonas aeruginosa* infection: Archives of Disease in Childhood, 1987, 62, 357-361
10. Brett M M, Ghoneim A T M, Littlewood J M: Prediction and Diagnosis of Early *Pseudomonas aeruginosa* Infection in Cystic Fibrosis: a Follow-up Study: J of Clin Microb, Aug 1988 p 1565-70
11. Brett M M, Simmonds E J, Ghoneim A T M, Littlewood J M: The value of serum IgG titres against *Pseudomonas aeruginosa* in the management of early peudomonal infection in cystic fibrosis. Archives of Diseases in Childhood 1992; 67: 1086-1088
12. Danielsen L, Westh H, Balsev E, Rosdahl VT, Döring G: *Pseudomonas aeruginosa* exotoxin A antibodies in rapidly deteriorating chronic leg ulcers. The Lancet Jan 1996 Vol. 347 No.8996 p 265
13. Dasgupta MK; Lam J, Döring G, Harley FL, Zuberbuhler P, Lam K, Reichert A, Costerton JW, Dossetor JB: Prognostic implications of circulating immune complex and *Pseudomonas aeruginosa*-specific antibodies in cystic fibrosis. J Clin Lab. Immunol. 1987, 23, 25-30
14. Döring G, Dalhoff A, Vogerl O, Brunner H, Dröge U, Botzenhard K. In vivo activity of Proteases of *Pseudomonas aeruginosa* in a rat model. The J of Infectious Diseases, Apri 1994, 4/Vol 149 p 532-537
15. Döring G, Buhl V, Høiby N, Schiøth P O, Botzenhart K: Detection of *Pseudomonas aeruginosa* in immune complexes isolated from sputum of cystic fibrosis patients. Acta path microb. immunol. scand Sect. C, 92: 307-312, 1984
16. Döring G, Goldstein W, Röhl A, Schiøtz P O, Høiby N, Botzenhard K: Role of *Pseudomonas aeruginosa* Exoenzymes in Lung Infections of Patients with Cystic Fibrosis. Infection and Immunity, Sept 1985, p. 557-562
17. M Kappler 1 , A Kraxner , D Reinhardt , B Ganster , M Griese , T Lang: Diagnostic and Prognostic Value of Serum Antibodies Against *Pseudomonas Aeruginosa* in Cystic Fibrosis Thorax 61. 2006 Jan 31. 8, 684-688
18. Tramper-Stranders GA, van der Ent CK, Slieker MG, Terheggen-Lagro SW, Teding van Berkhout F, Kimpen JL, Wolfs TF. Diagnostic value of serological tests against *Pseudomonas aeruginosa* in a large cystic fibrosis population. Thorax. 2006 Aug; 61,(8):689-93.
19. Ratjen F1, Walter H, Haug M, Meisner C, Grasmann H, Döring G: Diagnostic value of serum antibodies in early *Pseudomonas aeruginosa* infection in cystic fibrosis patients. Pediatr Pulmonol. 2007 Mar;42(3):249-55.
20. Doğru D., Pekcan S., Yalçın E. Özçelik U., Kiper N., Gürcan N., Şener B.: The role of serum *Pseudomonas aeruginosa* antibodies in the diagnosis and follow-up of cystic fibrosis. The Turkish Journal of Pediatrics 2013; 55: 50-57
21. Anstead M., Heltshe S., Khan U., Barbieri J., Langkamp M., Döring G., Dharia S., Gibson R., Treggiari M., Lymp J., Rosenfeld M., B. *Pseudomonas aeruginosa* serology and risk for re-isolation in the EPIC trial. Journal of Cystic Fibrosis 2013, 12 ,2 , 147 - 153
22. Daines C., VanDeVante D., Khan U., Emerson J., Heltshe S., McNamara S., Anstead M., Langkamp M., Döring G., Ratjen F., Ramsey B., Gibson R., Morgan W., Rosenfeld M.: Serology as a diagnostic tool for predicting initial *Pseudomonas aeruginosa* acquisition in children with cystic. Journal of Cystic Fibrosis 2014, 13 , 5 , 542 - 549

# 1 INTERNATIONALE ASSAY DESCRIPTION

## International Test Description

DET	100x	1:100 DIL → DET 1:100		
WB	20x	1:20 A. dest. → WB 1:20		
-	SPE 1:1000 DIL ↔			
-	°C 20°C - 25°C ⌚ max. 1 h			
		MTP AP	MTP ELA	MTP Exo A
100 µL	NC	A1/A2	A1/A2	A1/A2
100 µL	PC1	B1/B2	--	--
100 µL	PC2	--	B1/B2	--
100 µL	PC3	--	--	B1/B2
100 µL	CTR	C1/C2	C1/C2	C1/C2
100 µL	SPE 1:1000 DIL			
TAPE				
⌚ 2 h °C 37°C				
3x 300 µL	3x WB 1:20			
100 µL	DET 1:100			
TAPE				
⌚ 2 h °C 37°C				
3x 300 µL	3x WB 1:20			
100 µL	S			
⌚ 0.5 h °C 20°C - 25 °C ☀				
STP				
MEASURE				

## 1 ASSAY PROCEDURE

Preparation of Reagents		Dilution		
<b>DET</b>	<b>Antibody-HRP-Conjugate Concentrate (100x)</b>	1:100 in Dilution Buffer <b>DIL</b> → <b>DET 1:100</b>		
<b>WB</b>	<b>Washing Buffer Concentrate (20x)</b>	in <b>Aqua dest.</b> (z.B. 100 mL <b>WB</b> mit <b>A.dest</b> auf 2000 mL auffüllen) → <b>WB 1:20</b>		
Dilute samples <b>1:1000</b> with Dilution Buffer <b>DIL</b> (qualitative test). For quantitative antibody assays it is recommended to perform serial dilutions of <b>1:1000</b> , <b>1:10000</b> and <b>1:100000</b> with Dilution Buffer <b>DIL</b> .				
Before the assay procedure bring all reagents to room temperature ( <b>20°C -25°C</b> ).				
<b>Assay procedure in double determination</b>				
Pipette	Reagents	Plate 1 AP	Plate 2 ELA	Plate 3 Exo A
100 µL	Negative Control <b>NC</b>	A1/A2	A1/A2	A1/A2
100 µL	Positive Control <b>PC1</b>	B1/B2	--	--
100 µL	Positive Control <b>PC2</b>	--	B1/B2	--
100 µL	Positive Control <b>PC3</b>	--	--	B1/B2
100 µL	Control <b>CTR</b>	C1/C2	C1/C2	C1/C2
100 µL	Sample Dilution <b>SPE</b>	Pipette in the rest of the wells according to requirements.		
Seal the wells of the 3 microtiter plates with the sealing tape				
Incubation: <b>2 h at 37°C</b>				
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL each Washing Buffer <b>WB 1:20</b>		in each well	
100 µL	Diluted <b>Antibody-HRP-Conjugate</b> <b>DET 1:100</b>		in each well	
Seal the wells of the 3 microtiter plates with the sealing tape				
Incubation: <b>2 h at 37°C</b>				
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL each Washing Buffer <b>WB 1:20/ well</b>		in each well	
100 µL	Substrate <b>S</b>		in each well	
Incubation: <b>30 min</b> in the dark at <b>RT (20°C -25°C)</b> .				
100 µL	Stop Solution <b>STP</b>		in each well	
Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with <b>≥ 590 nm</b> as reference wavelength.				

This assay procedure is a duplicate of the page 24.