# m/rIGF-I ELISA

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<b>m/rIGF-I ELISA, E25</b>	<b>96 Determinations</b>
Regulatory Status	For Research Use Only.
Principle of the test	Enzyme Immunoassay
Duration (incubation period)	2 h
Sample Buffer	Ready for use
Washing Buffer	20fold concentrated
Standard	Single Standards <b>0.5 - 18 ng/ml</b>
Calibration	The assay is calibrated against recombinant IGF-I
Assay Range	<b>0.315 ng/mL – 1800 ng/mL</b>
Control	2 Serum Controls, low and high concentrations, lyophilised
Sample	mouse / rat serum and plasma
Required sample volume	10 µL recommended, at least 5 µL
Sample dilution	1:100
Analytical Sensitivity	Ø 0.315 ng/mL
Intra- / Interassay Variance	Ø < 10 %

## 1 INTENDED USE


Measurement of IGF-I in mouse / rat serum and plasma.

## 2 INTRODUCTION

Beside different cell culture models and studies with human material, mice and rats are suitable model organisms for basic research and pre-clinical studies. Thus, we developed this test system as a tool for IGF-I measurements in mice and rat for usage in research and pre-clinical studies. Even if the comparability of mice and humans is limited, we offer some background information on the human IGF-I system in the following section:

Insulin-like growth factors (IGF) I and II play a pivotal role in regulating the proliferation and differentiation of many cell types (1-3). IGF-I is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 daltons (5). Its major regulators are growth hormone (GH) and nutrition (6). In contrast to many other peptide hormones, IGFs are avidly bound to specific binding proteins (IGFBP). The seven IGFBPs which are known at present (7,8,22) either bind IGF-I and IGF-II with similar affinities or show a preference for IGF-II (9,10).

A major problem of IGF-I measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-I/IGFBP-3 complexes during the assay incubation. Depending on the ratio IGF-I to IGFBP in the sample interference comes up (see example Figure 1)


**Figure 1** Interference of IGFBP in IGF-I measurements. Known concentrations of IGF-I were assayed in the presence of 0.5 ng (left) or 5 ng (right) hIGFBP-3 by a conventional (□) and by the IGFBP-blocked assay (\*).

To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation, except acidification or dilution in a specially designed buffer system, is not required before measurement.

### 3 PRINCIPLE

The Mediagnost ELISA for m/r IGF-I E25 is a so-called Sandwich-Assay. It utilizes two specific

and high affinity antibodies for this protein. The IGF-I in the sample binds to the immobilized first antibody on the microtiter plate, the biotinylated and Streptavidin-Peroxidase conjugated second specific anti-IGF-I-Antibody binds in turn to the immobilized IGF-I. In the closing substrate reaction, the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-I-level of the samples. In order to dissociate IGF-I from the IGFBPs, the samples must be diluted in an acidic buffer (Figure 2). The diluted samples are then pipetted into the wells, by this the pH-value will be neutralized. After neutralization of the samples, the excess IGF-II occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of resulting free IGF-I. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the extremely low cross-reactivity of the IGF-I antibody with IGF-II, the excess of IGF-II does not disturb the interaction with IGF-I. The test runs like a conventional ELISA using a Streptavidin-Peroxidase-Enzyme Conjugate.

**Figure 2** Principle of the IGFBP blocked IGF-I ELISA

## 4 WARNINGS AND PRECAUTIONS

### For In Vitro Use only. For Professional use only. For Research Use Only.

The Mediagnost kit is suitable only for research use only and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obviously damaged or microbial contaminated or spilled material.

**Caution: This kit contains material of human and/or animal origin. Therefore, all components and specimens should be treated as potentially infectious.**

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations. Material Safety Data Sheet is available on request.

#### Mouse / Rat Serum

Following components contain Mouse or Rat serum: Control Serum KS1 and KS2

#### Reagents A-E, AK, EK, VP, WP

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

#### Substrate Solution (S)

TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

#### Stopping Solution (SL)

The Stopping solution contains 0.2 M acid sulphur acid (H<sub>2</sub>SO<sub>4</sub>)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

### 4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

## 5 SAMPLES

### 5.1 Sample type

Serum samples as well as Heparin-, EDTA- and Citrat-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Influence of Heparin (30IE/mL), EDTA (6,8mM) and NaCitrat (0,015M) on the measurement of IGF-I has been investigated in recovery experiments. Buffer solution was enriched with recombinant IGF-I and the above-mentioned substances. No significant influence on the recovery of IGF-I was detected, on average the recovery of recombinant material in comparison to enriched PBS was 108%.

Cell culture medium is suitable as sample matrix after predilution of 1:2 with **Sample Buffer PP**.

### 5.2 Specimen collection

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

**Required sample volume:** recommended 10 µL, min. 5 µL

### 5.3 Sample stability

In firmly closed sample vials:

- Storage at 20-25°C: max. 2 days
- Storage at -20° C: max. 2 years
- Freeze-thaw cycles max. 2

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading.

Therefore, it is recommended to keep sample refrigerated or frozen as soon as possible after separation of coagulated and corpuscular blood components and to avoid more than 2 freeze-thaw cycles.

### 5.4 Sample dilution

- Dilution: **1:100** with **Sample Buffer PP**
- Pipette **990 µL Sample Buffer PP** in PE-/PP-Tube (application of a multi-stepper is recommended in larger series); add **10 µL sample** (dilution 1:100). After mixing use 50 µL of this dilution in the assay.
- Attention: serum and plasma samples must be diluted at least 1:10 in **Sample Buffer PP** in order to achieve sufficient acidification of the samples.
- Depending on the expected IGF-I values the samples can be diluted higher or lower in **Sample Buffer PP**.

## 6 MATERIALS

### 6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the standard curve.

<b>MTP</b>	<b>Microtiter plate</b> , ready for use, coated with hamster-anti-m/r-IGF-I-antibody. Wells are	<b>(8x12) wells</b>
<b>A-E</b>	<b>Standards</b> , lyophilized, (recombinant IGF-I), concentrations are given on vial labels and on quality certificate.	<b>5 x 1 mL</b>
<b>KS1</b>	<b>Control Serum 1</b> , lyophilised, (Mouse/ Rat Serum), concentration is given on quality certificate in ng/mL.	<b>1 x 500 µL</b>
<b>KS2</b>	<b>Control Serum 2</b> , lyophilised, (Mouse/ Rat Serum), concentration is given on quality certificate in ng/mL.	<b>1 x 500 µL</b>
<b>AK</b>	<b>Antibody Conjugate</b> , ready for use, contains goat biotinylated anti-m/rIGF-I antibody.	<b>1 x 7 mL</b>
<b>EK</b>	<b>Enzyme Conjugate</b> , ready for use, contains HRP (Horseradish-Peroxidase)-labelled Streptavidin.	<b>1 x 12 mL</b>
<b>PP</b>	<b>Sample Buffer</b> , ready for use	<b>1 x 125 mL</b>
<b>WP</b>	<b>Washing Buffer</b> , 20-fold concentrated solution	<b>1 x 50 mL</b>
<b>S</b>	<b>Substrate</b> , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine.	<b>1 x 12 mL</b>
<b>SL</b>	<b>Stopping Solution</b> , ready for use, 0.2 M sulphuric acid.	<b>1 x 12 mL</b>
-	<b>Sealing Tape</b> , for covering the <b>microtiter plate</b> .	<b>2 x</b>
	<b>Instructions for use</b>	<b>1 x</b>
--	<b>Quality Certificate</b>	<b>1 x</b>

### 6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WP (A. dest.)**, **950 mL**.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and <sup>3</sup>590 nm

## 7 TECHNICAL NOTES

### Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at –20 °C after reconstitution. Avoid repeated thawing and freezing.

### Storage Life

Store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** standards **A-E** and Control Sera **KS1** and **KS2** must be stored at –20°C (max. 2 Months). **Attention:** Standards should be thawed only once- where required please store aliquoted in adequate volumens. For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). One of these freeze-thaw cycles showed no significant effect on the test

The 1:20 diluted Washing Buffer **WP** is 4 weeks stable at 2-8°C

### Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Reconstitution

The Standards **A – E** and Control **KS1** and **KS2** are reconstituted with the Sample Buffer **PP**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

### Dilution

After reconstitution dilute the Control **KS1** and **KS2** with the Sample Buffer **PP** in the same ratio as the sample.

The required volume of Washing Buffer **WP** is prepared by 1:20 dilution of the provided 20 fold concentrate with Aqua dest.

### Assay Procedure

When performing the assay, Blank, Standards **A-E**, Controls **KS1** and **KS2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **AK** and the Enzyme Conjugate **EK** as well as the succeeding Substrate Solution **S** should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution **SL** should be added to the plate in the same order as Substrate Solution **S**. All determinations (Blank, Standards **A-E**, Control **KS1** and **KS2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

### Incubation

**Incubation at room temperature means: Incubation at 20 - 25°C.** The Substrate Solution **S**, stabilised Tetramethylbencidine, is photosensitive—store and incubation in the dark.

### Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending approx. 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

### Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WP** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

**Manual washing** should be used. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. Decant contents into a biohazard bin, then blot plate on absorbent tissue. Wash the plate by adding 300 µL Washing Buffer **WP**/well, then decant and blot on absorbent tissue. Repeat this step 4 more times for total of 5 washes.

## 8 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
A-E	Standards	in 1 mL Sample Buffer PP	-
KS1	Control Serum 1	in 500 µL Sample Buffer PP	1:100 with Sample Buffer PP
KS2	Control Serum 2	in 500 µL Sample Buffer PP	1:100 with Sample Buffer PP
WP	Washing Buffer	-	1:20 with Aqua dest.
<p>Dilute samples 1:100 in Sample Buffer PP, mix immediately, incubate at least 15 min, max 2 h, incubate max. 2h.            Use 50 µl for each well in the assay.            Before assay procedure bring all reagents to room temperature 20-25°C.</p>			
<b>Assay procedure in double determination</b>			
Pipette	Reagents		Position
50 µL	Antibody Conjugate AK		in <u>all</u> wells used
50 µL	Sample Buffer PP (Blank)		A1/A2
50 µL	Standard A (0.5 ng/mL)		B1/B2
50 µL	Standard B (2.5 ng/mL)		C1/C2
50 µL	Standard C (6 ng/mL)		D1/D2
50 µL	Standard D (12 ng/mL)		E1/E2
50 µL	Standard E (18 ng/mL)		F1/F2
50 µL	Control Serum KS1	(1:100 diluted)	G1/G2
50 µL	Control Serum KS2	(1:100 diluted)	H1/H2
50 µL	Sample	(1:100 diluted)	in the rest of the wells according the requirements
Cover the wells with the sealing tape.			
<b>Sample-Incubation: 1 h at 20-25°C, 350 rpm</b>			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well		In each well
100 µL	Enzyme Conjugate EK		In each well
Cover the wells with the sealing tape.			
<b>Incubation: 30 Minutes at 20-25°C, 350 rpm</b>			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well		In each well
100 µL	Substrate Solution S		In each well
<b>Incubation: 30 Minutes in the Dark at 20-25°C</b>			
100 µL	Stopping Solution SL		In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			

## 9 EVALUATION OF RESULTS

For the evaluation of the assay it is required that the absorbance values of the should be below 0.25, and the absorbance of standard E should be above 1.00. Samples, which yield higher absorbance values than **Standard E**, should be re-tested with a higher dilution.

### 9.1 Establishing of the standard curve

Standards are provided in the following concentrations:


Standard	A	B	C	D	E
ng/mL	0.5	2.5	6	12	18

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples and standards
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Calculation of the standard curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. A four parametric logistic (4-PL) curve fit should be used for recalculation of IGF-I concentrations.
- 5) The IGF-I concentration of the diluted sample or the diluted control sera KS1&2 in ng/ml (or µg/ml according the chosen unit for the standards) is calculated in this way, the IGF-I concentration of the **undiluted sample** and of KS1 & KS2 is calculated **by multiplication** with the respective dilution factor.

### 9.2 Example of a typical standard curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

	Blank	A	B	C	D	E
ng/mL	0	0.5	2.5	6	12	18
OD (450-620 nm)	0.07	0.114	0.614	1.283	1.690	1.923

**Figure 3** Exemplary standard curve

The exemplary shown standard curve in Figure 3 **cannot** be used for calculation of your test results. You have to establish a standard curve for each test you conduct!

### 9.3 Exemplary calculation of IGF-I concentrations

Sample dilution: 1:100

Measured extinction of your sample: 1.3525

Measured extinction of the blank 0.07

Your measurement program will calculate the IGF-I concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit.

In this exemplary case the following equation is solved by the program to calculate the IGF-I concentration in the sample:

$$1.283 = 0.0004x^3 - 0.0184x^2 + 0.3029x \quad 6.107 = x$$

If the dilution factor (1:100) is taken into account the IGF-I concentration of the undiluted sample is  $6.107 \times 100 = 610.7 \text{ ng/mL}$

## 10 PERFORMANCE CHARACTERISTICS

### 10.1 Sensitivity

The sensitivity was determined by measuring the blank value and calculating the theoretical concentration of the blank value +2 SA. The analytical sensitivity of the Mediagnost E25 is as an average  $\emptyset$  0.315 ng/mL (range 0,262-0,405).

### 10.2 Precision

#### Intra-Assay Variance

Several samples were measured 16 in the same assay. Exemplary results are shown in Table 1. On average, the variation coefficient was  $\emptyset < 10\%$ .

**Table 1 The Intra-assay variability.** The IGF-I concentrations were determined in a duplicate assay and variability was calculated as the coefficient of variation (CV).

	Number of Determinations	Mean Value (µg/L)	Standard Deviation (µg/L)	VC (%)
Sample 1	16	246	13.09	5.32
Sample 2	16	684	52.27	7.64
Sample 3	16	679	93.61	13.79

#### Inter-Assay-Variability

Serum samples were measured in independent assays. On average, the coefficient of variation was  $\emptyset < 10\%$ . Exemplary results are detailed in Table 2.

**Table 2 Inter-Assay Variability** IGF-I concentrations were determined in duplicate in independent assays and variability was calculated as the coefficient of variation (CV).

	Number of Determinations	Mean Value (µg/L)	Standard Deviation (µg/L)	VC (%)
Sample 1	24	291	20	7
Sample 2	26	695	51	7
Sample 3	23	773	76	10
Sample 4	23	256	15	6
Sample 5	26	151	21	14
Sample 6	26	444	37	8
Sample 7	26	127	12	9
Sample 8	26	686	59	9
Sample 9	26	581	50	9
Sample 10	26	178	12	7

### 10.3 The Comparison of Methodes


Figure 4 Method comparison of the Meddiagnost E25 m/r IGF-I ELISA and the IDS m/r IGF-I ELISA.

### 10.4 Linearity


Figure 5 Linearity. Comparison of the theoretically expected and the measured IGF-I concentration in the respective dilutions. From a 1:10 dilution to a concentration of 0.5 ng/mL, the linearity of the sample dilution is given

## 11 LITREATURE / LITERATUR

- 1) Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem*.25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors*. New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology*. 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem*. 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab*. 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh)*. 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding pro-teins. *J. Clin Endocrinol Metab*. 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem*. 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J*. 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest*. 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med*. 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol*. 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab*. 51:781-788.
- 15) Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr*. 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med*. 301:1138-1142
- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest*. 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence*. Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab*. 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh)*.118:374-380.
- 22) Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. *J Clin Endocrinol Metab* Vol 82, 4:1301-1303
- 23) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Horm Res* 54:60-68
- 24) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Horm Res* 55:155-124

## 12 INTERNATIONAL ASSAY DESCRIPTION

A-E	STD	Rec in 1 mL BUF PP	-
KS1	Control	Rec in 500 µL BUF PP	1:100 DILU BUF PP
KS2	Control	Rec in 500 µL BUF PP	1:100 DILU BUF PP
WP	WASHBUF 20x	-	1:20 DILU A. dest.
-	SPE		1:100 DILU BUF PP
-	°C 20-25 °C ; A 0.25 h ; « max 350 rpm		
50 µL	Ab AK		A1 - End
50 µL	BUF PP		A1/A2
50 µL	STD A (0.5 ng/ml)		B1/B2
50 µL	STD B (2.5 ng/ml)		C1/C2
50 µL	STD C (6 ng/ml)		D1/D2
50 µL	STD D (12 ng/ml)		E1/E2
50 µL	STD E (18 ng/ml)		F1/F2
50 µL	CONTROL KS 1 1:100 DILU BUF PP		G1/G2
50 µL	CONTROL KS 2 1:100 DILU BUF PP		H1/H2
50 µL	SPE 1:100 DILU BUF PP		
TAPE			
A 1 h °C 20-25 « 350 rpm			
5x 300 µL	5x WASHBUF WP		
100 µL	CONJ EK		
TAPE			
A 30 min °C 20-25 « 350 rpm			
5x 300 µL	5x WASHBUF WP		
100 µL	SUBST TMB S		
A 30 min °C 20-25 			
H <sub>2</sub> SO <sub>4</sub> SL			
MEASURE			