# Resistin ELISA

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**Table of Contents**

<b>1</b>	<b>INTENDED USE</b>	<b>18</b>
<b>2</b>	<b>INTRODUCTION</b>	<b>18</b>
<b>3</b>	<b>METHOD</b>	<b>19</b>
<b>4</b>	<b>WARNINGS AND PRECAUTIONS</b>	<b>20</b>
<b>5</b>	<b>MATERIALS</b>	<b>21</b>
<b>6</b>	<b>SAMPLES</b>	<b>22</b>
<b>7</b>	<b>TECHNICAL NOTES</b>	<b>23</b>
<b>8</b>	<b>ASSAY PROCEDURE</b>	<b>24</b>
<b>9</b>	<b>QUALITY CONTROL</b>	<b>25</b>
<b>10</b>	<b>ESTABLISHING THE CALIBRATION CURVE</b>	<b>25</b>
<b>11</b>	<b>PERFORMANCE CHARACTERISTICS</b>	<b>27</b>
<b>12</b>	<b>EVALUATION OF RESULTS</b>	<b>28</b>
<b>13</b>	<b>SCIENTIFIC APPLICATION</b>	<b>29</b>
<b>14</b>	<b>LITERATUR / LITERATURE</b>	<b>30</b>

## ENGLISH Resistin ELISA E50 Instructions for Use

Resistin ELISA E50	96 Determination
CE	DE/CA40/00809/14/1
Principle of the test	Enzyme-linked Immunoassay
Duration (incubation periods)	4 h
Antibodies	Specific polyclonal Antibodies
Calibrators	5 single Calibrators: 0.02 - 1 ng/mL, recombinant Resistin
Assay Range	0.012 – 21 ng/mL
Controls	2 Controls, lyophilised
Samples	human serum / plasma
Required sample volume	15 µL
Sample dilution	1:21
Analytical sensitivity	∅ 0.012 µg/L
Inter- / Intra Assay	∅ < 10 %

### 1 INTENDED USE

Measurement of human Resistin in human Serum and Plasma Sample.

### 2 INTRODUCTION

Resistin, a cysteine-rich protein of 11.3 kDa (1), was firstly found in mice (2) and constitutes together with RELM $\alpha$ , RELM $\beta$  and RELM $\gamma$  the protein family of resistin-like molecules (RELM).

In humans, resistin and RELM $\beta$  (1) but no other proteins of the RELM family were found. The human form of resistin shows a homology of 53% to the murine protein (4). It has 11 cysteine-residues, is synthesized as a propeptide of 108 amino acids and secreted as a dimer, build by a disulfide bridge of cysteine residues (22). Beside this intermolecular disulfide bridge, 5 additional intramolecular ones exist (5,6).

Appearance of multi- and oligomer formation was proved by size exclusion chromatography. Thereby it was shown, that oligomer formation is SDS-insensitive but can be inhibited by  $\beta$ -mercaptoethanol and is therefore likely to be caused by disulfide bridges (1). Further on, the resistin structure seems to be dependent on its concentration, as circular dichroism analysis shows a concentration dependent shift of  $\alpha$ -helical to  $\beta$ -sheet structure (1).

Resistin expression was demonstrated in white adipose tissue (10), pituitary (11) and pancreatic islets (12) of mice as well as in brown adipose tissue of rats. In humans, resistin expression in adipocytes can be detected but only at a very low level. But in vitro, resistin expression of non-adipocytes in fatty tissue was shown (13). Human resistin gene is also expressed in pancreatic islets (12), pre-adipocytes (14) macrophages (15) and bone marrow (39). So, resistin is of relevance for inflammation processes as well as for lipid metabolism.

Most investigation refers to the mouse model. Here, the existence of trimeric and hexameric resistin in serum was demonstrated (7). In comparison to adiponectin biology, it is highly probable that different resistin oligomers have different biologic function (8, 9).

In mice, a correlation between adiposity, insulin resistance and resistin expression was found empirically. In humans, respective study results are not clear – several studies show an association of resistin serum concentration and adiposity or insulin resistance (17, 25-31). But

## Exemplary version, do not use to perform assays

others failed in confirming these results (14, 16-24). Therefore, there is requirement for valid and reproducible determination of resistin serum concentration.

Relevance of resistin in other physiologic processes than energy metabolism was investigated by several different approaches. Experiments with endothelial cells gave interesting results. Here, resistin was shown to enhance expression of VCAM-1 and ICAM-1 (33, 34). By this way, resistin is potentially able to influence endothelial inflammation (35, 36) and, thereby atherosclerosis. These results were confirmed by experiments in mice, where endothelin-1 was shown to regulate resistin secretion (37, 38).

In recent research human resistin was shown to increase pre-adipocyte proliferation and lipolysis of mature adipocytes (38). By the way of modulating MAPK-signalling pathways resistin exerts crucial influence on energy metabolism.

Present research demonstrates, that Resistin exerts influence on a broad variety of physiological processes, however a clear and defined biological role of resistin remains still unexisting.

This ELISA-kit enables the user to determine the exact concentration of Resistin in human serum/plasma as well as other body fluids and thereby assists investigation of Resistin biology.

With optimized special sample buffer: reliable quantification of resistin - independent of the absolute content.

### **3 METHOD**

The enzyme immunoassay for Resistin E50 is a so-called Sandwich-Assay. It utilizes a specific high affinity polyclonal rabbit antiserum coated on the wells of a microtiter plate. The Resistin in the samples binds quantitatively to the immobilized antiserum. In the following step, the biotinylated antiserum binds in turn to Resistin. After washing, Streptavidin-Peroxidase-Enzyme conjugate will be added, which will bind highly specific to the biotin of the antiserum and will catalyse in the closing substrate reaction the turn of the colour, quantitatively depending on the Resistin level of the samples.

### **4 WARNINGS AND PRECAUTIONS**

**For In Vitro Diagnostic Use only. For Professional use only.**

## Exemplary version, do not use to perform assays

The Mediagnost kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

**Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.**

Do not use obvious damaged or microbial contaminated or spilled material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

### Human Serum

Following components contain human serum: **Controls CTR1 and CTR2**

Source human serum for the Controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

### Reagents CAL A-E, DET, EC, SB, DIL, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

### Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

### Stop Solution STP

The Stop solution contains 0.2 M acid sulphur acid (H<sub>2</sub>SO<sub>4</sub>)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

## 4.1 General first aid procedures:

**Skin contact:** Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

**Eye contact:** In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

**Ingestion:** After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

## 5 MATERIALS

### 5.1 Materials provided

The reagents listed below are sufficient for 96 wells including the Calibration curve.

<b>MTP</b>	<b>Microtiter plate</b> , ready for use, coated with rabbit-anti-resistin-antibody. Wells are separately breakable.	<b>(8x12) wells</b>
<b>CAL A-E</b>	<b>Calibrators</b> , lyophilized, (recombinant human Resistin), concentrations are given on vial labels and on quality certificate. <b>Reconstitute in 750 µL Sample Buffer SB.</b>	<b>5 x 750 µL</b>
<b>CTR1</b>	<b>Control 1</b> , lyophilised, (human serum), concentration is given on quality certificate. <b>Reconstitute in 250 µL Dilution Buffer DIL.</b> <b>Dilute with Sample Buffer SB.</b>	<b>1 x 250µL</b>
<b>CTR2</b>	<b>Control 2</b> , lyophilised, (human serum), concentration is given on quality certificate. <b>Reconstitute in 250 µL Dilution Buffer DIL.</b> <b>Dilute with Sample Buffer SB.</b>	<b>1 x 250 µL</b>
<b>DET</b>	<b>Antibody Conjugate, 100fold concentrated solution</b> , contains rabbit biotinylated anti-resistin antibody. <b>Before use, dilute 1:100 with Dilution Buffer DIL.</b>	<b>1 x 120 µL</b>
<b>EC</b>	<b>Enzyme Conjugate, 100fold concentrated solution</b> , contains HRP (Horseradish-Peroxidase)-labelled Streptavidin. <b>Before use, dilute 1:100 with Dilution Buffer DIL.</b>	<b>1 x 120 µL</b>
<b>SB</b>	<b>Sample Buffer</b> , ready for use.	<b>1 x 120 mL</b>
<b>DIL</b>	<b>Dilution Buffer</b> , ready for use.	<b>1 x 25 mL</b>
<b>WB</b>	<b>Washing Buffer</b> , 20-fold concentrated solution. Dilute 1:20 with Aqua dest.	<b>1 x 50 mL</b>
<b>S</b>	<b>Substrate</b> , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine.	<b>1 x 12 mL</b>
<b>STP</b>	<b>Stop Solution</b> , ready for use, 0.2 M sulphuric acid.	<b>1 x 12 mL</b>
-	<b>Sealing Tape</b> , for covering the <b>microtiter plate</b> .	<b>3 x</b>
	<b>Instructions for use</b>	<b>1 x</b>
-	<b>Quality Certificate</b>	<b>1 x</b>

### 5.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, **950 mL**.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm

## 6 SAMPLES

### 6.1 Sample type

Serum and Plasma

Serum as well plasma samples are suitable (significant deviation of Resistin levels in corresponding serum-, Heparin-, EDTA-, Citrate-plasma-Samples were not found. A possible dilution of the sample by the anticoagulant must be taken into account).

Haemolytic samples appear to show falsely high Resistin levels, using such samples should be checked out critically

By means of the special sample buffer an external sample preparation prior to the assay is not required.

### 6.2 Specimen collection

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

### 6.3 Required sample volume: 15 µL

### 6.4 Sample stability

In firmly closable plastic tubes.

- Storage at 20 -25°C: max. 2 days
- Storage at -20°C min. 2 years

The storage of samples, over a period of 2 years at -20°C, showed no effect on the measured value. Freeze/thaw cycles of samples should be minimized.

### 6.5 Interference

Triglyceride, bilirubin and hemolysate in the sample do not interfere to a concentration of 100 mg/mL and 100 µg/mL or 1 mg/mL. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

### 6.6 Sample dilution

- Dilution: **1:21** with Sample Buffer **SB**
- Pipette 300 µL Sample Buffer SB in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add 15 µL Serum- or Plasma (dilution 1:21). After mixing use 2 x 100 µL of this dilution in the assay
- According to expected Resistin levels the dilution with SB can be higher or lower. The excellent linearity of this test system allows sample dilution of 1:5 to 1:400.
- Because the Sample Buffer **SB** is special composed for the correct determination of Resistin, the dilution should be **at least 1:5!**

## 7 TECHNICAL NOTES

### Storage Conditions

## Exemplary version, do not use to perform assays

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at –20 °C after reconstitution. Avoid repeated thawing and freezing.

### Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components**: Calibrators **A-E** and Controls **CTR1** and **CTR2** must be stored at –20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C.

### Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Reconstitution

The Calibrators **A – E** are reconstituted with the Sample Buffer **SB**.

The Controls **CTR1** and **CTR2** are reconstituted with the Dilution Buffer **DIL**.

It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The reconstituted calibrators and controls can be stored for 4 weeks at –20°C. Repeated freeze/thaw cycles have to be avoided.

### Dilution

After **reconstitution**, dilute the Controls **CTR1** and **CTR2** with the Sample Buffer **SB** in the same ratio (**1:21**) as the samples.

The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

The diluted Washing Buffer is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage!

### Antibody and Enzyme Conjugate

Use the Dilution Buffer **DIL** for the dilution of the Antibody Conjugate **DET** and Enzyme Conjugate **EC** 100fold concentrates. The diluted solutions are only limited stable at 2-8°C and should be prepared daily fresh.

### Assay Procedure

When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **DET** and the Enzyme Conjugate **EC** as well as the succeeding Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**.

All determinations (Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

### Incubation

**Incubation at room temperature means: Incubation at 20 - 25°C.** The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive—storage and incubation in the dark.

### Shaking

The incubation steps should be performed at mean rotation frequency of a microtiter plate shaker. We recommend 350 rpm. Depending on the design of the shaker, the shaking frequency should be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

### Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

## 8 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
<b>CAL A-E</b>	<b>Calibrators</b>	in <b>750 µL</b> Sample Buffer <b>SB</b>	-
<b>CTR1</b>	<b>Control 1</b>	in <b>250 µL</b> Dilution Buffer <b>DIL</b>	<b>1:21</b> with Sample Buffer <b>SB</b>
<b>CTR2</b>	<b>Control 2</b>	in <b>250 µL</b> Dilution Buffer <b>DIL</b>	<b>1:21</b> with Sample Buffer <b>SB</b>
<b>DET</b>	<b>Antibody Conjugate conc.</b>	-	<b>1:100</b> with Dilution Buffer <b>DIL</b> → <b>DET 1:100</b>
<b>EC</b>	<b>Enzyme Conjugate conc.</b>	-	<b>1:100</b> with Dilution Buffer <b>DIL</b> → <b>EC 1:100</b>
<b>WB</b>	<b>Washing Buffer conc.</b>	-	<b>1:20</b> with <b>Aqua dest.</b> → <b>WB 1:20</b>
<b>Samples + Controls CTR1 and CTR2: dilute 1:21 in Sample Buffer SB, mix immediately</b>			
Before assay procedure bring all reagents to room temperature ( <b>20-25°C</b> ).			
<b>Assay procedure in double determination</b>			
Pipette	Reagents	Position	
100 µL	Sample Buffer <b>SB</b> (blank value)	A1/A2	
100 µL	Calibrator <b>A (0.02 ng/mL)</b>	B1/B2	
100 µL	Calibrator <b>B (0.1 ng/mL)</b>	C1/C2	
100 µL	Calibrator <b>C (0.3 ng/mL)</b>	D1/D2	
100 µL	Calibrator <b>D (0.6 ng/mL)</b>	E1/E2	
100 µL	Calibrator <b>E (1.0 ng/mL)</b>	F1/F2	
100 µL	Control <b>CTR1</b> (1:21 diluted)	G1/G2	
100 µL	Control <b>CTR2</b> (1:21 diluted)	H1/H2	
100 µL	Sample <b>SPE</b> (1:21 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
<b>Sample Incubation: 2 h at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5x</b> with <b>300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Antibody Conjugate <b>DET 1:100</b>	In each well	
Cover the wells with the sealing tape.			
<b>Incubation: 1 h at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5x</b> with <b>300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Enzyme Conjugate <b>EC 1:100</b>	In each well	
Cover the wells with the sealing tape.			
<b>Incubation: 30 Minutes at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5x</b> with <b>300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Substrate <b>S</b>	In each well	
<b>Incubation: 30 Minutes at 20-25°C, in the dark at RT</b>			
100 µL	Stop Solution <b>STP</b>	In each well	

Measure the absorbance within 30 min at **450 nm** with  $\geq 590$  nm as reference wavelength.

## 9 QUALITY CONTROL

Good laboratory practice (GLP) requires that controls are included in each assay. The test results are only valid if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP or other applicable federal, state or local standards/laws. All Calibrators and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated.

### 9.1 Quality Criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below **0.25**, and the absorbance of Calibrator **E** should be above **1.0**.

Samples, which yield higher absorbance values than **Calibrator E**, should be re-tested with a higher dilution.

## 10 ESTABLISHING THE CALIBRATION CURVE


The calibrators provided contain the following concentrations of recombinant Resistin:

Calibrator	A	B	C	D	E
ng/mL	0.02	0.10	0.30	0.60	1.00
pg/mL	20	100	300	600	1000

- 1) Calculate the mean absorbance value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all samples, controls and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The **Resistin concentrations** of the diluted samples or the diluted control sera in ng/mL (or pg/mL according the chosen unit for the calibrators) are calculated in this way, the Resistin concentrations of the **undiluted samples** and of controls are calculated **by multiplication with the respective dilution factor**.

Exemplary version, do not use to perform assays

### E50 Calibration Curve

**Fig. 1. Exemplary Calibration Curve** with a polynomial 2<sup>nd</sup> degree as curve fit.

The exemplary shown calibration curve in Fig.1 **cannot** be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

Exemplary calculation of the Resistin concentration of a 1:21 diluted sample:

Measured extinction of your sample	0.85
Measured extinction of the blank	0.05

Your measurement program will calculate the Resistin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2<sup>nd</sup> degree).

In this exemplary case the following equation is solved by the program to calculate the Resistin concentration in the sample:

$$y = 0.0046207 + 3.3674 x - 1.5709x^2$$
$$0.2686 = x$$

if the dilution factor (1:21) is taken into account the Resistin concentration of the undiluted sample is

$$0.2686 \times 21 = 5,64 \text{ ng/mL} = 0.00564 \text{ } \mu\text{g/mL}$$

## 11 PERFORMANCE CHARACTERISTICS

### 11.1 Calibrators

The calibrators are prepared from recombinant human Resistin (19.5 kDa, 2 x 92 amino acids, expressed in *E. coli*) in concentrations of 20, 100, 300, 600 and 1000 pg/mL (pico Gramm / mL, equal to 0.02 ng/mL-1 ng/mL).

### 11.2 Sensitivity

The analytical sensitivity was determined by measuring the blank value. For this purpose, the signal intensity of blank value + 2SD was converted into a resistin concentration using the calibration curve. The analytical sensitivity of the assay is 0.012 ng / mL.

### 11.3 Interference

Interference of physiological appearing substance with the Resistin measurement was investigated. Serum samples have been enriched with different concentrations of possibly interfering substances and the amount of Resistin was measured and compared with the Resistin concentration in the same sample without any enrichment. In table 1 the relative results are shown. None of the tested substances interfered significantly with Resistin measurement.

**Table 1:** Interference: Three serum samples where enriched with indicated amount of the potentially interfering substance and measured. Shown is % of Resistin of the native, non-enriched serum sample

	Triglyceride 100 mg/mL	Bilirubin 100 µg/mL	Haemolysate 1 mg/mL
Serum 1	101	93	94
Serum 2	115	99	99
Serum 3	104	103	147

**Table 2:** Effects of coagulation inhibitors were investigating by adding indicated amounts of inhibitors to SB enriched with 0.3 ng/mL Resistin. Relative amounts of Resistin measured in inhibitor containing samples in comparison to 0.3 ng/mL Resistin containing Sample Buffer **SB** are shown.

		% of Resistin in SB	
		Mean (n=3)	SD
3.8 g/L	Citrate	94	7.67
0.0068 mol/L	EDTA	93	4.96
30.000 IU/L	Heparin	96	4.89

### 11.4 Reproducibility and Precision

The inter- and intra assay coefficients of variability are in average < 10%. Exemplary determinations are shown in table 3 and table 4.

**Table 3:** Inter-Assay-Variation

	Number of determinations	Mean value (ng/mL)	Standard deviation (ng/mL)	VC (%)
Sample 1	26	2.91	0.16	5.55
Sample 2	15	4.58	0.24	5.33
Sample 3	17	4.60	0.23	5.04
Sample 4	7	2.50	0.09	3.37
Sample 5	23	4.09	0.27	6.67

**Table 4:** Intra-Assay-Variation

	Number of determinations	Mean value (ng/mL)	Standard deviation (ng/mL)	VC (%)
Sample 1	16	2.81	0.13	4.49
Sample 2	15	4.79	0.24	4.97

### 11.5 Recovery and Linearity

The Mediagnost Resistin ELISA E50 is over a very wide range dilution authentic. The linearity of serum dilutions is over a very wide range excellent (s.Tab.5).

**Table 5:** Recovery and linearity of the Sample Dilution (characteristic results of two different sera)

Dilution	Sample 1 (native 5.5 ng/mL)		Sample 2 (native 2.25 ng/mL)	
	plus 5 ng/mL	Recovery (%)	plus 12.25 ng/mL	Recovery (%)
1:50	9.71	92.5	14.99	103.4
1:100	10.60	101.0	13.64	94.1
1:200	10.44	99.4	14.10	97.2
1:400	10.32	98.3	14.33	98.8

Different human sera were spiked with recombinant human Resistin in varying concentrations (e.g. in Table 5). The recovery of Resistin yielded on average 98 % of the theoretically expected amount.

## 12 EVALUATION OF RESULTS

**Table 6:** The expected values for Resistin were determined with the Mediagnost ELISA E50 in healthy probands and analysed by Prof. Dr. J. Kratzsch. Institute for Laboratory Medicine, University of Leipzig.

Female				Resistin (ng/mL):		
Age (Years):	n:	AV Age:	AV BMI:	AV ± SD:	25.- 75. Percentile:	Min. – Max.:
18 - 30	96	23.0	23.1	7.2 ± 2.6	5.4 – 8.8	3.1 – 14.7
31 - 40	63	36.5	24.3	8.1 ± 2.3	6.4 – 9.6	3.6 – 13.1
41 - 50	67	44.9	24.8	7.3 ± 2.5	5.7 – 8.1	4.0 – 16.1
51 - 60	29	54.7	25.0	7.2 ± 2.6	5.4 – 8.5	4.0 – 15.5
61 - 65	9	62.7	25.2	6.6 ± 1.1	6.0 – 6.7	5.4 – 9.3
Male				Resistin (ng/mL):		
Age (Years):	n:	AV Age:	AV BMI:	AV ± SD:	25.- 75. Percentile:	Min. – Max.:
18 - 30	107	23.9	24.1	6.4 ± 1.8	5.0 – 7.6	2.5 – 13.1
31 - 40	59	35.9	25.0	6.7 ± 3.2	4.8 – 7.4	3.8 – 26.9
41 - 50	66	45.0	25.2	6.5 ± 2.8	4.5 – 7.4	2.4 – 16.7
51 - 60	36	54.8	26.4	6.1 ± 2.1	4.7 – 7.2	3.2 – 13.3
61 - 68	20	63.2	25.6	7.2 ± 1.8	6.0 – 8.2	4.5 – 11.2

n=Number of Probands. AV=Average Value. BMI=Body Mass Index (kg/m<sup>2</sup>). SD=Standard Deviation

**Table 7:** Grouping of the expected values

Sex	Number	Mean [ng/mL]	Standard deviation	2.5. Percentile	9.5. Percentile
Male	288	6.48	2.44	3.32	11.68
Female	264	7.41	2.47	3.68	13.60
Total	552	6.93	2.49	3.58	13.12

## Instructions for use for scientific application

### 13 SCIENTIFIC APPLICATION

Resistin is relevant e.g. in research of: Adiposity Insulin Resistance, Diabetes Arteriosclerosis Inflammation. Due to Assays high sensitivity (12 pg/mL) this is also suitable for measurements in cell culture supernatants and in non-serum / plasma samples (e.g. in breast milk, urine, saliva).

#### 13.1 Samples suitable for scientific application

Serum and Plasma, common cell culture medium, saliva, breast milk and urine were found to be suitable specimens too.

Resistin concentrations may be completely different in body fluids of human origin other than serum or cell culture supernatants.

(See Table 8 for recovery and dilution recommendations). Due to the special Sample Buffer, external sample pre-treatment is not necessary.

**Table 8:** Samples were enriched with 0.3 ng/mL Resistin and measured in comparison to non-enriched sample. Relative recovery of added Resistin is shown.

Matrix	Dilution	% Recovery
Cerebrospinal fluid	1:2	129
Cerebrospinal fluid	1:10	93
Cerebrospinal fluid	1:40	103
Amnion fluid	1:10	85
Amnion fluid	1:40	91
Saliva	1:10	99
Saliva	1:21	86
Urine	1:10	79
Urine	1:21	85
Breast milk	1:2	97
Breast milk	1:10	58
Breast milk	1:21	63
Cell culture supernatant	1:2	100

#### 13.2 Cross reactivity

The cross-reactivities of several commercially available animal sera were tested in this assay.

Commercially available sera from cattle, cat, chicken, dog, donkey, goat, guinea pig, horse, mouse, pig, rabbit, rat and sheep were diluted 1:10 and used in the test and the signal intensity was measured. No cross-reactivity could be detected.

## 14 LITERATUR / LITERATURE







### Literature

1. Aruna B. Ghosh S. Singh AK. Mande SC. Srinivas V. Chauhan R. Ehtesham NZ. Human recombinant resistin protein displays a tendency to aggregate by forming intermolecular disulfide linkages. *Biochemistry* 2003;42(36):10554-9.
2. Steppan CM. Lazar MA. Resistin and obesity-associated insulin resistance. *Trends Endocrinol Metab* 2002;13(1):18-23.
3. Yang RZ. Huang Q. Xu A. McLenithan JC. Eisen JA. Shuldiner AR. Alkan S. Gong DW. Eison JA. Comparative studies of resistin expression and phylogenomics in human and mouse. *Biochem Biophys Res Commun* 2003;310(3):927-35.
4. Banerjee RR. Lazar MA. Resistin: molecular history and prognosis. *J Mol Med* 2003;81(4):218-26.
5. Banerjee RR. Lazar MA. Dimerization of resistin and resistin-like molecules is determined by a single cysteine. *J Biol Chem* 2001;276(28):25970-3.
6. Raghu P. Ghosh S. Soundarya K. Haseeb A. Aruna B. Ehtesham NZ. Dimerization of human recombinant resistin involves covalent and noncovalent interactions. *Biochem Biophys Res Commun* 2004;313(3):642-6.
7. Patel SD. Rajala MW. Rossetti L. Scherer PE. Shapiro L. Disulfide-dependent multimeric assembly of resistin family hormones. *Science* 2004;304(5674):1154-8.
8. Kim KH. Zhao L. Moon Y. Kang C. Sul HS. Dominant inhibitory adipocyte-specific secretory factor (ADSF)/resistin enhances adipogenesis and improves insulin sensitivity. *Proc Natl Acad Sci U S A* 2004;101(17):6780-5.
9. Chen J. Wang L. Boeg YS. Xia B. Wang J. Differential dimerization and association among resistin family proteins with implications for functional specificity. *J Endocrinol* 2002;175(2):499-504.
10. O. Beltowski J. Adiponectin and resistin--new hormones of white adipose tissue. *Med Sci Monit* 2003;9(2):RA55-61.
11. Morash BA. Ur E. Wiesner G. Roy J. Wilkinson M. Pituitary resistin gene expression: effects of age, gender and obesity. *Neuroendocrinology* 2004;79(3):149-56.
12. Minn AH. Patterson NB. Pack S. Hoffmann SC. Gavrilova O. Vinson C. Harlan DM. Shalev A. Resistin is expressed in pancreatic islets. *Biochem Biophys Res Commun* 2003;310(2):641-5.
13. Fain JN. Cheema PS. Bahouth SW. Lloyd Hiler M. Resistin release by human adipose tissue explants in primary culture. *Biochem Biophys Res Commun* 2003;300(3):674-8.
14. Janke J. Engeli S. Gorzelniak K. Luft FC. Sharma AM. Resistin gene expression in human adipocytes is not related to insulin resistance. *Obes Res* 2002;10(1):1-5.
15. Patel L. Buckels AC. Kinghorn IJ. Murdock PR. Holbrook JD. Plumpton C. Macphee CH. Smith SA. Resistin is expressed in human macrophages and directly regulated by PPAR gamma activators. *Biochem Biophys Res Commun* 2003;300(2):472-6.
16. Youn BS. Yu KY. Park HJ. Lee NS. Min SS. Youn MY. Cho YM. Park YJ. Kim SY. Lee HK. Park KS. Plasma resistin concentrations measured by enzyme-linked immunosorbent assay using a newly developed monoclonal antibody are elevated in individuals with type 2 diabetes mellitus. *J Clin Endocrinol Metab* 2004;89(1):150-6.
17. Chen XD. Lei T. Xia T. Gan L. Yang ZQ. Increased expression of resistin and tumour necrosis factor-alpha in pig adipose tissue as well as effect of feeding treatment on resistin and cAMP pathway. *Diabetes Obes Metab* 2004;6(4):271-9.
18. Volarova de Courten B. Degawa-Yamauchi M. Considine RV. Tataranni PA. High serum resistin is associated with an increase in adiposity but not a worsening of insulin resistance in Pima Indians. *Diabetes* 2004;53(5):1279-84.
19. Shalev A. Patterson NB. Hirshberg B. Rother KI. Harlan DM. Resistin serum levels in type 1 diabetes pre- and post-islet transplantation. *Metabolism* 2004;53(4):403-4.
20. Fujinami A. Obayashi H. Ohta K. Ichimura T. Nishimura M. Matsui H. Kawahara Y. Yamazaki M. Ogata M. Hasegawa G. Nakamura N. Yoshikawa T. Nakano K. Ohta M. Enzyme-linked immunosorbent assay for circulating human resistin: resistin concentrations in normal subjects and patients with type 2 diabetes. *Clin Chim Acta* 2004;339(1-2):57-63.
21. Degawa-Yamauchi M. Bovenkerk JE. Juliar BE. Watson W. Kerr K. Jones R. Zhu Q. Considine RV. Serum resistin (FIZZ3) protein is increased in obese humans. *J Clin Endocrinol Metab* 2003;88(11):5452-5.
22. Azuma K. Katsukawa F. Oguchi S. Murata M. Yamazaki H. Shimada A. Saruta T. Correlation between serum resistin level and adiposity in obese individuals. *Obes Res* 2003;11(8):997-1001.
23. Wang H. Chu WS. Hemphill C. Elbein SC. Human resistin gene: molecular scanning and evaluation of association with insulin sensitivity and type 2 diabetes in Caucasians. *J Clin Endocrinol Metab* 2002;87(6):2520-4.
24. Brichard SM. Delporte ML. Lambert M. Adipocytokines in anorexia nervosa: a review focusing on leptin and adiponectin. *Horm Metab Res* 2003;35(6):337-42.
25. Kawanami D. Maemura K. Takeda N. Harada T. Nojiri T. Imai Y. Manabe I. Utsunomiya K. Nagai R. Direct reciprocal effects of resistin and adiponectin on vascular endothelial cells: a new insight into adipocytokine-endothelial cell interactions. *Biochem Biophys Res Commun* 2004;314(2):415-9.
26. Verma S. Li SH. Wang CH. Fedak PW. Li RK. Weisel RD. Mickle DA. Resistin promotes endothelial cell activation: further evidence of adipokine-endothelial interaction. *Circulation* 2003;108(6):736-40.
27. Kaser S. Kaser A. Sandhofer A. Ebenbichler CF. Tilg H. Patsch JR. Resistin messenger-RNA expression is increased by proinflammatory cytokines in vitro. *Biochem Biophys Res Commun* 2003;309(2):286-90.
28. Lu SC. Shieh WY. Chen CY. Hsu SC. Chen HL. Lipopolysaccharide increases resistin gene expression in vivo and in vitro. *FEBS Lett* 2002;530(1-3):158-62.
29. Zhong Q. Lin CY. Clarke KJ. Kempainen RJ. Schwartz DD. Judd RL. Endothelin-1 inhibits resistin secretion in 3T3-L1 adipocytes. *Biochem Biophys Res Commun* 2002;296(2):383-7.
30. Barton M. Carmona R. Ortmann J. Krieger JE. Traupe T. Obesity-associated activation of angiotensin and endothelin in the cardiovascular system. *Int J Biochem Cell Biol* 2003;35(6):826-37.
31. Ort TA. A. A.; MacDougall. J.R.; Nelson. P.J.; Rothenberg. M.E.; Wu. F.; Eisen. A.; Halvorsen. Y.-D.C. Recombinant human FIZZ3/Resistin Stimulates Lipolysis in Cultured Human Adipocytes. Mouse Adipose Explants and in Normal Mice. *Endocrinology* 2005; doi:10.1210/en.2004-1421.

Exemplary version, do not use to perform assays

## INTERNATIONAL ASSAY PROCEDURE

### International Test Description

CAL A-E	Rec in 750 µL SB	-
CTR1	Rec in 250 µL DIL	1:21 SB
CTR2	Rec in 250 µL DIL	1:21 SB
DET	-	1:100 DIL → DET 1:100
EC	-	1:100 DIL → EC 1:100
WB 20x	-	1:20 A. dest. → WB 1:20
SPE		1:21 SB
-	°C 20-25°C,  	
100 µL	SB	A1/A2
100 µL	CAL A (0.02 ng/mL)	B1/B2
100 µL	CAL B (0.1 ng/mL)	C1/C2
100 µL	CAL C (0.3 ng/mL)	D1/D2
100 µL	CAL D (0.6 ng/mL)	E1/E2
100 µL	CAL E (1 ng/mL)	F1/F2
100 µL	CTR1 1:21 SB	G1/G2
100 µL	CTR2 1:21 SB	H1/H2
100 µL	SPE 1:21 SB	
<b>TAPE</b>		
 2 h °C 20-25°C  350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	DET 1:100	
<b>TAPE</b>		
 1 h °C 20-25°C  350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	EC 1:100	
<b>TAPE</b>		
 30 min °C 20-25°C  350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	S	
 30 min °C 20-25°C 		
100 µL	STP	
<b>MEASURE</b>		

## ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
<b>CAL A-E</b>	<b>Calibrators</b>	in <b>750 µL</b> Sample Buffer <b>SB</b>	-
<b>CTR1</b>	<b>Control 1</b>	in <b>250 µL</b> Dilution Buffer <b>DIL</b>	<b>1:21</b> with Sample Buffer <b>SB</b>
<b>CTR2</b>	<b>Control 2</b>	in <b>250 µL</b> Dilution Buffer <b>DIL</b>	<b>1:21</b> with Sample Buffer <b>SB</b>
<b>DET</b>	<b>Antibody Conjugate Conc.</b>	-	<b>1:100</b> with Dilution Buffer <b>DIL</b> → <b>DET 1:100</b>
<b>EC</b>	<b>Enzyme Conjugate Conc.</b>	-	<b>1:100</b> with Dilution Buffer <b>DIL</b> → <b>EC 1:100</b>
<b>WB</b>	<b>Washing Buffer conc.</b>	-	<b>1:20</b> with <b>Aqua dest.</b> → <b>WB 1:20</b>
<b>Samples + Controls CTR1 and CTR2: dilute 1:21 in Sample Buffer SB, mix immediately</b>			
Before assay procedure bring all reagents to room temperature ( <b>20-25°C</b> ).			
<b>Assay procedure in double determination</b>			
Pipette	Reagents	Position	
100 µL	Sample Buffer <b>SB</b> (blank value)	A1/A2	
100 µL	Calibrator <b>A (0.02 ng/mL)</b>	B1/B2	
100 µL	Calibrator <b>B (0.1 ng/mL)</b>	C1/C2	
100 µL	Calibrator <b>C (0.3 ng/mL)</b>	D1/D2	
100 µL	Calibrator <b>D (0.6 ng/mL)</b>	E1/E2	
100 µL	Calibrator <b>E (1.0 ng/mL)</b>	F1/F2	
100 µL	Control <b>CTR1</b> (1:21 diluted)	G1/G2	
100 µL	Control <b>CTR2</b> (1:21 diluted)	H1/H2	
100 µL	Sample <b>SPE</b> (1:21 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
<b>Sample Incubation: 2 h at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5 x with 300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Antibody Conjugate <b>DET 1:100</b>	In each well	
Cover the wells with the sealing tape.			
<b>Incubation: 1 h at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5 x with 300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Enzyme Conjugate <b>EC 1:100</b>	In each well	
Cover the wells with the sealing tape.			
<b>Incubation: 30 Minutes at 20-25°C, 350 rpm</b>			
5x 300 µL	Aspirate the contents of the wells and <b>wash 5 x with 300 µL</b> each Washing Buffer <b>WB 1:20</b> / well	In each well	
100 µL	Substrate <b>S</b>	In each well	
<b>Incubation: 30 Minutes at 20-25°C, in the dark at RT</b>			
100 µL	Stop Solution <b>STP</b>	In each well	
Measure the absorbance within 30 min at <b>450 nm</b> with $\geq 590$ nm as reference wavelength.			

Exemplary version, do not use to perform assays