# Leptin RIA-CT

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## PACKAGE INSERT ENGLISH

### TECHNICAL FEATURES+APPLICATIONS

- ◆ Measures total Leptin concentration in **serum, EDTA/Heparin plasma**
- ◆ Calibrated against the **WHO** International Standard: **NIBSC Code 97/594**
- ◆ High affinity antibodies, no influence of Leptin Bindingprotein
- ◆ ready-to-use calibrators with **0-64 ng/mL** are
- ◆ Small serum volumina sufficient, measurement in capillar blood possible
- ◆ Easy handling

### INTRODUCTION

Leptin, the product of the ob gene (1,2), is a single-chain proteohormone, discovered 1994, with a molecular weight of 16 kD, which is thought to play a key role in the regulation of body weight. Its amino acid sequence exhibits no major homologies with other proteins (1). Leptin is almost exclusively produced by differentiated adipocytes (3-5). It acts on the central nervous system, in particular the hypothalamus, thereby suppressing food intake and stimulating energy expenditure (2,6-9). Leptin receptors - alternatively spliced forms exist that differ in length - belong to the cytokine class I receptor family (10-12). They are found ubiquitously in the body (10, 11, 13, 14) indicating a general role of leptin, which is currently not fully understood. A circulating form of the leptin receptor exists which acts as one of several leptin binding proteins (15). Besides its metabolic effects, leptin was shown to have a strong influence on a number of endocrine axes. In male mice, it blunted the starvation-induced marked decline of LH, testosterone, thyroxine and the increase of ACTH and corticosterone. In female mice, leptin prevented the starvation-induced delay in ovulation (16). Ob/ob mice, which are leptin deficient due to an ob gene mutation, are infertile. This defect could be corrected by administration of leptin, but not through weight loss due to fasting (17), suggesting that leptin is pivotal for reproductive functions.

All these actions may, at least in part, be explained by the suppressive effect of leptin on neuropeptide Y (NPY) expression and secretion by neurons in the arcuate nucleus (6,18,19). NPY is a strong stimulator of appetite (20,21) and is known to be involved in the regulation of various pituitary hormones, e.g. suppression of GH through stimulation of somatostatin (22,23), suppression of gonadotropins (23) or stimulation of the pituitary-adrenal axis (21).

The most important variable that determines circulating leptin levels is body fat mass (24-26). Obviously, under conditions of regular eating cycles, leptin reflects the proportion of adipose tissue (27) showing an exponential relationship (37). This constitutive synthesis of leptin is modulated by a number of non-hormonal and hormonal variables. Stimulators in both rodents and humans are overfeeding (28,29), insulin (3,5,30-33) and glucocorticoids (5,34-36). Suppression has been shown for fasting (27), cAMP and beta-3-adrenoceptor agonists (35). From these findings it becomes clear that leptin is an integral component of various metabolic and endocrine feedback loops (38).

For clinical purposes, it is important to note that serum leptin levels show a moderate circadian variation with a peak during the night at about 2 a.m. (37). The leptin values at this time are about 30 to 100 % higher than the levels measured in the morning or early afternoon. This variation together with the influence of food intake needs to be taken into account, when blood samples are collected.

Under fairly standardized conditions, i.e. normal eating cycles and blood sampling in the morning or early afternoon, a single leptin measurement is informative.

For the appropriate interpretation of measured leptin levels, reference ranges are required. Because body fat mass is the major confounding variable, these ranges should be referred to measures of the percentage body fat such as body mass index (BMI) or percent body fat determined by, e.g., bioelectric impedance assessment (BIA). Leptin levels are higher in females than in males (38,39)

and an age dependence was shown in children and adolescents (40). Therefore, reference ranges referring to measures of body fat should be stratified according to gender and pubertal development.

Leptin levels are high in most obese patients suggesting the presence of leptin insensitivity (20,26,37,38,41,42). In a small percentage of patients, however, leptin levels have been found inappropriately low with respect to their fat mass. It remains for future studies to prove that these patients represent a new pathophysiologic entity: leptin deficiency. Since leptin has also been shown to be of great importance for reproductive functions, possible new pathophysiologic mechanisms may be discovered relating infertility to insufficient leptin production.

The discovery of leptin has released an avalanche of research activities seeking to understand the regulation and actions of this new hormone. Most importantly, it has provided a key to better understand the physiology of body weight regulation and to unveil possible pathophysiologic mechanisms in both obesity and eating disorders. Further, it may provide new insights into certain causes of infertility.

### **INTENDED USE**

This radio immunoassay kit is suited for measuring human leptin in serum or plasma EDTA-/Heparin-Plasma).

Under conditions of normal eating cycles, measurement in a single blood sample collected in the morning or early afternoon is sufficient.

The comparison with BMI-related reference ranges may be useful to detect conditions of relative **leptin deficiency** as a possible cause of obesity or provide an indication for **leptin resistance** respectively.

Due to its high correlation with body fat mass leptin measurements under standardized conditions may be used as a simple and inexpensive test for determination of **body fat**.

### **PRINCIPLE**

For the Radioimmunoassay for the determination of human Leptin a polyclonal rabbit-antibody of high specificity is used. Leptin is measured quantitatively and binding protein does not influence the test results. Tubes were coated with streptavidin and separation is facilitated by biotinylated anti-rabbit-antibody.

Calibrators are prepared of recombinant human leptin in concentrations between **1, 2, 4, 8, 16, 32 and 64 ng/mL**.

### **Calibration of the Assay**

The Assay is calibrated against the International Standard for human leptin. The standard preparation of the WHO with code 97/594 (y) is available from the NIBSC (z). One ampoule of the preparation, reconstituted in 1 mL Assay Buffer, will be quantified with the nominal content of 5000 ng human leptin.

## PERFORMANCE CHARACTERISTICS

### Sensitivity

The analytical **sensitivity** of the assay yields **0.1 ng/mL** (2x SD of zero calibrator)

### Intra-Assay-Variation

|          | Number of determinations | Mean value [ng/mL] | Standard Deviation [ng/mL] | VC% |
|----------|--------------------------|--------------------|----------------------------|-----|
| Sample 1 | 6                        | 19.67              | 0.87                       | 4.4 |
| Sample 2 | 6                        | 6.95               | 0.33                       | 4.8 |

### Inter-Assay-Variation

|          | Mean value (ng/mL) | Standard deviation (ng/mL) | VC% (%) |
|----------|--------------------|----------------------------|---------|
| Sample 1 | 8.975              | 0.48                       | 5.3     |
| Sample 2 | 6.3                | 0.32                       | 5.1     |
| Sample 3 | 20.9               | 1.05                       | 5.0     |

### Linearity

The Leptin RIA Lep R44 can be used for diluted samples. The linearity of serum dilutions is excellent (Fig.9).

### Recovery

Serum spiking experiments with recombinant human leptin yielded a recovery of 97% ( $\pm$  2%).

## SPECIMEN COLLECTION, PREPARATION AND STORAGE

**Serum** as well as EDTA/Heparin **plasma samples** are suitable (significant deviation of hLeptin levels in corresponding Serum or -Plasma samples were not found). An external sample preparation prior to assay is not required.

### Warning statement

Biotin concentrations of >125 ng/mL in R44 showed significant interference with the assay system and result in false high values. Therefore samples from patients under high dose biotin therapy (>10 mg/day) cannot be used for Leptin measurement by this test.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen at  $-20^{\circ}\text{C}$  or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although Leptin levels were found to be unaffected by few cycles (5x) in our experiments.

The high sensitivity of the assay allows measurement of hLeptin in small sample volumes. Because of the wide effective range of this RIA kit a preparative sample dilution is generally not necessary. For most of the determinations (serum or plasma samples, and no extreme values expected) **the use of undiluted samples 25  $\mu\text{L}$  per tube, should be appropriate.** In case Leptin levels of more than 64 ng/mL are expected, e.g. obese patients (BMI>35), the sample should and can be diluted, e.g. 1:10.

The hLeptin concentrations may be completely different in body fluids of human origin other than serum or cell culture supernatants.

## REAGENTS PROVIDED

|  |   |
|--|---|
| <b>T</b>   | <b>Tubes</b> , ready for use: coated with Streptavidin, 125 Tubes.  |
| <b>F-L</b>   | <b>Calibrators</b> , ready for use: <b>750 µL</b> , contain recombinant Leptin. Calibrator values are between 1 – 64 ng/mL (1; 2; 4; 8; 16; 32; 64 ng/mL) recombinant human leptin. Use <b>25 µL</b> / tube.              |
| <b>A</b>   | <b>Assay Buffer</b> , ready for use, <b>125 mL</b> .  |
| <b>N</b>   | <b>Control Serum</b> , lyophilized, <b>750 µL</b> , contains human serum, reconstitute with <b>750 µL Assay Buffer A</b> . Exact human leptin concentration and the <b>acceptable range</b> are given on the certificate. |
| <b>S</b>   | <b>Specific Antibody</b> , lyophilized, <b>7 mL</b> , contains anti-human leptin antibody. Reconstitute in <b>7 mL Assay Buffer A</b> .   |
| <b>R</b>   | <b>Capture-Antibody</b> , lyophilized, <b>7 mL</b> , biotinylated ant-rabbit-IgG, reconstitute in <b>7 mL Assay Buffer A</b> .  |
| <b>C</b>   | <b>Tracer</b> , lyophilized, <b>13 mL</b> , < 2.25 µCi or < 83 kBq, reconstitute in <b>13 mL Assay Buffer A</b> , red coloured.   |
|  | <b>Instructions for use</b>   |
|  | <b>Quality Certificate</b>  |

## MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes (100 and 200 µL) Micropipettes and multichannel pipettes with disposable plastic tips  
Vortex-mixer  
Shaking Device  
Gamma Counter

## REAGENT PREPARATION

In conducting the assay, follow strictly the test protocol. Room temperature incubation means: Incubation at 20 - 25°C.

Reagents with different lot numbers should not be mixed. All reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

**Control N** Reagent **R**, **S** and **C** have to be reconstituted in **Assay Buffer A**. It is recommended to keep the reconstituted reagents at room temperature for 30 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The shelf-life of the reagents after opening is in accordance with the **Tracer C** shelf life.

Reconstituted Components (**Control N** and Reagent **R**, **S**, and **C**) should be stored at -20°C (or below). Repeated freeze-thaw cycles have to be avoided.

Before use, all kit components should be brought to room temperature. **Precipitates, possible in buffers, should be dissolved before use through mixing and warming.**

## WARNINGS AND PRECAUTIONS

### For in-vitro diagnostic use only. For professional use only.

Possession and use of the kit are subject to the regulations of the national nuclear regulatory authorities.

Reagents with different lot numbers should not be mixed.

- Reagents contain sodium azide as preservative, however, highly diluted (0.02%).
- Safety Data Sheet available on request.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to room temperature at 20 - 25°C. Precipitates in buffers should be dissolved before use by thorough mixing and warming. Temperature WILL affect the the assay. However, Values for the patient samples will not be affected.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Source human serum for the Control provided in this kit was tested by recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total

## Exemplary version, do not use to perform assays

assurance of the absence of infectious agents; therefore, all components and patient's specimens should be treated as potentially infectious.

**Radioactivity** - Before ordering or using radioactive materials, it is necessary to take the appropriate actions to ensure compliance with national regulations governing their use. Local rules in each establishment, which define actions and behaviour in the radioactivity working areas, should also be adhered to. The advice given here does not replace any local rules, instructions or training in the establishment, or advice from the radiation protection advisers. It is important to follow the code of good laboratory practice in addition to the specific precautions relating to the radionuclide I-125 used.

Iodine-125 has a radioactive half-life  $T_{1/2}$  of 60 days and emits 35.5 keV gamma radiation, 27 – 32 keV x-rays and no beta radiation. Shielding is effectively done by lead, first half value layer is 0.02 mm lead, reduction to 10 % is made by 0.2 mm.

To reduce the radiation dose time spent handling radioactivity should be minimized (plan ahead), and distance from source of radiation should be maximized (doubling the distance from the source quarters the radiation dose).

Formation of aerosols, e. g. by improper opening and mixing of vials or pipetting of solutions which may cause minute droplets of radioactivity become airborne, is a hazard and should be avoided.

Solutions containing iodine should not be made acidic, because this might lead to the formation of volatile elemental iodine.

As some iodo-compounds can penetrate rubber gloves, it is advisable to wear two pairs, or polyethylene gloves over rubber.

For cleaning of contaminated areas or equipment, the Iodine-125 should be rendered chemically stable by using alkaline sodium thiosulphate solution together with paper or cellulose tissue.

### **General first aid procedures:**

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

### **The handling of radioactive and potentially infectious material must comply with the following guidelines:**

The material should be stored and used in a special designated area.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Avoid direct contact with these materials by wearing laboratory coats and disposable gloves.

Spilled material must be wiped off immediately. Clean contaminated areas and equipment with a suitable detergent.

Unused radioactive material and radioactive waste should be disposed according to the recommendations of the national regulatory authorities.

### **Warning statement**

Biotin concentrations of >125 ng/mL in R44 showed significant interference with the assay system and result in false high values. Therefore, samples from patients under high dose biotin therapy (>10 mg/day) cannot be used for Leptin measurement by this test.

## ASSAY PROCEDURE

NOTES: All determinations (Calibrators, Control and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Calibrators, Control and the samples should be pipette as fast as possible.

### Flow Chart of Assay Protocol

| #     | Tube           | A, F-L, N Samples | R   | S   | C   |
|-------|----------------|-------------------|-----|-----|-----|
| 1,2   | TC             | ---               | --- | --- | 100 |
| 3,4   | B <sub>0</sub> | A: 25             | 50  | 50  | 100 |
| 5-18  | Calibrators    | F-L: 25           | 50  | 50  | 100 |
| 19,20 | Control        | N: 25             | 50  | 50  | 100 |
| 21,22 | Sample 1       | 25                | 50  | 50  | 100 |
| 23,24 | Sample 2       | 25                | 50  | 50  | 100 |
| etc.  |                |                   |     |     |     |

(All volumes are given in µL)

- 1) Labelling of the assay tubes should be done in the following order:
  - 1, 2 total counts (TC),
  - 3, 4 zero calibrator (B<sub>0</sub>),
  - 5 - 18 duplicates of calibrators (F to L),
  - 19, 20 duplicates of control (N),
  - 21, 22 etc. duplicates of samples.
- 2) Add 25 µL of reagent A (assay buffer) to tubes 3 and 4.  
Add 25 µL of reagents F - L (calibrators) to tubes 5 - 18:
  - 5, 6 calibrator F
  - 7, 8 calibrator G, etc.
- 3) Add 25 µL of control (reagent N) to tubes 19 and 20.
- 4) Add 25 µL of the first sample to tubes 21 and 22, etc.
- 5) Add 50 µL reagent R (capture-antibody) beginning with tube 3.
- 6) Add 50 µL reagent S (spec. antibody) beginning with tube 3.
- 7) Add 100 µL reagent C (tracer) to all tubes.
- 8) Remove tubes 1 and 2 (total counts) or seal with a stopper.
- 9) Shake the tubes on a shaking device (350 rpm) overnight (at least 15 h) at room temperature.
- 10) Decant or aspirate the liquid (except tubes 1 and 2!) completely.  
Take care that the coating of the tubes remains intact.
- 11) Add 500 µL of assay buffer A to the tubes (except tubes 1 and 2!).
- 12) Decant or aspirate the liquid (see step 10).
- 13) Count the radioactivity of all tubes in a Gamma Counter for 1 to 3 min.

## CALCULATION OF RESULTS

### Establishing the Calibration Curve

| Calibrator | B <sub>0</sub> | F | G | H | I | J  | K  | L  |
|------------|----------------|---|---|---|---|----|----|----|
| ng/mL      | 0              | 1 | 2 | 4 | 8 | 16 | 32 | 64 |

1. Calculate the average counts of each pair of tubes. This equals to the binding value B.
2. The mean of the zero calibrator (tubes 3 and 4) equals B<sub>0</sub>.
3. Calculate the percent bound (% B/B<sub>0</sub>) by dividing the B values by B<sub>0</sub>:  

$$\%B/B_0 = B/B_0 \times 100\%$$
4. Plot % B/B<sub>0</sub> versus the calibrator concentrations on a semi-logarithmic or logit-log paper respectively. For convenience, it is recommended to use computer assisted data reduction programs.

Exemplary version, do not use to perform assays

5. Calculate the 'percentage bound of the zero calibrator':  
Average of tubes 3, 4 ( $B_0$ ) divided by average of tubes 1, 2 (total counts) multiplied by 100%.  
It should be  $\%B_0/TC > 20\%$ .

**Evaluation of sample concentrations:**

Read the concentration value (abscissa) corresponding to the  $\%B/B_0$  of the sample as in the example given below:

|  |          |
|--|----------|
| average counts of zero calibrator ( $B_0$ ): | 8648 cpm |
| average counts of sample tube:               | 4283 cpm |

$$\begin{aligned}\%B/B_0 &= \text{sample counts} / B_0 \times 100\% \\ &= 4283 / 8648 \times 100\% \\ &= 0.495 \times 100\% = 49.5\%\end{aligned}$$

For a 49.5% value on the y-axis (ordinate) a value of 7.28 ng/mL on the x-axis (abscissa) may be obtained.

In case of dilution of the samples, multiplication of this value determined graphically or by the aid of a computer program with the dilution factor gives the leptin concentration of the sample.

**Concentration of the control:**

The control N should fit with the labelled concentration range. The measured Leptin concentration are only valid, if the measured value of the control is in the labelled concentration range. Otherwise process analysis is required and depending on the results the measured values are accepted or not.

**EXPECTED NORMAL VALUES**

Serum leptin levels are mainly determined by body fat mass with low levels in lean individuals and high levels in obese subjects. In addition, there is a clear gender difference with higher levels in females at a given percentage body fat. Further, leptin levels are influenced by pubertal development. Any attempt, therefore, to give ranges of expected leptin levels must account for these relationships.

Various methods for the estimation of body fat are available such as calculation of body mass index (weight (kg) divided by the square of height (m)) (BMI), bioelectric impedance assessment (BIA) or total body dual energy x-ray absorptiometry (DXA). Although the accuracy of BMI with respect to reflecting true fat mass is inferior to other more sophisticated methods such as BIA or DXA, BMI provides a number of advantages:

- 1) It is independent of the regression models applied.
- 2) It is easy to determine, only weight and height measurements are required.
- 3) It is retrospectively mostly available.
- 4) It is the most precise measure during short-term changes of fat mass, e.g. during fasting.

Therefore, the following expectation ranges of serum leptin levels were referred to BMI as the major confounding independent variable and were stratified according to gender and pubertal development (45; see figures 1-8 and tables 4 - 11). After the age of 20 years, no significant age dependence was observed. These gender and age adjusted expectation ranges may be used to compare a measured leptin level at a given BMI with normal subjects to detect pathologic deviations.

The best-fit regression lines for the various subgroups are exponential curves of the form leptin =  $a \cdot e^{(b \cdot \text{BMI})}$ .

The 5th and 95th percentiles are given by the following equations:

$$\text{leptin} = a \cdot e^{(b \cdot \text{BMI} - c)} \text{ and } \text{leptin} = a \cdot e^{(b \cdot \text{BMI} + c)} \text{ respectively.}$$

## Exemplary version, do not use to perform assays

In a semi-logarithmic plot (y-axis = log leptin), these curves give straight lines. The values for a, b and c are given in table 3 according to gender and pubertal stage and also for adults. Using these values, the expectation ranges of leptin levels can be easily extended to lower or higher BMI ranges if required.

Example:

The 50th percentile for boys at Tanner stages 3 and 4 is given by the following curve:

$$\text{leptin} = 0.0181 \cdot e^{(0.2067 \cdot \text{BMI})}$$

The 5th percentile is given by:  $\text{leptin} = 0.0181 \cdot e^{(0.2067 \cdot \text{BMI} - 1.1919)}$

and the 95th percentile is given by:  $\text{leptin} = 0.0181 \cdot e^{(0.2067 \cdot \text{BMI} + 1.1919)}$

In a semi-logarithmic plot, these lines are parallel with an equal distance to the 50th percentile.

### Calculation of standard deviation scores (SDS; Z-scores)

A convenient method to detect any deviation of a measured leptin level from the corresponding reference range is to calculate its standard deviation score by relating the leptin level at the patient's BMI to the average leptin value of the corresponding sex and age group and expressing its deviation by the x-fold standard deviation. This method may be considered as normalization to the normal reference cohort. Thus, the leptin values can be adjusted for BMI, gender and pubertal stage/age (i.e., the influence of gender, age and BMI are removed) and may be pooled for further analysis.

Accounting for the logarithmic distribution of leptin levels, the leptin SDS can be calculated by the following equation:

$$\text{leptin SDS} = (\ln(\text{leptin}) - \ln(a) - b \cdot \text{BMI}) / d$$

In this equation, ln represents the natural logarithm (referring to the basis e). The constants a, b and d are given in table 3 according to gender and pubertal stage/age.

Example:

A boy at Tanner stage 3, BMI = 25 kg/m<sup>2</sup>, measured leptin concentration = 5 ng/mL.

$$\begin{aligned} \text{leptin SDS} &= (\ln(5) - \ln(0.0181) - 0.2067 \cdot 25) / 0.6850 = \\ &= (1.6094 - (-4.0118) - 5.1675) / 0.6850 = 0.66 \end{aligned}$$

### Estimation of optimal dilution of samples

Because serum leptin levels vary widely over several orders of magnitude, depending mainly on the body fat mass, adequate dilution might be a prerequisite for precise measurement. Therefore, samples should be diluted such, that the leptin concentration is near this value in order to take maximum advantage of the precision of the kit. The reference ranges provide a useful tool for a good estimation of the expected leptin value according to BMI, gender and age.

Example:

Adult woman, BMI = 45 kg/m<sup>2</sup> (130 kgs, 1.70 m height). From the reference range for adult women, the average leptin level at a BMI of 45 kg/m<sup>2</sup> is approximately 224 ng/mL. The optimal dilution would be 1:30.

**Appendix / Anhang**

**Table 3:** Constants a, b, c and d for calculation of leptin reference ranges and leptin SDS based on BMI. Groups of normal healthy individuals were stratified according to gender and pubertal stage/age. TS= Tanner stage, n= number of subjects, a,b,c and d = constants as defined in the text.

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

**Figure 1:** Reference ranges of human serum levels referring to BMI: Girls Tanner stage 1 & 2 (see text for details).

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**Figure 2:** Reference ranges of human serum levels referring to BMI: Boys Tanner stage 1 & 2 (see text for details).

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

**Figure 3:** Reference ranges of human serum levels referring to BMI: Girls Tanner stage 3 & 4 (see text for details).

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**Figure 4:** Reference ranges of human serum levels referring to BMI: Boys Tanner stage 3 & 4 (see text for details).

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

**Figure 5:** Reference ranges of human serum levels referring to BMI: Girls Tanner stage 5 (see text for details)

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**Figure 6:** Reference ranges of human serum levels referring to BMI: Boys Tanner stage 5 (see text for details).

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

**Figure 7:** Reference ranges of human serum levels referring to BMI: Adult women (see text for details)

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
**Figure 8:** Reference ranges of human serum levels referring to BMI: Adult men (see text for details).

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**Figure 9:** Serial dilution of leptin calibrator (filled circle) and different sera with varying concentrations of leptin. The leptin level of the undiluted calibrator was 16 ng/mL. Means of duplicate measurements are shown.

## LITREATURE / LITERATUR

### Literature

1. Zhang Y, Proenca R, Maffei M, Barone M, Leopold L, Friedman JM. 1994 Positional cloning of the mouse obese gene and its human homologue. *Nature*. 372:425-432.
2. Halaas JL, Gajiwala KS, Maffei M, et al. 1995 Weight-reducing effects of the plasma protein encoded by the obese gene. *Science*. 269:543-546.
3. MacDougald OA, Hwang CS, Fan H, Lane MD. 1995 Regulated expression of the obese gene product (leptin) in white adipose tissue and 3T3-L1 adipocytes. *Proc Natl Acad Sci U S A*. 92:9034-9037.
4. Rentsch J, Chiesi M. 1996 Regulation of ob gene mRNA levels in cultured adipocytes. *FEBS Lett*. 379:55-59.
5. Wabitsch M, Jensen PB, Blum WF, et al. 1996 Insulin and cortisol promote leptin production in cultured human fat cells. *Diabetes*. 45:1435-1438.
6. Stephens TW, Basinski M, Bristow PK, et al. 1995 The role of neuropeptide Y in the antiobesity action of the obese gene product. *Nature*. 377:530-532.
7. Campfield LA, Smith FJ, Guisez Y, Devos R, Burn P. 1995 Recombinant mouse OB protein: evidence for a peripheral signal linking adiposity and central neural networks. *Science*. 269:546-549.
8. Pelleymounter MA, Cullen MJ, Baker MB, et al. 1995 Effects of the obese gene product on body weight regulation in ob/ob mice. *Science*. 269:540-543.
9. Levin N, Nelson C, Gurney A, Vandlen R, de-Sauvage F. 1996 Decreased food intake does not completely account for adiposity reduction after ob protein infusion. *Proc Natl Acad Sci U S A*. 93:1726-1730.
10. Tartaglia LA, Dembski M, Weng X, et al. 1995 Identification and expression cloning of a leptin receptor, OB-R. *Cell*. 83:1263-1271.
11. Chen H, Charlat O, Tartaglia LA, et al. 1996 Evidence that the diabetes gene encodes the leptin receptor: identification of a mutation in the leptin receptor gene in db/db mice. *Cell*. 84:491-495.
12. Lee GH, Proenca R, Montez JM, et al. 1996 Abnormal splicing of the leptin receptor in diabetic mice. *Nature*. 379:632-635.
13. Lynn RB, Cao GY, Considine RV, Hyde TM, Caro JF. 1996 Autoradiographic localization of leptin binding in the choroid plexus of ob/ob and db/db mice. *Biochem Biophys Res Commun*. 219:884-889.
14. Considine RV, Considine EL, Williams CJ, Hyde TM, Caro JF. 1996 The hypothalamic leptin receptor in humans: identification of incidental sequence polymorphisms and absence of the db/db mouse and fa/fa rat mutations. *Diabetes*. 45:992-994.
15. Sinha MK, Opentanova I, Ohannesian JP, et al. 1996 Evidence of free and bound leptin in human circulation. *J Clin Invest*. 98:1277-1282.
16. Ahima RS, Prabakaran D, Mantzoros C, et al. 1996 Role of leptin in the neuroendocrine response to fasting. *Nature*. 382:250-252.
17. Chehab FF, Lim ME, Lu R. 1996 Correction of the sterility defect in homozygous obese female mice by treatment with the human recombinant leptin. *Nat Genet*. 12:318-320.
18. Schwartz MW, Baskin DG, Bukowski TR, et al. 1996 Specificity of leptin action on elevated blood glucose levels and hypothalamic neuropeptide Y gene expression in ob/ob mice. *Diabetes*. 45:531-535.
19. Schwartz MW, Seeley RJ, Campfield LA, Burn P, Baskin DG. 1996 Identification of targets of leptin action in rat hypothalamus. *J Clin Invest*. 98:1101-1106.
20. Campfield LA, Smith FJ, Burn P. 1996 The OB protein (leptin) pathway - a link between adipose tissue mass and central neural networks. *Horm Metab Res*. 28:619-632.
21. Rohner-Jeanrenaud F, Cusin I, Sainsbury A, Zakrzewska KE, Jeanrenaud B. 1996 The loop system between neuropeptide Y and leptin in normal and obese rodents. *Horm Metab Res*. 28:642-648.
22. Chan YY, Steiner RA, Clifton DK. 1996 Regulation of hypothalamic neuropeptide-Y neurons by growth hormone in the rat. *Endocrinol*. 137:1319-1325.
23. Pierroz DD, Catzeflis C, Aebi AC, Rivier JE, Aubert ML. 1996 Chronic administration of neuropeptide Y into the lateral ventricle inhibits both the pituitary-testicular axis and growth hormone and insulin-like growth factor I secretion in intact adult male rats. *Endocrinol*. 137:3-12.
24. Frederich RC, Hamann A, Anderson S, Lollmann B, Lowell BB, Flier JS. 1995 Leptin levels reflect body lipid content in mice: evidence for diet-induced resistance to leptin action. *Nat Med*. 1:1311-1314.
25. Maffei M, Halaas J, Ravussin E, et al. 1995 Leptin levels in human and rodent: measurement of plasma leptin and ob RNA in obese and weight-reduced subjects. *Nat Med*. 1:1155-1161.
26. Considine RV, Sinha MK, Heiman ML, et al. 1996 Serum immunoreactive-leptin concentrations in normal-weight and obese humans. *N Engl J Med*. 334:292-295.
27. Kolaczynski JW, Considine RV, Ohannesian J, et al. 1996 Responses of leptin to short-term fasting and refeeding in humans: A link with keto-genesis but not ketones themselves. *Diabetes*. 45:1511-1515.
28. Harris RB, Ramsay TG, Smith SR, Bruch RC. 1996 Early and late stimulation of ob mRNA expression in meal-fed and overfed rats. *J Clin Invest*. 97:2020-2026.
29. Kolaczynski JW, Ohannesian J, Considine RV, Marco C, Caro JF. 1996 Response of leptin to short-term and prolonged overfeeding in humans. *J Clin Endocrinol Metab*. 91:4162-4165.

Exemplary version, do not use to perform assays

30. Saladin R, De-Vos P, Guerre-Millo M, et al. 1995 Transient increase in obese gene expression after food intake or insulin administration. *Nature*. 377:527-529.
31. Cusin I, Sainsbury A, Doyle P, Rohner-Jeanrenaud F, Jeanrenaud B. 1995 The ob gene and insulin. A relationship leading to clues to the understanding of obesity. *Diabetes*. 44:1467-1470.
32. Kolaczynski JW, Nyce MR, Considine RV, et al. 1996 Acute and chronic effects of insulin on leptin production in humans: studies in vivo and in vitro. *Diabetes*. 45:699-701.
33. Malström R, Taskinen M-R, Karonen S-L, Yki-Järvinen H. 1996 Insulin increases plasma leptin concentrations in normal subjects and patients with NIDDM. *Diabetologia*. 39:993-996.
34. De-Vos P, Saladin R, Auwerx J, Staels B. 1995 Induction of ob gene expression by corticosteroids is accompanied by body weight loss and reduced food intake. *J Biol Chem*. 270:15958-15961.
35. Slieker LJ, Sloop KW, Surface PL, et al. 1996 Regulation of expression of ob mRNA and protein by glucocorticoids and cAMP. *J Biol Chem*. 271:5301-5304.
36. Miell JP, Englaro P, Blum WF. 1996 Dexamethasone induces an acute and sustained rise in circulating leptin levels in normal human subjects. *Horm Metab Res*. 28:704-707.
37. Sinha MK, Ohannesian JP, Heiman ML, et al. 1996 Nocturnal rise of leptin in lean, obese, and non-insulin-dependent diabetes mellitus subjects. *J Clin Invest*. 97:1344-1347.
38. Havel PJ, Kasim-Karakas S, Dubuc GR, Mueller W, Phinney SD. 1996 Gender difference in plasma leptin concentrations. *Nature Med*. 2:949-950.
39. Rosenbaum M, Nicolson M, Hirsch J, et al. 1996 Effects of gender, body composition, and menopause on plasma concentrations of leptin. *J Clin Endocrinol Metab*. 81:3424-3427.
40. Hassink SG, Sheslow DV, de Lancey E, Opentanova I, Considine RV, Caro JF 1996 Serum leptin in children with obesity: relationship to gender and development. *Pediatrics*. 98:201-203.
41. Scholz GH, Englaro P, Thiele I, et al. 1996 Dissociation of serum leptin concentration and body fat content during long term dietary intervention in obese individuals. *Horm Metab Res*. 28:718-723.
42. Caro JF, Kolaczynski JW, Nyce MR, et al. 1996 Decreased cerebrospinal-fluid/serum leptin ratio in obesity: a possible mechanism for leptin resistance. *Lancet*. 348:159-161.
43. Robinson CJ, Gaines-Das R, Woollacott D, et al. 2001 The first international standard for human leptin and the first international standard for mouse leptin: comparison of candidate preparations by in vitro bioassays and immunoassays. *J Molecular Endocrinol*. 27: 69-76.
44. Adresse NIBSC: Blanche Lane, South Mimms, Potters Bar, Hertfordshire EN6 3QG, Great Britain.
45. Blum WF, Juul A; Reference ranges of serum leptin, s.318-326. In: *Leptin -the voice of adipose tissue*, Blum WF et al, eds., Johann Ambrosius Verlag, Heidelberg, 1997.
46. Holtkamp K, Blum WF, et al., 2006, Elevated Physical Activity and Low Leptin Levels Co-Occur in Patients with Anorexia Nervosa, *The Journal of Clinical Endocrinology Metabolism* 88(11) 5169-5174
47. Guebre-Egziabher, Bernhard J., et al. 2005, Adiponectin in chronic kidney disease is related more to metabolic disturbances than to decline in renal function, *Nephrol Dial Transplant* 20: 129-134

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**SUMMARY – Mediagnost LEPTIN RIA-CT R44**

|  |                                   |
|--|-----------------------------------|
| <b>Reagent Preparation:</b>                        | <b>Reconstitution:</b>            |
| <b>Capture Antibody (R)</b>                        | in <b>7 mL Assay Buffer (A)</b>   |
| <b>Specific Antibody (S)</b>                       | in <b>7 mL Assay Buffer (A)</b>   |
| <b>Tracer (C)</b>                                  | in <b>13 mL Assay Buffer (A)</b>  |
| <b>Control (N)</b>                                 | in <b>750 µL Assay Buffer (A)</b> |
| <b>Plasma/Serum Samples</b> can be used undiluted. |                                   |

**Proposed Assay Procedure for Double Determination**

| <b>Addition of Reagent [µL]</b> |                          |   |                                   |                                    |                     |
|---------------------------------|--------------------------|---|-----------------------------------|------------------------------------|---------------------|
| <b>Nr. of Tubes</b>             | <b>Contents of Tubes</b> | <b>A Assay Buffer<br/>F-L Calibrators<br/>N Control<br/>Samples</b> | <b>R<br/>Capture<br/>Antibody</b> | <b>S<br/>Specific<br/>Antibody</b> | <b>C<br/>Tracer</b> |
| <b>1,2</b>                      | TC                       | –   | –                                 | –                                  | 100                 |
| <b>3,4</b>                      | B <sub>0</sub>           | <b>A:</b> 25  | 50                                | 50                                 | 100                 |
| <b>5-18</b>                     | Calibrators              | <b>F-L:</b> 25  | 50                                | 50                                 | 100                 |
| <b>19,20</b>                    | Control                  | <b>N:</b> 25  | 50                                | 50                                 | 100                 |
| <b>21,22</b>                    | Sample 1                 | 25  | 50                                | 50                                 | 100                 |
| <b>23,24 (etc)</b>              | Sample 2 (etc)           | 25  | 50                                | 50                                 | 100                 |

**Tubes Nr.:1,2 remove until counting the activity**

**Incubation, overnight, at least. 15 hours at RT, 350 rpm**

Decant or aspirate the liquid completely. Take care that the coating of the tubes remains intact.

Add **500 µL** of reagent **A (Assay Buffer)** to the tubes.

Aspirate the liquid completely (see above).

Count the radioactivity of **all** tubes.

Exemplary version do not use to perform assays



**REF R44**  **International Test Description**

International Test Description

|              |                        |  |
|--------------|------------------------|--|
| <b>F - L</b> |                        |  |
| <b>N</b>     | <b>Rec in</b> 750 µL A |  |
| <b>R</b>     | <b>Rec in</b> 7 mL A   |  |
| <b>S</b>     | <b>Rec in</b> 7 mL A   |  |
| <b>C</b>     | <b>Rec in</b> 13 mL A  |  |

|                    |       |
|--------------------|-------|
| <b>SPE</b>         | 25 µL |
| <b>°C</b> 20-25 °C |       |

| T       |                      |       | A     | F-L   | R     | S     | C      |
|---------|----------------------|-------|-------|-------|-------|-------|--------|
| 1 / 2   | <b>C = TC</b>        | -     | -     | -     | -     | -     | 100 µL |
| 3 / 4   | <b>B<sub>0</sub></b> | -     | 25 µL | -     | 50 µL | 50 µL | 100 µL |
| 5 / 6   | <b>F (1 ng/mL)</b>   | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 7 / 8   | <b>G (2 ng/mL)</b>   | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 9 / 10  | <b>H (4 ng/mL)</b>   | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 11 / 12 | <b>I (8 ng/mL)</b>   | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 13 / 14 | <b>J (16 ng/mL)</b>  | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 15 / 16 | <b>K (32 ng/mL)</b>  | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 17 / 18 | <b>L (64 ng/mL)</b>  | -     | -     | 25 µL | 50 µL | 50 µL | 100 µL |
| 19 / 20 | <b>N</b>             | 25 µL | -     | -     | 50 µL | 50 µL | 100 µL |
| 21 ...  | <b>SPE</b>           | 25 µL | -     | -     | 50 µL | 50 µL | 100 µL |

 >15 h  20-25  350 rpm

|  |          |
|--|----------|
| <b>ASP</b> <b>ALL</b> <b>NOT</b> <b>TC</b> |          |
| 500 µL                                     | <b>A</b> |
| <b>ASP</b>                                 |          |
| <b>MEASURE</b>                             |          |