


Leptin Sensitive ELISA

Enzyme Immunoassay for Sensitive Quantitative Determination of

human Leptin
English


For Research Use Only.
Not for use in diagnostic procedures.













REF **E077**



Gesellschaft für Forschung und Herstellung von Diagnostika GmbH

 : Aspenhastr. 25 • D-72770 Reutlingen / Germany
Phone: + 49 - (0) 7121 51484-0 • Fax: + 49 - (0) 7121 51484-10
E-mail: contact@mediagnost.de • <http://www.mediagnost.de>

	Expiry date/ Verfallsdatum/ Date de péremption/ Data di scadenza/ Fecha de caducidad/ Data de validade/ Uiterste gebruiksdatum/ Udløbsdato/ Bäst före-datum/ Termin ważności/ Lejárati idő/ Čas expirácie/ Doba expirace/ Срок на годност/ Αεγυμίσκυυρπäv/ Ημερομηνία λήξης/ Data de expirare/ Rok uporabe/ Viimeinen käyttöpäivä
	Consider instructions for use/ Bitte Gebrauchsanweisung beachten/ Consultez la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consulte las instrucciones de uso/ Respeitar as instruções de utilização/ A.u.b de gebruiksaanwijzing volgen/ Se brugsanvisningen/ Läs anvisningarna före användning/ Proszę przeczytać instrukcję obsługi/ Vegye figyelembe a használati utasításban foglaltakat/ Postupujte podľa pokynov na použitie/ Dodržujte návod k použití/ Моля, спазвайте инструкцията за употреба/ Palun järgige kasutusjuhendit/ Λάβετε υπόψη σας τις οδηγίες χρήσης/ Vă rugăm să respectați instrucțiunile de utilizare/ Upoštečajte navodila za uporabo/ Lue käyttöohje huolellisesti!
	Lot-Batch Number/ Charge-Chargennummer/ Lot-Code du lot/ Lotto-Numero di lotto/ Lote-Código de lote/ Lote-Código do Lote/Lot-Partijnummer/ Lot-Batchkode/ Partisatskod/ Numer serii/ Tétel-sarzs szám/ Číslo šarže/ Číslo šarže/ Партиден номер/Partii number/ Παρτίδα-αριθμός παρτίδας/ Lot-număr lot/ Številka serije/ Erä
	Manufactured by/ Hergestellt von/ Fabriqué par/ Prodotto da/ Fabricado por/ Fabricado por/ Vervaardigd door/ Fabrikation af/ Tillverkad av/ Wyprodukowane przez/ Gyártotta/ Vyrobene/ Vyrobeno v/ Производител/ Τοῦτjα/ Κατασκευάζεται από/ Produs de/ Proizvajalec/ Valmistaja
	Catalogue Number/ Bestellnummer/ Numéro de référence/ Numero di riferimento/ Número de referencia/ Número de Referência/ Referentienummer/ Referencenummer/ Bestellningsnummer/ Numer katalogowy/ Rendelési szám/ Katalógové číslo/ Objednací číslo/ Καταλογен номер/ Tellimisnumber/ Αρ. παραγγελίας/Număr comandă/ Številka naročila/ Viite tai tilausnumero
	Store at between/ Lagerung bei zwischen/ Conserver à entre/ Conservare a tra/ Conservar a temp. entre/ Armazena entre/ Bewaar bij tussen/ Opbevares mellem/ Förvaras vid/ Przechowywać w/ Tárolási tartomány/ Skladujte v rozsahu / Skladujte v rozmezí/ Температурно ограничение/ Säilitada temperatuuridel/ Φύλαξη σε θερμοκρασία/ Depozitare între/ Skladiščenje med/ Säilytys x-y Celsiusasteen lämpötilassa
	Contains sufficient for x tests/ Inhalt ausreichend für x Tests/ Contient suffisant pour x tests/ Contenuto sufficiente per x test/ Contenido suficiente para x pruebas/ Conteúdo suficiente para x testes/ Bevat voldoende voor x bepalingen/ Ineholder tilstrækkeligt til x prøver/ Innehållet räcker till x analyser/ Zawartość na x testów/ Tartalma x teszt elvégzésére elegendő/ Obsahuje materiál pre x testov/ Obsah dostatočný pro x testů/ Съдържание достатъчно за x тестове/ Sisust jätkub x katse jaoks/ Το περιεχόμενο επαρκεί για x δοκιμές/ Conținut suficient pentru x teste/ Vsebina zadostuje za x preizkusov/ Sisältö riittää x testille
	Incubation time/ Inkubationszeit/ Temps d'incubation/ Tempo d'incubazione/ Tiempo de incubación/ Tempo de incubação/ incubatietijd/ Inkubationstid/ inkubationstid/ Czas inkubacji/ Inkubációs idő/ Inkubačná lehota/ Inkubační doba/ Инкубационен период/ Inkubatsiooniaeg/ Χρόνος επώασης/ Timp de incubare/ Inkubacijska doba/ inkubaatioaika
	Incubate at/ Inkubation bei/ Incuber à/ Incubare a/ incubar a/ Incubar a/ incubatietemperatuur/ Inkubation ved/ inkubation vid/ Inkubacja przy/ Inkubáció hőmérséklete/ Inkubácia pri/ Inkubace při/ Инкубира се при/ Inkubatsioon temperatuuril/ Επώαση στους/ Incubare la/ Inkubacija pri/ Inkubaatiolämpötila
	Mix tubes with a Vortex mixer/ Mix Röhrchen mit Vortex Mixer/ Mélanger à l'aide d'un vortex/ Miscelare la provetta con agitatore Vortex/ Tubos de mezcla con mezclador de vortex/ Misturar os tubos com um agitador Vortex/ buisjes mengen met een Vortex/ Blanderør med Vortex-mixer/ Blanda rören med en vortexblandare/ Miksowanie rurek w mikserze Vortex/ Csővecskék keverése örvénykeverővel/ Premiešat' pomocou prístroja Vortex/ Promíchat pomocí přístroje Vortex/ Разбъркване на епруветките с миксер Vortex/ Segada torukesi Vortexi mikseriga/ Αναμίξτε τους σωληνίσκους με αναδευτήρα Vortex/ Amestecați erubetele cu ajutorul unui agitator vortex/ Mešanje cevčic z mešalnikom Vortex/ Sekoita putket Vortex sekoittajalla
MTP	Mikrotiterplate/ Mikrotiterplatte/ plaque de microtitrage/ Piastra di microtitolazione/ Placa de microtitulación/ Placa de Microtitulação/ Microtiterplaat/ Mikrotiterplade/ mikrotiterplatta/ microtiterplaat/ Płytko microtiter/ Mikrotiter lap/ Mikrotitračná podložka/ Mikrotitrační podložka/ Микротитърна плака/ Mikrotiterplaat/ Τρυβλίο μικροτιτλοδότησης/ Microplacă/ Mikrotitrska plošča/ Mikrotitrauslevy
Rec in	Reconstitute in/ Rekonstituieren in/ Reconstituer dans/ Ricostituire nel/ Reconstituir en/ Reconstituir em/ Rekonstituieren in/ Rekonstituier i/ Rekonstituera/ Rekonstituować w/ Helyreállítás/ Znovu připravit za/ Znovu připravit za/ Разтваряне в/ Moodustada uuesti/ Ανασυστήστε σε/ Reconstituire în/ Predelava v/ Rekonstituoi
SPE	Sample/ Probe/ Echantillon/ Campione/ Muestra/ Amostra/ Monster/ Prøve/ prov/ Próbká/ Minta/ Vzorka/ Vzorek/ Проба/ Proov/ Δείγμα/ Probă/ Vzorec/ Näyte
DET	Antibody and Enzyme Conjugate/ Antikörper- und Enzymkonjugat/ anticorps conjugué et conjugué enzymatique/ Coniugato di anticorpo ed enzima/ Conjugado de anticuerpos y enzimas/ Conjugado Anticorpo-Enzima/ antilichaam- en enzymconjugaat/ Antistoffer og enzym-konjugat/ antikropp- och enzymkonjugat (antikropp och enzym, konjugat)/ Koniugat antyciał i enzymów/ Antitest és enzim páros/ Protílátkový a enzymatický konjugát/ Protílátkový a enzymatický konjugát/ Антитяло и ензим конюгат/ Antikehad ja ensüümi konjugaat/ Σύμπλοκο αντισώματος-ενζύμου/ Compuși din anticorpi și enzime/ Antitelesa in konjugat encima/ Vasta-aine ja entsümi konjugaatti
DIL	Dilution Buffer/ Verdünnungspuffer/ Tampon de dilution/ Tampone di diluizione/ Tampón de dilución/ / Tampão de diluição/ Verdunningsbuffer/ / Fortyndingsbuffer/ Utspänningsbuffert / Bufor rozcieńczający/ / Hígító puffer/ Riediaci pufor/ Ředící pufr / Буфер за разреждане/ Lahjenduspuhver/ Ρυθμιστικό διάλυμα αραιώσης / Tampon de diluare/ Pufer za redčenje/ Laimennuspuskuri

X:X	Dilute / Verdünnen / Diluer / Diluire / Diluir / Diluir / Verdunnen / Fortyndes / Späd / Rozcieńczanie / Hígítás / Riedit' / Ředit / Разреждане / Lahjendada / Αραιώστε / Diluaŋi / Razredčiti / Laimennetaan
CAL A-E	Calibrator X/ Kalibrator X/ calibreur X/ calibratore X/ calibrador X/ calibrador X/ kalibrator X/ kalibrator X/ kalibrator X/ kalibrator X/ kalibrátor X/ kalibrátor X/ kalibrátor X/ kalibrátor X/ калибратор X/ kalibraator X/ Βαθμονομητής X/ calibrator X/ kalibrator X/ kalibraattori X
CTR1 / CTR2	Control X/ Kontrolle X/ Contrôle X/ controllo X/ control X/ Controle X/ controle X/ Kontrol X/ Kontroll X/ kontrolne X/ Ellenőrző X/ Kontrolné X/ Kontrolní X/ Контролен X/ Kontroll X/ ελέγχου X/ control X/ Kontrolni X/ Kontrolli X
WB	Washing Buffer Concentrate/ Waschpufferkonzentrat/ Tampon de lavage conc./ Tampone di lavaggio concentrato/ Tampón de lavado concentrado/ Tampão de Lavagem Concentrado/ wasbuffer, geconcentreerd/ Vaskebufferkonzentrat/ Vaskebufferkonzentrat/ tvättbuffertkonzentrat/ Bufor płukania koncentrat/ Mosópufer koncentrátum/ Koncentrát vymývacieho pufru/ Концентрат на промивен буфер/ Pesurpuhvi kontsentraat/ Συμπύκνωμα ρυθμιστικού διαλύματος έκπλυσης/ Concentrat pentru tamponul de spălare/ Koncentrat izpiralnega pufra/ Pesuliuositiivist
WB 1:20	Washing Buffer/ Waschpuffer/ Tampon de lavage/ Tampone di lavaggio/ Tampón de lavado/Tampão de Lavagem/ wasbuffer/ Vaskebuffer/ tvättbuffert/ Bufor płukania/ Mosópufer/ Vymývací pufer/ Vymývací pufr/ Промивен буфер/ Pesurpuhver/ Ρυθμιστικό διάλυμα έκπλυσης/ Tampon pentru spălare /Izpiralni pufer/ Pesuliuos
S	Substrate/ Substrat/ Substrat/ Substrato/ Substrato/ Substrato/ Substrato/ Substrat/ Substrat/ Substrat/ Substrat/ Substrat/ Substrat/ Substrat/ Szubsztrátum/ Substrát/ Substrát/ Substraat/ Υπόστρωμα/ Substrat/ Substrat/ Substraattiliuos
STP	Stop Solution/ Stopplösung/ Stop Solution/ Soluzione di stop/ Stop Solución/ Solução Stop/ stopoplossing/ Stopopløsning/ Stopplösning/ Stop roztwór/ Megállító oldat/ Roztok na ukončenie/ Roztok pro ukončení/ Стопирец разтвор/ Stopp-lahus/ Διάλυμα διακοπής/ Soluție de oprire/ Stop raztopina/ Pysäytysliuos
TAPE	Cover Plate with sealing tape/ Platte abkleben/ Recouvrir la microplaque avec bande adhésive/ Coprire la piastra con nastro adesivo/ Cubrir la placa con una cinta adhesiva/ Cobrir a Placa com fita adesiva/ plaatje met tape afdekken/ Afdækningsplade med tape/ maskera platta/ Odkleić płytke/ Tányér leragasztása/ Oblepit' podložku lepiacou páskou/ Olepit podložku lepicí páskou/ Плака с лента за запечатване/ Katta plaat isoleerikleerlindiga/ Κολλήστε το πλακίδιο με κολλητική ταινία/ Aoperiti placa cu o bandă adezivă/ Prelepiti ploščo/ Peitä mikrotitrauslevy oheisella teipillä
MEASURE	Measure plate within 30 min at 450 nm (Referencefilter ≥590 nm)/ Ausmessung innerhalb von 30 min bei 450 nm (Referenzfilter ≥ 590 nm)/ Mesure l'absorbance en l'espace de 30 min à 450 nm avec ≥590 nm longueur d'onde pour référence/ Misurazione entro 30 min. a 450 nm (filtro di riferimento ≥ 590 nm)/ Medición de la placa dentro de los siguientes 30 min a 450 nm (filtro de referencia ≥ 590nm)/ Medir a placa dentro de 30 min a 450 nm (Filtro de referência ≥ 590 nm)/ Binnen 30 minuten bij 450 nm meten (referentiefilter ≥ 590 nm)./ Mål plade i løbet af 30 min ved 450 nm (referencefilter ≥590 nm)/ Mät inom 30 min vid 450 nm (referensfilter ≥ 590 nm)./ Pomiar w ciągu 30 min przy 450 nm (filtr odniesienia ≥ 590 nm)/ Ki mérés 30 percen belül 450 nm-nél (referenciaszűrő ≥ 590 nm)/ Merat 30 minut pri 450 nm (Referenčných filtrov ≥590 nm)/ Měřit 30 minut při 450 nm (Referenční filtr ≥ 590 nm)/ Отчитане в рамките на 30 min при 450 nm (референтен филтър ≥ 590 nm)/ Mõõtmise 30 min jooksul 450 nm korral (võrdlusfilter ≥ 590 nm)/ Μέτρηση εντός 30 min στα 450 nm (φίλτρο αναφοράς ≥ 590 nm)/ Măsurare în decurs de 30 min la 450 nm (filtru de referință ≥ 590 nm)/ Izmerite ploščico v 30 min pri 450 nm (referenčni filter ≥590nm)/ Mittaa 30 minuutin aikana 450 nm:ssä (referenssi suodatin ≥ 590 nm)
Literature	Literature/ Literatur/ Bibliographie/ Letterario/ Bibliografia/ Literatura documentação/ literatuur/ Litteratur/ litteratur/ Literatura/ Irodalom/ Literatúra/ Literatura/ Литература// Kirjandus/ Βιβλιογραφία/ Bibliografie/ literatura/ Lähdeluettelo
International Test Description	International test description/ internationale Testanleitung/ description internationale de test/ Istruzioni per il test internazionali/ Descripción de ensayo internacional/ Descrição internacional do teste/ internationale testbeschrijving/ internationell testbeskrivning/ Opis testu międzynarodowego/ nemzetközi teszt-útmutató/ Medzinárodný návod k testu/ Mezinárodní návod k testu/ rahvusvaheline katse kirjeldus/ Διεθνείς οδηγίες για εργαστηριακές δοκιμές/ instrucțiuni internaționale pentru testare/ mednarodna navodila za preizkus/ Kansainvälinen käyttöohje
End	in all required wells/ in allen benötigten Vertiefungen/ dans tous les godets requis/ in tutti i pozzetti richiesti/ en todos los pozos requeridos/ em todos os tubos necessários/ in alle nodige putjes/ i alle nødvendige brønde/ i alla nödvändiga brunnar/ we wszystkich potrzebnych wgłębieniach/ minden szükséges forrásban/ vo všetkých potrebných miestach/ ve všech potřebných místech/ във всички необходими ямки/ kõigis vajalikes süvendites/ σε όλες τις απαραίτητες κοιλότητες/ în toate cavitățile necesare/ v vseh zahtevanih vdolbinah/ kaikkiin tarvittaviin mikrotitrauslevyn syvennyyksiin

Read entire protocol before use!

For Research Use Only.

Not for use in diagnostic procedures.

CAUTION: Not for human or animal therapeutic or diagnostic use.

For in vitro use only.

For professional use only.

Table of Contents

1	INTENDED USE	5
2	INTRODUCTION	5
3	PRINCIPLE.....	6
4	WARNINGS AND PRECAUTIONS	7
5	SAMPLES.....	8
6	REAGENTS PROVIDED	9
7	TECHNICAL NOTES.....	10
8	ASSAY PROCEDURE.....	11
9	CALCULATION OF RESULTS	12
10	EXEMPLARY VALUES	14
11	PERFORMANCE CHARACTERISTICS	19
12	COMPARISON STUDIES.....	20
13	LITERATURE	21
14	International Assay Description	23
15	ASSAY PROCEDURE.....	24

Instructions for use

Leptin Sensitive ELISA E077	96 Determinations
Principle of the test	Enzyme-linked Immunoassay
Duration (incubation period)	1.75 h
Antibodies	monoclonal antibodies, ready for use
Buffer	Ready for use and 20fold concentrate
Calibrators	5 single Calibrators: 0.05 – 5 µg/L, lyophilized recombinant human Leptin
Reference material	International Standard WHO/ NIBSC 97/594, recombinant Leptin
Assay Range	0.01 – 50 µg/L
Control	2 Controls, freeze-dried
Sample	human serum / plasma, other human body fluids e.g. urine, saliva Precise measurement also in lean subjects
Required sample volume	25 µL
Sample dilution	1:10
Analytical sensitivity	ø 0.01 µg/L
Intra- / Interassay Variance	ø < 10 %
Exemplary values	Blum WF, Juul A; Reference ranges of serum leptin, In: Leptin- the voice of adipose tissue, Blum WF et al. eds. Johann Ambrosius Verlag, Heidelberg, 1997

1 INTENDED USE

Measurement of human Leptin in human serum, plasma and in other human body fluids e.g. urine, saliva, breast milk, for research use only!

2 INTRODUCTION

Leptin, the product of the ob gene (1,2), is a recently discovered single-chain proteohormone with a molecular weight of 16 kD, which is thought to play a key role in the regulation of body weight. Its amino acid sequence exhibits no major homologies with other proteins (1). Leptin is almost exclusively produced by differentiated adipocytes (3-5). It acts on the central nervous system, in particular the hypothalamus, thereby suppressing food intake and stimulating energy expenditure (2,6-9). Leptin receptors - alternatively spliced forms exist that differ in length - belong to the cytokine class I receptor family (10-12). They are found ubiquitously in the body (10, 11, 13, 14) indicating a general role of leptin, which is currently not fully understood. A circulating form of the leptin receptor exists which acts as one of several leptin binding proteins (15). Besides its metabolic effects, leptin was shown to have a strong influence on a number of endocrine axes. In male mice, it blunted the starvation-induced marked decline of LH, testosterone, thyroxine and the increase of ACTH and corticosterone. In female mice, leptin prevented the starvation-induced delay in ovulation (16). Ob/ob mice, which are leptin deficient due to an ob gene mutation, are infertile. This defect could be corrected by administration of leptin, but not through weight loss due to fasting (17), suggesting that leptin is pivotal for reproductive functions.

All these actions may, at least in part, be explained by the suppressive effect of leptin on neuropeptide Y (NPY) expression and secretion by neurons in the arcuate nucleus (6,18,19). NPY is a strong stimulator of appetite (20,21) and is known to be involved in the regulation of

various pituitary hormones, e.g. suppression of GH through stimulation of somatostatin (22,23), suppression of gonadotropins (23) or stimulation of the pituitary-adrenal axis (21).

The most important variable that determines circulating leptin levels is body fat mass (24-26). Obviously, under conditions of regular eating cycles, leptin reflects the proportion of adipose tissue (27) showing an exponential relationship (37). This constitutive synthesis of leptin is modulated by a number of non-hormonal and hormonal variables. Stimulators in both rodents and humans are overfeeding (28,29), insulin (3,5,30-33) and glucocorticoids (5,34-36). Suppression has been shown for fasting (27), cAMP and beta-3-adrenoceptor agonists (35). From these findings it becomes clear that leptin is an integral component of various metabolic and endocrine feedback loops (38).

Serum leptin levels show a moderate **circadian variation** with a peak during the night at about 2 a.m. (37). The leptin values at this time are about 30 to 100 % higher than the levels measured in the morning or early afternoon. This variation together with the influence of food intake needs to be taken into account, when blood samples are collected.

This ELISA kit is suited for measuring human leptin in serum or plasma, other biological fluids (e.g. saliva, urine or cerebrospinal fluid), and conditioned adipocyte culture media for scientific purposes.

3 PRINCIPLE

The Mediagnost Sensitive ELISA for Leptin E077 is a so-called Sandwich-Assay using two specific and high affinity antibodies. The Leptin in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-Leptin-Antibody binds in turn to the immobilised Leptin. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction the turn of the colour will be catalysed quantitatively depending on the Leptin-level of the samples.

4 WARNINGS AND PRECAUTIONS

For research and professional use only not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided. A Material Safety Data sheet is available on request. Do not use obviously damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Controls CTR1 and CTR2**

Source human serum for the Controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore, all components and specimens should be treated as potentially infectious.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

Reagents CAL A-E, DET, DIL, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate S

The TMB-Substrate S contains 3,3',5,5' Tetramethylbenzidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution STP

The Stop solution contains 0.2 M acid sulphur acid (H₂SO₄)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1 General first aid procedures

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

Serum and plasma samples.

With this immunoassay kit, human leptin can also be determined for the measurement of leptin with low concentrations, e.g. in samples in body fluids other than serum (e.g. urine, liquor, breast milk) as well as in conditioned adipocyte culture media **for research use only**, for **scientific purposes**.

5.1 Blood Sample Types

Significant deviation of Leptin levels in corresponding Serum- or EDTA-Plasma samples were not found. If commercially available tubes for Citrate-plasma are used for preparation, samples will be diluted, resulting Leptin values will be therefore reduced.

5.2 Blood Samples

The blood samples for serum preparation should be gained according to standardized venipuncture procedure. Hemolytic reactions have to be avoided.

The human serum or plasma samples should be obtained in the morning or in the early afternoon (2:00 p.m.) if the nutritional status is normal.

Leptin levels show a circadian variation with a peak during the night at about 2 a.m. (37). This variation together with the influence of food intake needs to be taken into account, when blood samples are collected.

Required Serum/ Plasma sample volume: 25 µL

5.3 Serum- / Plasma- Sample stability

In firmly closable sample vials

- Storage at 20-25°C: 2 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles max. 5

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized, 5 Freezing-Thawing showed no effect on samples.

5.4 Interference

Hemoglobin, triglyceride and bilirubin in the sample do not interfere to a concentration of **1 mg/mL**, **100 mg/mL** and **100 µg/mL**. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

5.5 Sample dilution

Samples have to be diluted in Dilution Buffer **DIL**.

For **Serum and Plasma** samples we recommend a dilution of 1:10.

Suggestion for dilution protocol:

Dilute e.g. 225 µL Dilution Buffer **DIL** in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add 25 µL Serum-or Plasma (Dilution 1:10). After mixing use 2 x 100 µL from the dilution in the assay.


If Leptin levels over 50 ng/mL are expected, the sample should be diluted higher in Dilution Buffer **DIL**, e.g. 1:20.

The hLeptin concentrations may be **completely different in cell culture supernatants and body fluids of human origin other than serum**. The optimal dilution must be found out by the customer. E.g. in samples of urine or breast milk very low levels were determined. Furthermore also this kind of samples must be diluted at least 1:5. If samples are used undiluted, the reaction will be interfered and the values are reduced.

6 REAGENTS PROVIDED

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the calibration curve.

MTP	Microtiter plate , ready for use, coated with mouse-anti-Leptin-antibody. Wells are separately breakable.	(8x12) wells
CAL A-E	Calibrators , lyophilized, (recombinant human Leptin), concentrations are given on vial labels and on quality certificate ng/mL.	5 x 1 mL
CTR1	Control 1 , lyophilised, (human serum) concentration is given on quality certificate	1 x 250 µL
CTR2	Control 2 , lyophilised, (human serum) concentration is given on quality certificate	1 x 250 µL
DET	Antibody-HRP-Conjugate , ready for use, mouse-anti-hLeptin-antibody biotinylated + streptavidin horseradish peroxidase conjugate	1 x 12 mL
DIL	Dilution Buffer , ready for use	1 x 60 mL
WB	Washing Buffer , 20-fold concentrated solution	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised tetramethylbencidine.	1 x 12 mL
STP	Stop Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x
	Instructions for use	1 x
--	Quality Certificate	1 x

6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm
- Polyethylene PE / Polypropylene PP tubes for diluting of samples

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**. Store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components**: Calibrators **A-E** and Controls **CTR1** and **CTR2** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C.

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Calibrators **A – E** and Controls **CTR1** and **CTR2** are reconstituted with the Dilution Buffer **DIL**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest. After reconstitution dilute the Controls **CTR1** and **CTR2** with the Dilution Buffer **DIL** in the same ratio (1:10) as the sample.

Assay Procedure

When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1**, **CTR2** and the samples should be pipette as fast as possible. To avoid distortions due to differences in incubation times, Antibody-HRP-Conjugate **DET** as well as the succeeding Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**.

All determinations (Blank, Calibrators **A-E**, Controls **CTR1**, **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive—storage and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending **350 rpm**. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8 ASSAY PROCEDURE

All determinations (Calibrators, Control **CTR1 & CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Calibrators, Controls and the samples should be pipetted as fast as possible (e.g., <15 minutes).

All incubations have to be conducted at room temperature (20-25°C)!

To avoid distortions due to differences in incubation times, Antibody-POD-Conjugate **DET** as well as the following Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as the Substrate **S**.

- 1) Add **100 µL** Dilution Buffer **DIL** to the first wells (blank). Subsequently, add 100 µL of each Calibrator or diluted Control (CTR1&CTR2) or diluted Sample to the following wells.
- 2) Cover the wells with sealing tape and incubate the plate for **1 hour** shaking with **350 rpm**.
- 3) After incubation aspirate the contents of the wells into a disinfectant (possible theoretically risk of infection!) and wash the wells **3 times** with **300 µL** of **Washing Buffer WB** / well respectively. The washing buffer WB should incubate at least for 15 seconds/cycle
- 4) Pipette **100 µL** of the **Antibody-HRP-Conjugate DET** in each well and incubate **30 minutes** shaking with **350 rpm**.
- 5) After incubation wash the wells 3 times with Washing Buffer WB as described in step 3
- 6) Pipette **100 µL of the Substrate S** in each well.
- 7) Incubate the plate for **15 minutes in the dark at room temperature (20 - 25°C)**.
- 8) Stop the reaction by adding **100 µL of Stop Solution STP**.
- 9) Measure the colour reaction within 30 minutes at **450 nm** (reference filter **≥590 nm**).

9 CALCULATION OF RESULTS

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of Calibrator E should be above 1.00. Samples, which yield higher absorbance values than Calibrator E, should be re-tested with a dilution.

9.1 Establishing of the calibration curve

Calibrator	A	B	C	D	E
ng/mL	0.05	0.5	1.5	3.5	5

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples, controls and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The leptin concentration of the diluted controls **CTR1** and **CTR2** and diluted **samples** is obtained from the calibration curve. The **concentration** of the undiluted samples / controls in **ng/mL** can be calculated **by multiplication** with the respective dilution factor.

9.2 Example of a typical Calibration Curve

The following Calibration Curve is for demonstration only and cannot be used in place of data generation at the time of assay.




Figure 1: Exemplary Standard Curve with a polynomial 2nd degree as curve fit.

Exemplary calculation of the leptin concentration of the 1:10 diluted sample:

Measured extinction of your sample 0.293

Measured extinction of the blank 0.015

Your **software program** will calculate the leptin concentration of the diluted sample automatically by using the difference of sample and blank (0.278) for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2nd degree).

In this exemplary case the following equation is solved by the program to calculate the Leptin concentration in the sample:

$$0.278 = -0.0089658 + 0.64743x - 0.034025x^2$$
$$0.4524 = x$$

if the dilution factor (1:10) is taken into account the leptin concentration of the undiluted sample is

$$0.4524 \text{ ng/mL} \times 10 = 4.4524 \text{ ng/mL}$$

10 EXEMPLARY VALUES

Serum leptin levels are mainly determined by body fat mass with low levels in lean individuals and high levels in obese subjects. In addition, there is a clear gender difference with higher levels in females at a given percentage body fat. Further, leptin levels are influenced by pubertal development. Any attempt, therefore, to give ranges of expected leptin levels must account for these relationships.

The following exemplary serum leptin levels were referred to BMI as the major confounding independent variable and were stratified according to gender and pubertal development (45; see figures 2-9 and tables 2-9). After the age of 20 years, no significant age dependence was observed. These. The best-fit regression lines for the various subgroups are exponential curves of the form:

$$Leptin = a \cdot e^{(b \cdot BMI)}$$

Table 1 Constants a, b, c and d for calculation of leptin ranges and leptin SDS based on BMI. Groups of normal healthy individuals were stratified according to gender and pubertal stage/age. TS= Tanner stage, n= number of subjects, a,b,c and d = constants as defined in the text (see Chapter Expected Normal Values).

Cohort	n	a	b	c	d
Males:					
TS 1&2	136	0.0146	0.2706	0.8821	0.5379
TS 3&4	50	0.0181	0.2067	1.1919	0.6850
TS 5	112	0.0316	0.1462	1.0821	0.6558
Adults	380	0.0130	0.2200	1.1053	0.6740
Females					
TS 1&2	136	0.0422	0.2499	0.7849	0.4786
TS 3&4	43	0.0543	0.2357	0.5745	0.3379
TS 5	157	0.2550	0.1508	0.7053	0.4301
Adults	587	0.3042	0.1467	0.8548	0.5212

Table 2 Girls Tanner stages 1 and 2

BMI (kg/m ²)	Percentile (µg/L)				
	1	5	50	95	99
11	0.22	0.30	0.66	1.45	1.99
12	0.28	0.39	0.85	1.86	2.56
13	0.36	0.50	1.09	2.38	3.29
14	0.46	0.64	1.40	3.06	4.22
15	0.60	0.82	1.79	3.93	5.42
16	0.76	1.05	2.30	5.04	6.96
17	0.98	1.35	2.95	6.47	8.93
18	1.25	1.73	3.79	8.31	11.5
19	1.61	2.22	4.87	10.7	14.7
20	2.07	2.85	6.25	13.7	18.9
21	2.65	3.66	8.03	17.6	24.3
22	3.41	4.70	10.3	22.6	31.2
23	4.37	6.03	13.2	29.0	40.0
24	5.62	7.75	17.0	37.2	51.4
25	7.21	9.95	21.8	47.8	65.9
26	9.26	12.8	28.0	61.4	84.7
27	11.9	16.4	35.9	78.8	109.0
28	15.3	21.1	46.1	101.0	140.0
29	19.6	27.0	59.2	130.0	
30	15.2	34.7	76.1		
31	32.3	44.6	97.7		
32	41.5	57.2	125.		
33	53.2	73.4			
34	68.4	94.3			
35	87.8	121.			
36	113				
37	145				

Table 3 Boys Tanner stages 1 and 2

BMI	Percentile (µg/L)				
	1	5	50	95	99
11	0.08	0.12	0.29	0.69	0.99
12	0.01	0.16	0.38	0.91	1.30
13	0.14	0.20	0.49	1.19	1.71
14	0.19	0.26	0.65	1.56	2.24
15	0.24	0.35	0.85	2.04	2.93
16	0.32	0.46	1.11	2.68	3.84
17	0.41	0.60	1.45	3.51	5.04
18	0.55	0.79	1.90	4.60	6.60
19	0.72	1.03	2.50	6.03	8.66
20	0.94	1.35	3.27	7.90	11.3
21	1.24	1.77	4.29	10.4	14.9
22	1.62	2.33	5.62	13.6	19.5
23	2.12	3.05	7.37	17.8	25.5
24	2.78	3.99	9.66	23.3	33.5
25	3.65	5.24	12.7	30.6	43.9
26	7.78	6.87	16.9	40.1	57.5
27	6.27	9.0	21.7	52.5	75.4
28	8.22	11.8	28.5	68.9	98.8
29	10.7	15.5	37.4	90.3	129.0
30	14.1	20.3	48.9	118.0	
31	18.5	26.6	64.2		
32	24.3	34.8	84.1		
33	31.8	45.6	110.0		
34	41.7	59.8	144.0		
35	54.6	78.4			
36	71.6	102.0			
37	93.9	134.0			
38	123.0				




Figure 2 Reference ranges of human serum levels referring to BMI: Girls Tanner stage 1 & 2 (see text for details)




Figure 3 Reference ranges of human serum levels referring to BMI: Boys Tanner stage 1 & 2 (see text for details)

Table 4 Girls Tanner stages 3 and 4.

BMI (kg/m ²)	Percentile (µg/L)				
	1	5	50	95	99
11	0.32	0.41	0.73	1.29	1.63
12	0.41	0.52	0.92	1.63	2.06
13	0.52	0.66	1.16	2.07	2.61
14	0.65	0.83	1.47	2.61	3.31
15	0.83	1.05	1.87	3.31	4.19
16	1.05	1.33	2.36	4.19	5.30
17	1.33	1.68	2.99	5.30	6.71
18	1.68	2.13	3.78	6.71	8.49
19	2.13	2.69	4.79	8.5	10.8
20	2.69	3.41	6.06	10.7	13.6
21	3.41	4.31	7.67	13.61	17.2
22	4.32	5.46	9.71	17.2	21.8
23	5.46	6.91	12.3	21.8	27.6
24	6.91	8.75	15.6	27.6	34.9
25	8.75	11.1	19.7	34.9	44.2
26	11.1	14.0	24.9	44.2	56.0
27	14.0	17.7	31.6	56.0	70.9
28	17.8	22.5	39.9	70.9	89.7
29	22.5	28.4	50.5	89.7	114.0
30	28.4	36.0	63.9	114.0	144.0
31	36.0	45.6	80.9	144.0	
32	45.6	57.7	102.0	144.0	
33	57.7	73.0	130.0		
34	73.0	92.4			
35	92.4	117.0			
36	117.0	148.0			
37	148.0				

Table 5 Boys Tanner stage 3 and 4

BMI (kg/m ²)	Percentile (µg/L)				
	1	5	50	95	99
11	0.03	0.05	0.18	0.58	0.94
12	0.04	0.07	0.22	0.71	1.16
13	0.49	0.08	0.27	0.88	1.43
14	0.06	0.10	0.33	1.08	1.75
15	0.07	0.12	0.40	1.32	2.16
16	0.09	0.15	0.49	1.63	2.65
17	0.11	0.18	0.61	2.00	3.26
18	0.14	0.23	0.75	2.46	4.01
19	0.17	0.28	0.92	3.03	4.93
20	0.21	0.34	1.13	3.72	6.06
21	0.26	0.42	1.39	4.58	7.46
22	0.32	0.52	1.71	5.63	9.17
23	0.39	0.64	2.10	6.92	11.3
24	0.48	0.78	2.58	8.51	13.9
25	0.59	0.96	3.18	10.5	17.0
26	0.73	1.19	3.91	12.9	21.0
27	0.89	1.46	4.80	15.8	25.8
28	1.10	1.79	5.90	19.4	31.7
29	1.35	2.20	7.26	23.9	39.0
30	1.66	2.71	8.93	29.4	48.0
31	2.05	3.33	11.0	36.2	58.9
32	2.51	4.09	13.5	44.5	72.4
33	3.09	5.04	16.6	54.7	89.1
34	3.80	6.20	20.4	67.2	109.0
35	4.68	7.62	25.1	82.6	134.0
36	5.75	9.37	30.9	101.0	
37	7.07	11.5	37.	124.0	
38	8.7	14.2	46.7		
39	10.7	17.4	57.4		
40	13.1	21.4	70.5		




Figure 4 Reference ranges of human serum levels referring to BMI: Girls Tanner stage 3 & 4 (see text for details).




Figure 5 Reference ranges of human serum levels referring to BMI: Boys Tanner stage 3 & 4 (see text for details).

Table 6 Girls Tanner stage 5

Percentile (µg/L)					
BMI (kg/m ²)	1	5	50	95	99
11	0.50	0.66	1.34	2.71	3.62
12	0.58	0.77	1.56	3.15	4.21
13	0.67	0.89	1.81	3.67	4.89
14	0.78	1.04	2.11	4.26	5.69
15	0.91	1.21	2.45	4.96	6.62
16	1.05	1.41	2.85	5.76	7.70
17	1.22	1.64	3.31	6.70	8.95
18	1.42	1.90	3.85	7.79	10.4
19	1.66	2.21	4.48	9.06	12.1
20	1.93	2.57	5.20	10.5	14.1
21	2.24	2.99	6.05	12.3	16.4
22	2.60	3.48	7.03	14.2	19.0
23	3.03	4.04	8.18	16.6	22.1
24	3.52	4.70	9.51	19.3	25.7
25	4.09	5.46	11.0	22.4	29.9
26	4.76	6.35	12.9	26.0	34.8
27	5.53	7.39	15.0	30.3	40.4
28	6.43	8.59	17.39	35.2	47.0
29	7.48	9.99	20.2	40.9	54.7
30	8.70	11.6	23.5	47.6	63.5
31	10.1	13.5	27.3	55.3	73.9
32	11.8	15.7	31.8	64.4	85.9
33	13.7	18.3	37.0	74.9	99.9
34	15.9	21.2	43.0	87.0	116.0
35	18.5	24.7	50.0	101.0	135.0
36	21.5	28.7	58.1	118.0	
37	25.0	33.4	67.6	137.0	
38	29.1	38.8	78.6		
39	33.8	45.1	91.4		
40	39.4	52.5	106.0		

Table 7 Boys Tanner stage 5

Percentile (µg/L)					
BMI (kg/m ²)	1	5	50	95	99
11	0.03	0.05	0.16	0.47	0.73
12	0.04	0.06	0.18	0.54	0.84
13	0.05	0.07	0.21	0.62	0.97
14	0.05	0.08	0.24	0.72	1.12
15	0.06	0.10	0.28	0.84	1.30
16	0.07	0.11	0.33	0.97	1.51
17	0.08	0.13	0.38	1.12	1.74
18	0.1	0.15	0.44	1.3	2.02
19	0.11	0.17	0.51	1.50	2.34
20	0.13	0.2	0.59	1.74	2.7
21	0.15	0.23	0.68	2.01	3.13
22	0.17	0.27	0.79	2.33	3.62
23	0.20	0.31	0.91	2.69	4.19
24	0.23	0.36	1.05	3.12	4.85
25	0.27	0.41	1.22	3.61	5.62
26	0.31	0.48	1.41	4.17	6.5
27	0.36	0.55	1.63	4.83	7.52
28	0.41	0.64	1.89	5.59	8.71
29	0.48	0.74	2.19	6.47	10.1
30	0.55	0.86	2.54	7.49	11.7
31	0.64	1.00	2.94	8.67	13.5
32	0.74	1.15	3.4	10.0	15.6
33	0.86	1.33	3.94	11.6	18.1
34	0.99	1.54	4.55	13.4	20.9
35	1.15	1.79	5.27	15.6	24.2
36	1.33	2.07	6.10	18.0	28.1
37	1.54	2.39	7.06	20.8	32.5
38	1.78	2.77	8.17	24.1	37.6
39	2.06	3.21	9.46	27.9	43.5
40	2.38	3.71	10.9	32.3	50.3




Figure 6 Reference ranges of human serum levels referring to BMI: Girls Tanner stage 5 (see text for details)




Figure 7 Reference ranges of human serum levels referring to BMI: Boys Tanner stage 5 (see text for details).

Table 8 Adult women

BMI (kg/m ²)	Percentile (µg/L)				
	1	5	50	95	99
11	0.46	0.65	1.53	3.59	5.10
12	0.53	0.75	1.77	4.16	5.90
13	0.61	0.87	2.05	4.82	6.83
14	0.71	1.01	2.37	5.58	7.91
15	7.82	1.17	2.75	6.46	9.17
16	0.95	1.35	3.18	7.48	10.61
17	1.10	1.57	3.68	8.66	12.3
18	1.28	1.81	4.27	10.0	14.2
19	1.48	2.10	4.94	11.6	16.5
20	1.71	2.43	5.72	13.4	19.1
21	1.99	2.82	6.62	15.6	22.1
22	2.30	3.26	7.67	18.0	25.6
23	2.66	3.78	8.88	20.9	29.3
24	3.08	4.38	10.3	24.2	34.3
25	3.57	5.07	11.9	28.0	39.7
26	4.13	5.87	13.8	32.4	46.0
27	4.79	6.79	16.0	37.5	53.3
28	5.54	7.87	18.5	43.5	61.7
29	6.42	9.11	21.4	50.4	71.5
30	7.43	10.6	24.8	58.3	82.8
31	8.61	12.2	28.7	67.5	95.8
32	9.97	14.1	33.3	78.2	111.0
33	11.5	16.4	38.5	90.5	129.0
34	13.4	19.0	44.6	105.0	149.0
35	15.5	22.0	51.6	121.0	
36	17.9	25.4	59.8	141.0	
37	20.8	29.5	69.3		
38	24.0	34.1	80.2		
39	27.8	39.5	92.9		
40	32.2	45.7	108.0		

Table 9 Adult men

BMI (kg/m ²)	Percentile (µg/L)				
	1	5	50	95	99
11	0.03	0.05	0.15	0.44	0.69
12	0.04	0.06	0.18	0.55	0.87
13	0.05	0.08	0.23	0.69	1.08
14	0.06	0.09	0.28	0.85	1.34
15	0.07	0.12	0.35	1.06	1.67
16	0.09	0.15	0.44	1.33	2.09
17	0.12	0.18	0.55	1.65	2.60
18	0.14	0.23	0.68	2.06	3.24
19	0.18	0.28	0.85	2.57	4.04
20	0.22	0.35	1.06	3.20	5.03
21	0.23	0.44	1.32	3.98	6.27
22	0.35	0.54	1.64	4.97	7.81
23	0.43	0.78	2.05	6.19	9.73
24	0.54	0.85	2.55	7.71	12.1
25	0.67	1.05	3.18	9.61	15.1
26	0.83	1.31	3.96	12.0	18.8
27	1.04	1.64	4.94	14.9	23.5
28	1.30	2.04	6.15	18.6	29.2
29	1.61	2.54	7.67	23.2	36.4
30	2.01	3.16	9.56	28.9	45.4
31	2.51	3.94	11.9	36.0	56.6
32	3.12	4.91	14.8	44.9	70.5
33	3.89	6.12	18.5	55.8	87.8
34	4.85	7.63	23.0	69.6	109.0
35	6.04	9.51	28.7	86.7	136.0
36	7.53	11.8	35.8	108.0	
37	9.38	14.8	44.6	135.0	
38	11.7	18.4	55.5		
39	14.6	22.9	69.2		
40	18.2	28.6	86.2		




Figure 8 Reference ranges of human serum levels referring to BMI: Adult women (see text for details)




Figure 9 Reference ranges of human serum levels referring to BMI: Adult men (see text for details).

11 PERFORMANCE CHARACTERISTICS

11.1 Sensitivity

The analytical sensitivity of the assay yields 0.01 ng/mL (2x SD of zero standards in 12-fold determinations).

11.2 Precision

Intra- and Inter-Assay Variance

The inter- and intra assay coefficients of variability are below 10%. Exemplary determinations are shown in table 10 and 11.

Table 10 Inter-Assay-Variation

	Mean value (ng/mL)	Standard Deviation (ng/mL)	VC (%)
Sample 1	2.04	0.147	7.2
Sample 2	6.93	0.423	6.1
Sample 3	14.86	1.11	7.5

Table 11 Intra-Assay-Variation

	Mean value (ng/mL)	Standard Deviation (ng/mL)	VC (%)
Sample 1	22.44	1.28	4.35
Sample 2	4.1	0.108	2.63

11.3 Linearity

Dilutions of 1:5 up to 1:320 have been tested with four serum samples and the amount of Leptin measured is shown in table 12. The data show that dilutions of 1:5 to 1:320 seem to give reasonable results (no deviation of more than 20% of the mean value of all dilutions tested). A dilution of 1:10 of serum samples is recommended in general.

Table 12 Dilution linearity of 1:5 to 1:320

Dilution	Sample A ng/mL	Sample B ng/mL	Sample C ng/mL	Sample D ng/mL
1:5	10	26		23
1:8	10	25	39	25
1:10	11	26	36	26
1:20	11	28	37	28
1:40	11	29	39	29
1:80	12	29	41	30
1:160	10	30	40	31
1:320	11	27	37	30

11.4 Recovery

Recombinant Leptin (NIBSC) was added to Dilution Buffer (DIL) and serum samples. The leptin content of the so enriched samples was measured and relative recovery was calculated. Results are shown in table 13.

Table 13 Relative recovery of NIBSC leptin in serum samples [%]

Leptin added [µg/L]	5	10	15
Sample 1	109	108	108
Sample 2	108	104	109
Sample 3	105	102	85

11.5 Interference

Interference of bilirubin, hemoglobin and triglycerides was tested by adding different amounts of these substances to human serum containing Leptin. For comparison the same amount of buffer without any substance was also added to the serum. Table 14 demonstrates that neither bilirubin nor triglycerides or haemoglobin exert any influence on the measurement of Leptin in human serum.

Table 14 Interference of physiological substances

	Triglyceride 100 mg/mL	Bilirubin 100 µg/mL	Hemoglobin 1 mg/mL
Serum 1	95	101	94
Serum 2	107	105	105
Serum 3	111	101	101

11.6 Species Cross-Reactivity

The cross-reactivities of several commercially available animal sera were tested in this assay. For this purpose, the sera in various dilutions were used as a sample in the Mediagnost Leptin ELISA E077.

NO SIGNAL was detected in the sera of the following species

Horse, cattle, chicken, rabbit, dog, guinea pig, sheep, mouse, goat, donkey, rat, cat, pig.

12 COMPARISON STUDIES

The correlation between ELISA E077 and the Mediagnost LEP-R40 radioimmunoassay is shown in Figure 10.




Figure 10 Assay comparison
– Mediagnost Radioimmunoassay LEP-R40 and ELISA E077

13 LITERATURE







LITERATURE

1. Zhang Y, Proenca R, Maffei M, Barone M, Leopold L, Friedman JM. 1994 Positional cloning of the mouse obese gene and its human homologue. *Nature*. 372:425-432.
2. Halaas JL, Gajiwala KS, Maffei M, et al. 1995 Weight-reducing effects of the plasma protein encoded by the obese gene. *Science*. 269:543-546.
3. MacDougald OA, Hwang CS, Fan H, Lane MD. 1995 Regulated expression of the obese gene product (leptin) in white adipose tissue and 3T3-L1 adipocytes. *Proc Natl Acad Sci U S A*. 92:9034-9037.
4. Rentsch J, Chiesi M. 1996 Regulation of ob gene mRNA levels in cultured adipocytes. *FEBS Lett*. 379:55-59.
5. Wabitsch M, Jensen PB, Blum WF, et al. 1996 Insulin and cortisol promote leptin production in cultured human fat cells. *Diabetes*. 45:1435-1438.
6. Stephens TW, Basinski M, Bristow PK, et al. 1995 The role of neuropeptide Y in the antiobesity action of the obese gene product. *Nature*. 377:530-532.
7. Campfield LA, Smith FJ, Guisez Y, Devos R, Burn P. 1995 Recombinant mouse OB protein: evidence for a peripheral signal linking adiposity and central neural networks. *Science*. 269:546-549.
8. Pellemounter MA, Cullen MJ, Baker MB, et al. 1995 Effects of the obese gene product on body weight regulation in ob/ob mice. *Science*. 269:540-543.
9. Levin N, Nelson C, Gurney A, Vandlen R, de-Sauvage F. 1996 Decreased food intake does not completely account for adiposity reduction after ob protein infusion. *Proc Natl Acad Sci U S A*. 93:1726-1730.
10. Tartaglia LA, Dembski M, Weng X, et al. 1995 Identification and expression cloning of a leptin receptor, OB-R. *Cell*. 83:1263-1271.
11. Chen H, Charlat O, Tartaglia LA, et al. 1996 Evidence that the diabetes gene encodes the leptin receptor: identification of a mutation in the leptin receptor gene in db/db mice. *Cell*. 84:491-495.
12. Lee GH, Proenca R, Montez JM, et al. 1996 Abnormal splicing of the leptin receptor in diabetic mice. *Nature*. 379:632-635.
13. Lynn RB, Cao GY, Considine RV, Hyde TM, Caro JF. 1996 Autoradiographic localization of leptin binding in the choroid plexus of ob/ob and db/db mice. *Biochem Biophys Res Commun*. 219:884-889.
14. Considine RV, Considine EL, Williams CJ, Hyde TM, Caro JF. 1996 The hypothalamic leptin receptor in humans: identification of incidental sequence polymorphisms and absence of the db/db mouse and fa/fa rat mutations. *Diabetes*. 45:992-994.
15. Sinha MK, Opentanova I, Ohannesian JP, et al. 1996 Evidence of free and bound leptin in human circulation. *J Clin Invest*. 98:1277-1282.
16. Ahima RS, Prabakaran D, Mantzoros C, et al. 1996 Role of leptin in the neuroendocrine response to fasting. *Nature*. 382:250-252.
17. Chehab FF, Lim ME, Lu R. 1996 Correction of the sterility defect in homozygous obese female mice by treatment with the human recombinant leptin. *Nat Genet*. 12:318-320.
18. Schwartz MW, Baskin DG, Bukowski TR, et al. 1996 Specificity of leptin action on elevated blood glucose levels and hypothalamic neuropeptide Y gene expression in ob/ob mice. *Diabetes*. 45:531-535.
19. Schwartz MW, Seeley RJ, Campfield LA, Burn P, Baskin DG. 1996 Identification of targets of leptin action in rat hypothalamus. *J Clin Invest*. 98:1101-1106.
20. Campfield LA, Smith FJ, Burn P. 1996 The OB protein (leptin) pathway - a link between adipose tissue mass and central neural networks. *Horm Metab Res*. 28:619-632.
21. Rohner-Jeanrenaud F, Cusin I, Sainsbury A, Zakrzewska KE, Jeanrenaud B. 1996 The loop system between neuropeptide Y and leptin in normal and obese rodents. *Horm Metab Res*. 28:642-648.
22. Chan YY, Steiner RA, Clifton DK. 1996 Regulation of hypothalamic neuropeptide-Y neurons by growth hormone in the rat. *Endocrinol*. 137:1319-1325.
23. Pierroz DD, Catzeflis C, Aebi AC, Rivier JE, Aubert ML. 1996 Chronic administration of neuropeptide Y into the lateral ventricle inhibits both the pituitary-testicular axis and growth hormone and insulin-like growth factor I secretion in intact adult male rats. *Endocrinol*. 137:3-12.
24. Frederich RC, Hamann A, Anderson S, Lollmann B, Lowell BB, Flier JS. 1995 Leptin levels reflect body lipid content in mice: evidence for diet-induced resistance to leptin action. *Nat Med*. 1:1311-1314.
25. Maffei M, Halaas J, Ravussin E, et al. 1995 Leptin levels in human and rodent: measurement of plasma leptin and ob RNA in obese and weight-reduced subjects. *Nat Med*. 1:1155-1161.
26. Considine RV, Sinha MK, Heiman ML, et al. 1996 Serum immunoreactive-leptin concentrations in normal-weight and obese humans. *N Engl J Med*. 334:292-295.
27. Kolaczynski JW, Considine RV, Ohannesian J, et al. 1996 Responses of leptin to short-term fasting and refeeding in humans: A link with keto-genesis but not ketones themselves. *Diabetes*. 45:1511-1515.
28. Harris RB, Ramsay TG, Smith SR, Bruch RC. 1996 Early and late stimulation of ob mRNA expression in meal-fed and overfed rats. *J Clin Invest*. 97:2020-2026.
29. Kolaczynski JW, Ohannesian J, Considine RV, Marco C, Caro JF. 1996 Response of leptin to short-term and prolonged overfeeding in humans. *J Clin Endocrinol Metab*. 91:4162-4165.
30. Saladin R, De-Vos P, Guerre-Millo M, et al. 1995 Transient increase in obese gene expression after food intake or insulin administration. *Nature*. 377:527-529.
31. Cusin I, Sainsbury A, Doyle P, Rohner-Jeanrenaud F, Jeanrenaud B. 1995 The ob gene and insulin. A relationship leading to clues to the understanding of obesity. *Diabetes*. 44:1467-1470.

32. Kolaczynski JW, Nyce MR, Considine RV, et al. 1996 Acute and chronic effects of insulin on leptin production in humans: studies in vivo and in vitro. *Diabetes*. 45:699-701.
33. Malström R, Taskinen M-R, Karonen S-L, Yki-Järvinen H. 1996 Insulin increases plasma leptin concentrations in normal subjects and patients with NIDDM. *Diabetologia*. 39:993-996.
34. De-Vos P, Saladin R, Auwerx J, Staels B. 1995 Induction of ob gene expression by corticosteroids is accompanied by body weight loss and reduced food intake. *J Biol Chem*. 270:15958-15961.
35. Sliker LJ, Sloop KW, Surface PL, et al. 1996 Regulation of expression of ob mRNA and protein by glucocorticoids and cAMP. *J Biol Chem*. 271:5301-5304.
36. Miell JP, Englaro P, Blum WF. 1996 Dexamethasone induces an acute and sustained rise in circulating leptin levels in normal human subjects. *Horm Metab Res*. 28:704-707.
37. Sinha MK, Ohannesian JP, Heiman ML, et al. 1996 Nocturnal rise of leptin in lean, obese, and non-insulin-dependent diabetes mellitus subjects. *J Clin Invest*. 97:1344-1347.
38. Havel PJ, Kasim-Karakas S, Dubuc GR, Mueller W, Phinney SD. 1996 Gender difference in plasma leptin concentrations. *Nature Med*. 2:949-950.
39. Rosenbaum M, Nicolson M, Hirsch J, et al. 1996 Effects of gender, body composition, and menopause on plasma concentrations of leptin. *J Clin Endocrinol Metab*. 81:3424-3427.
40. Hassink SG, Sheslow DV, de Lancey E, Opentanova I, Considine RV, Caro JF 1996 Serum leptin in children with obesity: relationship to gender and development. *Pediatrics*. 98:201-203.
41. Scholz GH, Englaro P, Thiele I, et al. 1996 Dissociation of serum leptin concentration and body fat content during long term dietary intervention in obese individuals. *Horm Metab Res*. 28:718-723.
42. Caro JF, Kolaczynski JW, Nyce MR, et al. 1996 Decreased cerebrospinal-fluid/serum leptin ratio in obesity: a possible mechanism for leptin resistance. *Lancet*. 348:159-161.
43. Robinson CJ, Gaines-Das R, Woollacott D, et al. 2001 The first international standard for human leptin and the first international standard for mouse leptin: comparison of candidate preparations by in vitro bioassays and immunoassays. *J Molecular Endocrinol*. 27: 69-76.
44. Adresse NIBSC: Blanche Lane, South Mimms, Potters Bar, Hertfordshire EN6 3QG, Great Britain.
45. Blum WF, Juul A; Reference ranges of serum leptin,s.318-326. In:Leptin -the voice of adipose tissue, Blum WF et al, eds., Johann Ambrosius Verlag, Heidelberg, 1997.

14 International Assay Description

International Test Description

CAL A-E	Rec in 1 mL DIL	-
CTR1	Rec in 250 µL DIL	1:10 DIL → CTR1 1:10
CTR2	Rec in 250 µL DIL	1:10 DIL → CTR2 1:10
WB	-	1:20 A. dest. → WB 1:20
SPE	-	1:10 DIL → SPE 1:10
°C 20-25°C		
100 µL	DIL	A1/A2
100 µL	CAL A (0.05 ng/mL)	B1/B2
100 µL	CAL B (0.5 ng/mL)	C1/C2
100 µL	CAL C (1.5 ng/mL)	D1/D2
100 µL	CAL D (3.5 ng/mL)	E1/E2
100 µL	CAL E (5 ng/mL)	F1/F2
100 µL	CTR1 1:10	G1/G2
100 µL	CTR2 1:10	H1/H2
100 µL	SPE 1:10	
TAPE		
 1 h °C 20-25°C  350 rpm		
3x 300 µL	3x WB 1:20	
100 µL	DET	
TAPE		
 0.5 h °C 20-25°C  350 rpm		
3x 300 µL	3x WB 1:20	
100 µL	S	
 0.25 h °C 20-25°C 		
STP		
MEASURE		

15 ASSAY PROCEDURE

Preparation of reagents		Reconstitution	Dilution
CAL A-E	Calibrators	in 1 mL Dilution Buffer DIL	-
CTR1 CTR2	Controls	in 250 µL Dilution Buffer DIL	1:10 in Dilution Buffer DIL → CTR1 1:10 → CTR2 1:10
WB	Washing Buffer 20-fold conc.	-	1:20 with Aqua dest. → WB 1:20
Sample dilution in Dilution Buffer DIL, e.g. 1:10 → SPE 1:10			
Before assay procedure bring all reagents to room temperature 20-25°C .			
Assay Procedure in Double Determination:			
Pipette	Reagents	Position	
100 µL	Dilution Buffer DIL (Blank)	A1/A2	
100 µL	Calibrator A (0.05 ng/mL)	B1/B2	
100 µL	Calibrator B (0.5 ng/mL)	C1/C2	
100 µL	Calibrator C (1.5 ng/mL)	D1/D2	
100 µL	Calibrator D (3.5 ng/mL)	E1/E2	
100 µL	Calibrator E (5 ng/mL)	F1/F2	
100 µL	Control CTR1 1:10	G1/G2	
100 µL	Control CTR2 1:10	H1/H2	
100 µL	Sample SPE 1:10	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
Sample Incubation: 1 h at 20-25°C, 350 rpm			
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WB 1:20 / well	In each well	
100 µL	Antibody-HRP-Conjugate DET	In each well	
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at 20-25°C, 350 rpm			
3 x 300 µL	Aspirate the contents of the wells and wash 3 x with 300 µL each Washing Buffer WB 1:20 / well	In each well	
100 µL	Substrate S	In each well	
Incubation: 15 Minutes in the Dark at 20-25°C			
100 µL	Stop Solution STP	In each well	
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			