Anti-HAV IgM ELISA

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Table of Contents

1.	INTENDED USE	15
2.	INTRODUCTION	15
3.	ASSAY PRINCIPLE	15
4.	WARNINGS AND PRECAUTIONS	16
5.	SAMPLES	17
6.	MATERIALS	18
7.	TECHNICAL NOTES	19
8.	ASSAY PROCEDURE, Qualitative Test Procedure	20
9.	QUALITY CONTROL	21
10.	QUALITATIVE RESULTS	21
11.	QUANTITATIVE TEST PROTOCOL	21
12.	QUANTITATIVE RESULTS	21
13.	LIMITATION OF PROCEDURE	22
14.	INTERNATIONAL TEST DESCRIPTION, QUALITATIVE Assay	23
15.	QUALITATIVE ASSAY PROCEDURE	24

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ENGLISH Instructions for use

Anti-HAV ELISA E11	96 Determinations
Regulatory Status	For Research Use Only. Not for diagnostic purposes
Principle of the test	Sandwich ELISA
Duration (incubation period)	3.5 h
Conjugate Concentrate	100-fold concentrated
Antigen solution	Ready for use
Dilution Buffer and Substrate	Ready for Use
Washing buffer	20-fold concentrated
Positive and Negative Control	Ready for use
Sample	human Serum / Plasma
Required sample volume	10 µL
Sample dilution	1:2000

1. INTENDED USE

The mediagnost anti-HAV IgM ELISA E11, is an enzyme immunoassay for research use for the qualitative and quantitative detection of IgM antibodies in human serum or plasma directed against the Hepatitis A-Virus.

2. INTRODUCTION

The detection of anti-HAV IgM (IgM antibodies specific for Hepatitis A virus) in the serum or plasma of patients indicates a fresh infection with the Hepatitis-A virus. Depending on the anti-HAV IgM titre it is possible to differentiate between an acute infection (1 - 3 months after the start of clinical symptoms) and early convalescence (3 - 6 months). In rare cases (approx. 10% of all clinical cases) the anti-HAV IgM response remains positive (6 - 12 months).

3. ASSAY PRINCIPLE

The Mediagnost anti-HAV IgM ELISA is a "class-capture" enzyme immunoassay. Serum or plasma samples are diluted 1:2000 and added to the wells of a microtiter plate, which have been previously coated with an antibody directed against human IgM antibodies (μ -chain specific). All IgM antibodies present in the sample are bound during an incubation step. After washing, the Hepatitis A-virus is added, which binds only to the HAV-specific IgM antibodies. The bound antigen is detected by the addition of the conjugate (monoclonal anti-HAV antibody, peroxidase conjugated). After washing and incubation with a colorimetric substrate the reaction is terminated by addition of stop solution and the color turns to yellow. The absorbance of the coloured reaction product is measured on a microtiter plate reader. The colour intensity of the reaction corresponds to the concentration of antibodies in the sample.

4. WARNINGS AND PRECAUTIONS

For In Vitro Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro use and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human material: **PK, NK, KK**

The antigen has been inactivated with formaldehyde. Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious. The antigen AG has been inactivated with formaldehyde.

Reagents KK, VP, WP

Contain as preservative 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (<0.015%)

H317 May cause an allergic skin reaction.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P272 Contaminated work clothing should not be allowed out of the workplace.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/ attention.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P501 Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate Solution S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H335 May cause respiratory irritation.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P305+P351+ IF IN EYES: Rinse cautiously with water for several minutes.

P338 Remove contact lenses, if present and easy to do. Continue rinsing.

Stopping Solution SL

The Stopping solution contains 0.2 M acid sulphur acid (H₂SO₄)

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301+P330+ IF SWALLOWED: rinse mouth.

P331 Do NOT induce vomiting.

P305+P351+ IF IN EYES: Rinse cautiously with water for several minutes.

P338 Remove contact lenses, if present and easy to do. Continue rinsing.

P309+P310 IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1. General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5. SAMPLES

5.1. Sample type

Serum and EDTA-plasma

5.2. Specimen collections

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

5.3. Required sample volume:

10 μ L

5.4. Sample dilution

Qualitative Test: Dilution

Samples should be diluted in PE/PP-tubes. For larger sample numbers usage of a multi-stepper is recommended.

Recommended dilution for Qualitative Test: **1:2000** with Dilution Buffer **V**

Quantitative Test: Dilution

For quantitative results the samples are diluted 1:2000, 1:20000, 1:200000, 1:2000000 with Dilution Buffer **VP**.

Exemplary Dilution Protocol

1:20 Pre-Dilution

190 μ L Dilution Buffer **VP** + **10 μ L** Serum- or EDTA-Plasma

The pre-diluted sample should then be diluted as follows:

1:2000 Dilution

Add **10 μ L** of the 1:20 pre-diluted sample to **990 μ L** Dilution Buffer **VP**

1:20000 Dilution

Add **100 μ L** of the 1:2000 dilution to **900 μ L** Dilution Buffer **VP**

1:200000 Dilution

Add **100 μ L** of the 1:20000 dilution to **900 μ L** Dilution Buffer **VP**

1:2000000 Dilution

Add **100 μ L** of the 1:200000 dilution to **900 μ L** Dilution Buffer **VP**

After mixing use **100 μ L** each of the needed sample dilution per well

It is recommended to use the diluted sample within one hour.

6. MATERIALS

6.1. Materials provided

The reagents listed below are sufficient for 96 tests. The microtiter plate is divided up in 12 strips with separately breakable wells.

MTP	Microtiter plate , ready for use, coated with an antibody directed against human IgM antibodies (anti- μ -chain). Wells are separately breakable.	(8x12) wells
KK	Conjugate Concentrate , 100fold concentrated (peroxidase conjugated anti-HAV IgG) (monoclonal, mouse).	1 x 250 μL
AG	Hepatitis A Virus antigen , inactivated, ready for use.	1 x 13.5 mL
PK	Positive Control , ready for use, recalcificated human plasma, reactive for anti-HAV-IgM, Titre > 1:10000.	1 x 2 mL
NK	Negative Control , ready for use, recalcificated human plasma, not reactive for anti-HAV-IgM/IgG. Anti-HBsAg, anti-HIV and anti HCV negative.	1 x 1 mL
VP	Dilution Buffer , ready for use. Please shake before use.	1x 120 mL
WP	Washing Buffer , 20-fold concentrated solution.	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine (TMB).	1 x 12 mL
SL	Stopping Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate.	2 x
	Instructions for use	1 x
--	Quality Control Certificate	1 x

6.2. Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WP (A.dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Incubator (suitable for incubations at **37°C**)
- Vortex-mixer
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

7. TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components after initial opening is warranted for 4 weeks, store the unused strips and microtiter wells airtight together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The 1:20 diluted Washing Buffer WP is 4 weeks stable at 2-8°C.

Preparation of reagents

Bring all reagents to room temperature (20 – 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Dilution

Use the Dilution Buffer **VP** for the dilution of the Conjugate Concentrate **KK**. Prepare only the amount required. The diluted conjugate can be stored for at least one week at 4°C.

The required volume of Washing Buffer **WP** is prepared by 1:20 dilution of the provided 20fold concentrate with **Aqua dest.** Please dilute only according to daily requirements.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate Solution **S**, is photosensitive—store and incubation in the dark.

Assay Procedure

When performing the assay reagents and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, diluted Conjugate Concentrate **KK**, as well as the Substrate Solution **S** should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution **SL** should be added to the plate in the same order as Substrate Solution **S**. All determinations should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Washing

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WP** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter plate washer**, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8. ASSAY PROCEDURE, Qualitative Test Procedure

Reagent Preparation			
Before use bring all reagents to room temperature: 20-25 °C			
KK	Conjugate Concentrate	in Dilution Buffer VP (e.g. 10 µL KK + 990 µL VP)	1:100
WP	Washing Buffer	in Aqua dest. (e.g 50 mL WP + 950 mL A.dest)	1:20
-	Samples	In Dilution Buffer VP , see Page 17	1:2000
Assay Procedure			
Pipette	Reagents		Position
-	The first well (A1, blank) will be left empty		A1
100 µL	Positive Control		PK B1/C1/D1
100 µL	Negative Control		NK E1/F1
100 µL	Samples e.g. 1:2000 diluted, we recommend double determination.		according to demand
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	HAV-Antigen AG Solution		from B1 in each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	1:100 diluted Conjugate Concentrate		from B1 in each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	Substrate Solution S		in each well, also A1
Substrate S Incubation: 30 Minutes in the dark at 20-25°C			
100 µL	Stopping Solution SL		in each well, also A1
Measure the absorbance within 30 min at 450 nm (reference filter ≥ 590 nm).			

The absorbance of the blank wells are subsequently subtracted from the control and sera to obtain the absolute values

9. QUALITY CONTROL

9.1. Quality Criteria

The average absorbance values of the Positive Controls **PK**, the Negative Controls **NK** and the samples are calculated.

The Negative Control **NK** extinction value should **not exceed 0.1**. The difference between extinctions of the Positive Control **PK** and Negative Control **NK** must be at least **0.4**, otherwise the test is considered invalid.

10. QUALITATIVE RESULTS

Cut-off calculation:

The sum of: 30% of the mean absorbance value of the Positive Control **PK** + mean absorbance value of the Negative Control **NK** is used as cut-off. Interpretation of positive: Samples with extinction equal or higher than the cut-off (30% of the average of the positive control + extinction average of the negative control) are regarded as positive. Samples with an extinction of $\pm 10\%$ of the cut-off value are regarded as border line samples. In this case repeat the test.

Calculation (example)

Positive Control PK

	<u>Extinction:</u>
1. value	1.120
2. value	1.205
3. value	1.196
<u>total:</u>	<u>3.521</u>
mean value:	$3.521 / 3 = 1,174$

Negative Control NK:

1. value	0.021
2. value	0.025
<u>total:</u>	<u>0.046</u>
mean value:	$0.046 / 2 = 0.023$

The exemplary cut-off: $1.174 \times 0.3 + 0.023 = 0.375$

Exemplary Sera with extinctions > 0.375 are regarded as positive.

Exemplary Border line samples: $0.375 \times 1.1 = 0.412$ (+ 10%); $0.375 \times 0.9 = 0.337$ (-10%)

Extinctions between 0.337 and 0.412

11. QUANTITATIVE TEST PROTOCOL

To determine the serum titer for example to control the course of disease the quantitative protocol is recommended.

11.1. Samples preparation

For quantitative results the patients serum is diluted 1:2000, 1:20000, 1: 200000, 1:2000000 with dilution buffer (See page 17).

11.2. Test Preparation

Pipette from the above-mentioned serum dilutions, 100 μ L per well (we recommend double determinations). The test protocol is identical to the "Qualitative" method.

12. QUANTITATIVE RESULTS

The evaluation is carried out graphically or via an appropriate evaluation program.

A graph is used to determine the titre in the quantitative calculation of anti-HAV IgM. The graph is plotted on semi-logarithmic paper. The extinction values are plotted on the Y-axis and the serum dilution values on the X-axis.

Cut-off calculation

The sum of: 30% of the mean absorbance value of the Positive Control **PK** + mean absorbance value of the Negative Control **NK** is used as cut-off. Interpretation of positive:

Titer calculation (example)

Positive Control **PK** average: 1.240
Negative Control **NK** average: 0.026
cut-off: $1.240 \times 0.3 + 0.026 = 0.398$


Sample dilution	Extinction
1:2000	1.302
1:20000	1.098
1:200000	0.630
1:2000000	0.060

The measured extinction values of the individual serum dilutions are plotted onto the graph including the cut-off value (drawn parallel to the Y-axis). The individual dots are joined by a line; where this line crosses the limiting value line it is possible to read, from the X-axis, the anti-HAV titre. The serum titre in this example is 1: 500000.

13. LIMITATION OF PROCEDURE

The influence of the heterophilic antibodies, rheuma factors, anti-species antibodies is reduced by the assay design, but cannot be completely excluded. Higher concentrations of physiological and pharmaceutical or other substances may interfere with the measurement.

Example for the graphical analysis

14. INTERNATIONAL TEST DESCRIPTION, QUALITATIVE Assay

KK	CONJ	1:100 DILU BUF VP
WP	WASHBUF 20x	1:20 DILU A. dest.
SPE		1:2000 DILU BUF VP
°C 20-25 °C		

-	-	A1
100 µL	Control+ PK	B1/C1/D1
100 µL	Control- NK	E1/F1
100 µL	SPE	
TAPE		
🕒 1 h °C 37°C		
3x 300 µL	3x WASHBUF WP	3x 300 µL
100 µL	AG HAV-Antigen Solution	B1 → End
TAPE		
🕒 1 h °C 37°C		
3x 300 µL	3x WASHBUF WP	3x 300 µL
100 µL	CONJ KK	B1 → End
TAPE		
🕒 1 h °C 37°C		
3x 300 µL	3x WASHBUF WP	
100 µL	SUBST TMB S	
🕒 0.5 h °C 20-25°C 🌞		
100 µL	H ₂ SO ₄ SL	
MEASURE		

15. QUALITATIVE ASSAY PROCEDURE

Reagent Preparation			
Before use bring all reagents to room temperature: 20-25 °C			
KK	Conjugate Concentrate	in Dilution Buffer VP (e.g. 10 µL KK + 990 µL VP)	1:100
WP	Washing Buffer	in Aqua dest. (e.g. 50 mL WP + 950 mL A.dest)	1:20
-	Samples	In Dilution Buffer VP , see Page 17	1:2000
Assay Procedure in Double Determination			
Pipette	Reagents		Position
-	-		A1
100 µL	Positive Control		PK B1/C1/D1
100 µL	Negative Control		NK E1/F1
100 µL	Samples e.g. 1:2000 diluted		according to demand
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	HAV-Antigen AG Solution		from B1 in each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	1:100 diluted Conjugate Concentrate KK		from B1 in each well
Cover the wells with the sealing tape.			
Incubation: 1 h at 37°C			
3x 300 µL	Aspirate the contents of the wells and wash 3x with 300 µL Washing Buffer WP / well		in each well
100 µL	Substrate Solution S		in each well, also A1
Substrate S Incubation: 30 Minutes in the dark at 20-25°C			
100 µL	Stopping Solution SL		in each well, also A1
Measure the absorbance within 30 min at 450 nm (reference filter ≥ 590 nm).			