

Aspenhastraße 25
D - 72770 Reutlingen

Tel +49 - (0)7121 – 51484 - 0
Fax +49 - (0)7121 – 51484 -10
E-mail contact@mediagnost.de
Web http://www.mediagnost.de

PRODUCT SPECIFICATION MONOCLONAL ANTIBODY

REF. : M42
HYBRIDOMA CELL LINE: Clone 2B1, anti-Phosphatidylethanol (anti-PEth)
LOT: 270808,R4

PROTEIN CONCENTRATION : 1.94 mg/ml

Protein concentration has been determined by BCA-protein assay (PIERCE) with BSA as standard.

FORMULATION :

The product is supplied as a solution of antibody in phosphate-buffered saline without any preservative. The antibodies were derived from cell culture which showed no signs of contamination.

PURIFICATION :

The product has been purified by Ion-Exchange-Chromatography.

ISOTYPING: mouse IgM

Isotyping was performed by ELISA with Ab sets from DIANOVA, Germany.

PURITY: >90 %

Purity was estimated by silver stained SDS PAGE.

EXPIRY:

Store at $\leq -20^{\circ}\text{C}$, avoid repeated freeze/thaw cycles.

Product is stable for at least 1 year after receipt when stored appropriately.

DESCRIPTION:

Nissinen, A.E. et. al.: Immunological Detection of In Vitro Formed Phosphatidylethanol - An Alcohol Biomarker – With Monoclonal Antibodies; Alcohol Clin Exp Res, Vol. 32, No. 6, pp. 921-928, 2008

For in-vitro Research Use Only. Not for use in or on humans or animals. We are not responsible for violations or patent infringements which may occur with the use of our products



Geschäftsführung:
Dr. Angelika Haage, Dr. Lutz Pridzun
Sitz der Gesellschaft: Kusterdingen
Amtsgericht Stuttgart HRB 380863
UST.-ID. Nr.: DE146 893 439
Steuer-Nr.: 86115/70125

Baden-Württembergische Bank
Konto 7 477 501 403
BLZ 600 501 01
Swift/BIC SOLA DE ST
IBAN DE66 600501017477501403

Kreissparkasse Tübingen
Konto 217 288
BLZ 641 500 20
Swift/BIC SOLA DE S1 TUB
IBAN DE56 641500200000217288

Deutsche Bank 24
Konto 1 726 900
BLZ 640 700 24
Swift/BIC DEUTDE33
IBAN DE73640700240172690000