

Aspenhastraße 25
D - 72770 Reutlingen

Tel +49 - (0)7121 - 51484 - 0
Fax +49 - (0)7121 - 51484 -10

E-mail contact@mediagnost.de
<http://www.mediagnost.de>

PRODUCT CERTIFICATE

ANTI-HAV HORSERADISH PEROXIDASE CONJUGATE (K50)

Lot-Number 050511/3

PREPARATION OF CONJUGATE :

The conjugate was prepared from affinity purified anti-HAV-antibody 7E7. The monoclonal antibody 7E7 was tested in neutralization experiments. From this experiments, it may be concluded, that there is an immunodominant site, that is recognized by 7E7. In competition experiments 7E7 showed a 77 % competition versus a polyclonal serum (Ping and Lemon, 1992).

POD Activity: 270 U/mL

Activity was determined with Mediagnost POD test by usages of pyrogallol and hydrogen peroxide.

FORMULATION :

Horseradish peroxidase-labelled antibody is supplied in phosphate buffered saline containing 10 mg/ ml BSA for stabilisation.

PROTEIN CONCENTRATION : 4,54 mg/ml

Determined with BCA-protein assay reagent (Pierce, Europe), prior to addition of stabilizing proteins.

STORAGE :

Store at - 20 °C. After thawing store at 4°C. Avoid repeated freeze-thawing cycles. The conjugate remains stable for at least 6 month after delivery at appropriate storage conditions.

REFERENCE :

Li-Hua Ping and Stanley M. Lemon (1992): Antigenic Structure of Human Hepatitis A Virus Defined by Analysis of Escape Mutants Selected against Murine Monoclonal Antibodies. Journal of Virology Vol 66 No 4, 2208-2216.

Product is for in-vitro and research use only. Not for use in diagnostic procedures.

Certificate is valid without signature.

Version 4 Exemplary

Geschäftsführung:
Dr. Angelika Haage, Dr. Lutz Pridzun
Prof. Dr. Bertram Flehmig
Sitz der Gesellschaft: Kusterdingen
Amtsgericht Stuttgart HRB 380863
UST.-ID. Nr.: DE146 893 439
Steuer-Nr.: 86115/70125

Baden-Württembergische Bank
Konto 7 477 501 403
BLZ 600 501 01
Swift/BIC SOLA DE ST
IBAN DE66 600501017477501403

Kreissparkasse Tübingen
Konto 217 288
BLZ 641 500 20
Swift/BIC SOLA DE S1 TUB
IBAN DE56 641500200000217288

Deutsche Bank 24
Konto 1 726 900
BLZ 640 700 24
Swift/BIC DEUTDE33
IBAN DE73640700240172690000